MicroScan

HNID Procedural Manual

REF B1012-10B

TABLE OF CONTENTS

English	2
Français (FR)	9
Deutsch (DE)	16
Italiano (IT)	23
Español (ES)	31
Português Portugal (PT-PT)	38
Norsk (NO)	45
Ελληνικά (EL)	52
中文 (ZH-CN)	60
Lietuviškai (LT)	66
Magyar (HU)	73
Polski (PL)	80
Čeština (CZ)	88
한국어 (KO)	95
Türkçe (TR)	102
Русский (RU)	109
Română (RO)	117
Србија (SR)	124
Português Brasil (PT-BR)	131
Nederlands (NL)	138
Tiếng Việt (VI)	145
Македонски (MK)	152
REFERENCES	160

MicroScan

HNID Procedural Manual

REF B1012-10B

INTENDED USE

The HNID system is designed for rapid identification of and detection of β -lactamase production in *Neisseria* species, *Haemophilus* species, *Branhamella catarrhalis* and *Gardnerella vaginalis*. Biotyping of *Haemophilus influenzae* and *Haemophilus parainfluenzae* may also be performed.

SUMMARY AND PRINCIPLE

Modified conventional and chromogenic tests are used for identification of fastidious Gram negative cocci and bacilli.^{1,2,3,4,5,6,7,8,9,10,11,12} A 3 McFarland suspension of the organism to be tested is made in a special inoculum broth. The HNID panel is rehydrated and inoculated from the broth and incubated for 4 hours. After incubation, changes in pH and utilization of certain substrates or other biochemical reactions are interpreted for identification. A rapid acidometric method is used for the detection of the enzyme β -lactamase.^{13,14,15}

WARNING AND PRECAUTIONS

HNID panels and reagents are for *in vitro* diagnostic use only. Do not allow Peptidase Reagent or HNID Indole Reagent to come in contact with eyes, skin or clothing. Observe aseptic techniques and established precautions against microbiological hazards throughout the procedure, with special awareness that the inoculated panels contain potentially pathogenic organisms.

STORAGE

Store between 2-8°C. After rehydration with HNID Inoculum Broth, discard the panel if any substrate appears to have a positive reaction. Refer to the Biochemical Interpretation Chart for positive color reactions.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION AND PREPARATION

Specimens are collected and placed on primary isolation media according to laboratory procedures as described in the *Manual of Clinical Microbiology*.¹⁶

MATERIALS PROVIDED

HNID Panel (B1012-10B)

Procedural Manual

HNID Biotype Codebook

Worksheet

MATERIALS REQUIRED, BUT NOT PROVIDED

HNID Inoculum Broth (B1015-16)

Peptidase Reagent (B1012-30B)

Cover Tray (B1010-56B)

Mineral Oil, 60 mL (B1010-40)

Mineral Oil, 250 mL (B1010-40A) - ONLY for use with WalkAway *SI* instruments (and WalkAway instruments upgraded with the automated oil overlay feature)

0.5% N, N-Dimethyl- α -Naphthylamine (B1010-45A)

HNID Indole Reagent (B1015-15)

0.8% Sulfanilic Acid (B1010-44A)

50 μ L Pipettor with Sterile Tips

Sterile Wooden Applicator Sticks or Sterile Inoculating Loop

3 McFarland Barium Sulfate Turbidity Standard

4 McFarland Barium Sulfate Turbidity Standard

Vortex mixer

35-37°C non-CO₂ Incubator

Quality Control Organisms:

Branhamella catarrhalis ATCC 49143*
Gardnerella vaginalis ATCC 49145
Haemophilus influenzae ATCC 49144
Haemophilus paraphrophilus ATCC 49146
Neisseria lactamica ATCC 49142
Neisseria mucosa ATCC 49233
Stenotrophomonas maltophilia ATCC 49130

PROCEDURE

1. Panel Preparation

Remove the panel from refrigeration and *immediately* remove from foil pouch. **Panels should not be used if the dessicant is not present or is broken or if the integrity of the packaging is compromised (unsealed, punctured or torn).** Allow the panel and the HNID Inoculum Broth to come to room temperature before use. Panels may be stacked with a clean cover tray on top.

2. Organism Preparation

Using a sterile applicator stick or a sterile inoculating loop, transfer sufficient morphologically similar colonies from an 18-24 hour primary isolation plate (Blood, Chocolate, GV or Thayer Martin Agar) to the HNID Inoculum Broth comparable to a 3 McFarland Barium Sulfate Turbidity Standard. For *Neisseria* species, a 4 McFarland Barium Sulfate Turbidity Standard is desirable. Vortex the Inoculum Broth before inoculating the panel. Do not use less than a 3 McFarland Barium Sulfate Turbidity Standard to inoculate any HNID panels.

3. Inoculation

Transfer 50 µL (0.05 mL) of Inoculum Broth to each labeled well of the HNID panel. A 50 µL pipettor with sterile tip should be used. After each panel inoculation, MicroScan recommends subculturing 1 drop of the excess Inoculum Broth onto a suitable medium as a purity check.

NOTE: For use with autoSCAN-4 and WalkAway Systems the LOC well must also be rehydrated.

4. Biochemical Oil Overlays

Using a dropper bottle, overlay the URE well with 2-3 drops of Mineral Oil. The media in the well should be completely covered, but the well should not be overfilled.

NOTE: WalkAway *SI* instruments (and WalkAway instruments upgraded with the automated oil overlay feature) automatically add oil to the URE well.

5. Incubation

Incubate the inoculated panel for 4 hours at 35-37°C in a non-CO₂ incubator. Stack panels in groups of 4 to 5 with a Cover Tray over the top panel to prevent evaporation.

6. Visual or touchSCAN Panel Reading

A. After 4 hours incubation at 35-37°C, those tests which do not require the addition of reagents (IDX, GLU, SUC, MAL, FRU, LAC, GAL, URE, ORN, AGL, STA and BL) should be read first.

B. Add reagents as follows:

a. PRO, NGL, ZAR

Add 1 drop of Peptidase Reagent to the PRO, NGL, and ZAR wells. Allow 3-5 minutes for the color reaction to develop.

b. NO₃, NO₂

Add 1 drop of 0.8% Sulfanilic Acid and 1 drop of 0.5% N, N-Dimethyl-α-Naphthylamine to the NO₃ and NO₂ wells. Allow 3-5 minutes for the color reaction to develop.

c. IND

Add 1-2 drops of HNID Indole Reagent to the IND well. Color will develop immediately.

C. Interpret the reactions according to the Biochemical Interpretation Chart and record results on the worksheet provided, or input the positive reactions via the touchSCAN reader.

Table 1, Substrates

Substrates	Abbr.	Substrates	Abbr.
Indoxyl Phosphate	IDX	Urea	URE
Nitrate	NO ₃	Ornithine	ORN

*American Type Culture Collection, Manassas, VA USA

Substrates	Abbr.	Substrates	Abbr.
Nitrite	NO ₂	Indole	IND
Glucose	GLU	L-Prolyl-β-Naphthylamide	PRO
Sucrose	SUC	N-γ-L-Glutamyl-β-Naphthylamide	NGL
Maltose	MAL	N-α-Benzoyl-DL-Arginine-β-Naphthylamide	ZAR
Fructose	FRU	p-Nitrophenyl-α-D-glucoside	AGL
Lactose	LAC	Starch	STA
o-Nitrophenyl-β-D-Galactoside	GAL	β-Lactamase (Penicillin)	BL

7. autoSCAN-4 Panel Reading

- After 4 hours incubation at 35-37°C, overlay the GLU, SUC, MAL, FRU, LAC, STA and BL with mineral oil before adding reagents.
- Add reagents as follows:
 - PRO, NGL, ZAR
Add 1 drop of Peptidase Reagent to the PRO, NGL, and ZAR wells. Allow 3-5 minutes for the color reaction to develop.
 - NO₃, NO₂
Add 1 drop of 0.8% Sulfanilic Acid and 1 drop of 0.5% N, N-Dimethyl-α-Naphthylamine to the NO₃ and NO₂ wells. Allow 3-5 minutes for the color reaction to develop.
 - IND
Add 1-2 drops of HNID Indole Reagent to the IND well. Color will develop immediately.
- After each reaction has developed for the appropriate amount of time, insert the panel into the autoSCAN-4 for reading.

8. WalkAway Systems Panel Reading

- After 4 hours incubation, those tests which do not require the addition of reagents (IDX, GLU, SUC, MAL, FRU, LAC, GAL, URE, ORN, AGL, STA, and BL) are read first.
- The WalkAway then adds the appropriate reagents to the following: PRO, NGL, ZAR, NO₃, NO₂, and IND. These reactions are then interpreted.

NOTE: The WalkAway reads wells that do not require reagents first to prevent altered results due to reagent fumes that may be trapped under the tray lid. Therefore, biochemicals cannot be verified manually.

RESULTS

Table 2

Biochemical Interpretation Chart			
Well	Positive	Negative	Notes
IDX	Blue to Gray Precipitate	Clear	
NO ₃	Pink, Red to Magenta	Clear to Yellow	Results recorded after 3-5 minutes.
NO ₂	Clear to Yellow	Pink, Red to Magenta	Results recorded after 3-5 minutes.
GLU SUC LAC	Orange to Yellow	Red to Red/Orange	Any red color should be recorded as negative.
FRU	Yellow to Yellow/Orange	Red to Orange	Any yellow color should be recorded as positive.
MAL	Yellow	Red to Yellow/Orange	
GAL	Yellow	Clear	A pale yellow is positive. Use IND well before addition of reagent as a negative control.
URE	Red to Magenta	Yellow to Light Pink	
ORN	Green to Blue	Yellow to Yellow/Green	Any yellow color should be recorded as negative.
IND	Red Ring	Yellow	Any red color should be recorded as positive.
PRO NGL ZAR	Red to Magenta	Yellow to Orange	Results should be recorded within 3-5 minutes. A weak ZAR reaction is indicated by a cloudy appearance.
AGL	Yellow	Clear	Any yellow color should be recorded as positive. Use the IND well before the addition of reagent as a negative control.
STA BL	Yellow to Yellow/Orange	Red to Orange	Any yellow color should be considered as positive.
LOC	For autoSCAN-4 and WalkAway System Panel Alignment		

Reactions

Indoxyl Phosphatase (IDX) - Alkaline phosphatase splits indoxyl phosphate to free indoxyl which combines with molecular oxygen to form an insoluble blue compound.

Nitrate (NO₃) - The ability of an organism to reduce nitrate to nitrite is detected by the addition of Sulfanilic Acid and N, N-Dimethyl- α -Naphthylamine which produce a red color in the presence of nitrite.

Nitrite (NO₂) - The ability of an organism to reduce nitrite to free nitrogen gas is detected by the addition of Sulfanilic Acid and N, N-Dimethyl- α -Naphthylamine. In the presence of nitrite, a red color develops which is a negative test. If no color develops, the nitrite has been reduced and the test is a positive test.

Carbohydrate Utilization (GLU, SUC, MAL, FRU, LAC, STA) - The utilization of a specific carbohydrate results in acid formation with the consequent pH drop indicated by a phenol red indicator.

β -Galactosidase (GAL) - β -galactosidase hydrolyzes o-nitrophenyl- β -D-galactopyranoside which liberates the yellow-colored o-nitrophenol.

α -Glucosidase (AGL) - α -glucosidase hydrolyzes p-nitrophenyl- α -D-glucoside which liberates the yellow-colored p-nitrophenol.

Urea (URE) - The enzyme urease splits urea forming ammonia with the resulting pH increase detected by a phenol red indicator.

Ornithine (ORN) - Decarboxylation of ornithine results in an alkalinization of the media which is detected by the pH indicator bromthymol blue.

Indole (IND) - The metabolism of tryptophan results in the formation of indole which is detected by the addition of HNID Indole Reagent. If indole is present, a red ring forms.

Peptidases (PRO, NGL, ZAR) - The ability of an organism to produce a specific peptidase enzyme is determined by splitting the appropriate amino acid from its linkage to β -naphthylamide. The free β -naphthylamide is then detected by complexing with p-dimethylaminocinnamaldehyde producing a red to magenta color.

β -Lactamase (BL) - Penicillin G is degraded by β -lactamase to Penicilloic Acid. The pH shift is detected by a phenol red indicator.

Organism Identification

The MicroScan HNID Biotype Codebook is used for the identification of unknown test organisms. This Codebook has been generated from MicroScan's data base of information compiled from tests on the HNID panel.

The Biotype Codebook provides for computer analysis of 17 tests to identify unknown organisms. These tests are translated into a 6-digit "biotype number," for which the Codebook lists the organism identification and its relative probability. Refer to the HNID Biotype Codebook for the procedure for generating a biotype number. When only one organism is listed for a biotype number, the relative probability is 99.9%. If more than one species identification is possible for a particular biotype number, all possibilities are printed in the order of the highest probability up to a cumulative total probability of 99.9%.

Should a biotype number occur that is not in the Codebook, we suggest repeating the test from a fresh culture. If retesting yields the same results, consult the Biotype Lookup Service on the Beckman Coulter website or consult your local distributor.

LIMITATION OF THE PROCEDURE

1. Should more than one species identification be listed for a particular biotype number, supplemental tests may be necessary for a correct identification. Testing for other criteria, such as X and V factors and aerotolerance should also be performed.
2. Only *Haemophilus* species, *Neisseria* species, *Branhamella catarrhalis* and *Gardnerella vaginalis* can be tested - do not use this product for other coccobacillary organisms.
3. The interpretation of test results requires trained clinical personnel who should use judgement, knowledge and additional confirmatory testing where required prior to accepting the identification of an organism. Please refer to the footnotes listed in the HNID Biotype Codebook or the Microbiology Information Manual.
4. Biotype numbers should not be used for phenotypic identification of strains isolated from various specimens from the same patient.
5. CLSI guidelines M7-A7, Vol. 26, No. 2, recommends nitrocefin-based β -Lactamase test for testing *Branhamella catarrhalis*.¹⁷
6. Biotype groups for *Haemophilus* sp. do not correlate with serotypes, i.e., *Haemophilus influenzae* II is not definitively group B. Serotyping is necessary.
7. Occasionally, some strains of *N. meningitidis* do not produce acid from maltose. These strains may give an incorrect identification of *N. gonorrhoeae* if NGL is also negative. Please refer to the footnotes listed in the HNID Biotype Codebook or the Microbiology Information Manual.
8. *H. haemolyticus* isolates grown on chocolate agar may be incorrectly identified as *H. influenzae*. Hemolysis on horse blood agar may be used to distinguish *H. haemolyticus* from *H. influenzae*.

EXPECTED VALUES

Refer to the HNID Percent Positive Reaction chart for the expected results for each test and each species. (See Table 3.)

	IDX	NO ₃	NO ₂	GLU	SUC	MAL	FRU	LAC	GAL	URE	ORN	IND	PRO	NGL	ZAR	AGL	STA
<i>Haemophilus influenzae I</i>	99	99	10	95	1	5	90	1	1	99	99	99	50	1	1	1	1
<i>H. influenzae II</i>	99	99	5	95	1	10	90	1	1	99	1	99	50	1	1	1	1
<i>H. influenzae III</i>	99	99	5	95	1	25	75	1	1	99	1	1	25	1	1	1	1
<i>H. influenzae IV</i>	99	99	1	95	1	10	95	1	1	99	99	1	10	1	1	1	1
<i>H. influenzae V</i>	99	99	10	95	1	5	90	1	1	1	99	99	10	1	1	1	1
<i>H. influenzae VI-1²</i>	99	99	1	95	1	1	75	1	1	1	99	1	1	1	1	1	1
<i>H. influenzae VI-2²</i>	99	1	1	95	1	1	75	1	1	1	1	99	25	1	1	1	1
<i>H. parainfluenzae I</i>	99	99	1	95	95	95	95	1	75	1	99	1	1	1	1	1	1
<i>H. parainfluenzae II</i>	99	99	50	95	95	95	95	1	50	99	99	1	1	1	1	1	1
<i>H. parainfluenzae III</i>	99	99	25	95	95	95	95	10	50	99	1	1	25	1	1	1	1
<i>H. parainfluenzae IV</i>	99	99	1	95	95	95	95	1	25	1	1	1	1	1	1	1	1
<i>H. aphrophilus/paraphrophilus</i>	99	99	25	95	95	95	95	95	95	1	1	1	10	99	1	75	1
<i>H. haemolyticus</i>	99	99	1	95	1	95	95	1	1	95	50	95	1	1	1	1	1
<i>Neisseria gonorrhoeae</i>	1	1	25	95	1	1	1	1	1	1	1	1	99	1	1	1	1
<i>N. meningitidis</i>	1	1	50	95	1	95	1	1	1	1	1	1	99	99	1	1	1
<i>N. lactamica</i>	1	1	50	95	1	95	1	95	99	1	1	1	99	1	1	1	1
<i>N. mucosa (NIT +)</i>	1	99	99	75	95	95	75	1	1	1	1	1	99	1	1	1	1
<i>Neisseria species</i>	1	1	95	90	90	95	90	10	1	1	1	1	99	10	10	1	1
<i>Branhamella catarrhalis</i>	1	75	95	1	1	1	1	1	1	1	1	1	99	1	1	1	1
<i>Gardnerella vaginalis</i>	1	1	1	95	75	95	95	50	50	1	1	1	99	1	1	99	99

1. Positive Reactions Interpretations Values shown in this table are not actual percentages but representation based on published data and MicroScan results.
2. *Haemophilus influenzae VI-1* is the Biotype of Oberhofer biotype classification. *Haemophilus influenzae VI-2* is the Biotype VI of the Kilian biotype classification.¹⁸

Expected reactions have been expressed by nine values:					
Value	Meaning	Value	Meaning	Value	Meaning
1	Pos. Reaction Very Rare	25	Majority Reaction Negative	90	Neg. Reaction Uncommon
5	Pos. Reaction Rare	50	Pos. or Neg. Reaction	95	Neg. Reaction Rare
10	Pos. Reaction Uncommon	75	Majority Reaction Positive	99	Neg. Reaction Very Rare

QUALITY CONTROL

Culture purity should be checked by subculturing from the inoculum suspension onto suitable agar media. Test results are invalid if a mixed culture is detected. The acceptability of the identification media may be checked by testing with the following strains: (Table 4.)

Table 4

Quality Control Chart

	<i>N. lactamica</i> ATCC 49142	<i>H. influenzae</i> ATCC 49144	<i>H. paraphrophilus</i> ATCC 49146	<i>B. catarrhalis</i> ATCC 49143
IDX	-	+	+	-
NO ₃	-	+	+	V
NO ₂ ¹	V	V	-	V
GLU	+	+	+	-
SUC	-	-	+	-
MAL	+	-	+	-
FRU	-	V	+	-
LAC	+	-	+	-
GAL	+	-	+	-
URE	-	+	-	-
ORN	-	+	-	-
IND	-	+	-	-
PRO	+	V	-	V
NGL	-	-	+	-
ZAR ²	-	-	-	-
AGL ³	-	-	V	-
STA ³	-	-	-	-
BL	-	+	-	-

1. A positive result may be obtained with *Neisseria mucosa* ATCC 49233.
2. A positive result may be obtained with *Stenotrophomonas maltophilia* ATCC 49130.
3. A positive result may be obtained with *Gardnerella vaginalis* ATCC 49145.


WARRANTY STATEMENT

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

Symbols Key

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1; 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1; 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1; 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1; 5.3.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1; 5.5.5
	<i>In Vitro</i> Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1; 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1; 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a

Symbols Key

Symbol	Symbol Title	Description of Standard	Standard
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1; 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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REFERENCES

1. Baron, E.S., and A. K. Saz. 1976. Effect of types of media on the production of acid from glucose by so-called glucose-negative strains of *N. gonorrhoeae*. *J. Clin. Microbiol.* 3:330.
2. Brooker, D. C., M. S. Lund, and D. J. Blazevic. 1973. Rapid test for lysine decarboxylase activity in Enterobacteriaceae. *Appl. Microbiol.* 26:622.
3. D'Amato, R.F., L.A. Eriquez, K.M. Tomfohrde, and E. Singerman. 1978. Rapid Identification of *Neisseria gonorrhoeae* and *Neisseria meningitidis* by using enzymatic profiles. *J. Clin. Microbiol.* 7:77.
4. Davis, J.A., J.R. Mitzel, and W.E. Beam, Jr. 1971. Carbohydrate fermentation patterns of *Neisseria meningitidis* determined by a microtiter method. *Appl. Microbiol.* 21:1072.
5. Djurhuus, R. 1981. Ornithine decarboxylase (E.C.4.1.1.17) assay based upon the retention of putrescine by a strong cation-exchange paper. *Anal. Biochem.* 113:352.
6. Enriquez, L.A., and G.B. Knight. 1980. Aminopeptidases highly specific to glutamyl residues from *Neisseria meningitidis* and *Moraxella urethralis*. *J. Clin. Microbiol.* 12:667.
7. Fay, G.D. and A.L. Barry. 1972. Rapid ornithine decarboxylase test for the identification of Enterobacteriaceae. *Appl. Microbiol.* 23:710.
8. Kellogg, Jr., D.S., and E.M. Turner. 1973. Rapid fermentation confirmation of *Neisseria gonorrhoeae*. *Appl. Microbiol.* 25:550.
9. Lampe, A.S. 1981. Nonliquid reagent for detecting nitrate reduction. *J. Clin. Microbiol.* 14:452.
10. Plot, P., E. Van Dyck, P.A. Totten, and K.K. Holmes. 1982. Identification of *Gardnerella vaginalis*. *J. Clin. Microbiol.* 15:19.
11. Qadri, S.M. Hussain, and N. Parker. 1981. Clinical evaluation of the rapid carbohydrate degradation microtube method for identification of *Neisseria* species. *J. Clin. Microbiol.* 14:318.
12. Tyagi, A.K., C.W. Tabor, and H. Tabor. 1981. Ornithine decarboxylase from *Saccharomyces cerevisiae*. *J. Biol. Chem.* 256:12156.
13. Jorgensen, J.H., J.C. Lee, and G.A. Alexander. 1977. Rapid penicillinase paper strip test for the detection of beta-lactamase producing *Haemophilus influenzae* and *Neisseria gonorrhoeae*. *Antimicrob. Agent. Chemother.* 11:1087.
14. Park, C.H., J.S. Lopez, and C.B. Cook. 1978. Acidometric agar plate method for ampicillin susceptibility testing of *Haemophilus influenzae*. *Antimicrob. Agent. Chemother.* 13:318.
15. Weissfeld, A.S., G.D. Sanner, J.R. Childress, J.D. Dyckman, T.W. Huber, and R.P. Williams. 1977. Feasibility of screening for penicillinase producing *Neisseria gonorrhoeae* from primary culture plates by using a rapid microacidometric test. *Antimicrob. Agent. Chemother.* 12:703.
16. Balows, A., W. J. Hausler, K. L. Herrmann, H. D. Isenberg, and H.J. Shadomy (eds.) 1991. *Manual of Clinical Microbiology*, 5th ed. American Society for Microbiology, Washington, D.C.
17. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically. 2006. Approved Standard M7-A7. Clinical and Laboratory Standards Institute, Wayne, PA.
18. Oberhofer, T.R., and A. E. Bach. 1979. Biotypes of *Haemophilus* encountered in clinical laboratories. *J. Clin. Microbiol.* 10:168.

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.