MicroScan

MICroSTREP *plus* 6 Procedural Manual, Multi-regional

REF B1016-170

TABLE OF CONTENTS

English	2
Français (FR)	12
Deutsch (DE)	22
Italiano (IT)	33
Español (ES)	43
Português Portugal (PT-PT)	53
Norsk (NO)	63
Ελληνικά (EL)	73
中文 (ZH-CN)	84
Lietuviškai (LT)	93
Magyar (HU)	103
Polski (PL)	113
Čeština (CZ)	124
Slovák (SK)	134
한국어 (KO)	144
Türkçe (TR)	153
Русский (RU)	163
Hrvatski (HR)	174
Български (BG)	184
中文 (ZH-TW)	195
Română (RO)	204
Србија (SR)	214
Português Brasil (PT-BR)	224
Nederlands (NL)	234
Tiếng Việt (VI)	244
Македонски (MK)	254
REFERENCES	264

MicroScan

MICroSTREP *plus* 6 Procedural Manual, Multi-regional

REF B1016-170

INTENDED USE

MicroScan MICroSTREP *plus* 6 Panel is an *in vitro* diagnostic medical device intended for quantitative and qualitative determination of antimicrobial agent susceptibility of streptococci including *Streptococcus pneumoniae* from isolated colonies. The methods for incubation and results determination include non-automated offline incubation and manual reading, as well as automated incubation and reading the WalkAway instrumentation and LabPro data management system.

INTENDED USER

This product is intended for laboratory professional use.

SUMMARY AND PRINCIPLE

The antimicrobial susceptibility tests are miniaturizations of the broth dilution susceptibility test that have been dehydrated. Various antimicrobial agents are diluted in water or buffer to concentrations bridging the range of clinical interest. Panels are rehydrated with Mueller-Hinton broth supplemented with 2-5% lysed horse blood (LHB), after inoculation of the broth with a standardized suspension of the organism.

After incubation in a non-CO₂ incubator for 20-24 hours, the minimum inhibitory concentration (MIC) for the test organism is determined by observing the lowest antimicrobial concentration showing inhibition of growth.^{1,2,3,4,5,6,7,8,9,10,11,12,13} Alternatively, the panel can be incubated and read by the MicroScan WalkAway System.^{1,2,4,5,6,7,8,10,11,12}

Antimicrobial Agents	Abbr.	MICroSTREP <i>plus</i> Type 6 B1016-170 Dilutions µg/mL
Amoxicillin-K Clavulanate	Aug	1/0.5-4/2
Ampicillin	Am	0.06-4
Azithromycin	Azi	0.12-2
Cefepime	Cpe	0.25-2
Cefotaxime	Cft	0.25-2
Ceftriaxone	Cax	0.25-2
Cefuroxime	Crm	0.25-2
Chloramphenicol	C	2-8
Clarithromycin	Cla	0.12-2
Clindamycin	Cd	0.06-0.5
Daptomycin	Dap	0.25-2
Erythromycin	E	0.06-0.5
Levofloxacin	Lvx	0.5-4
Linezolid	Lzd	1-4
Meropenem	Mer	0.25-2
Minocycline	Min	0.5-2
Moxifloxacin	Mxf	0.25-2
Penicillin	P	0.03-4
Pristinamycin	Prs	1-2
Rifampin	Rif	0.5-2
Tetracycline	Te	1-4
Trimethoprim-Sulfamethoxazole	T/S	0.25/4.75-2/38
Vancomycin	Va	0.25-4

CLINICAL RELEVANCE

Antimicrobial susceptibility test results are used in conjunction with other laboratory results and/or clinical indicators to refine the selection of antimicrobial therapy in symptomatic patients.

WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with special awareness that inoculated panels contain potentially pathogenic organisms.
3. This material contains infectious agents and should be discarded as appropriate for biohazardous waste.
4. Broth used to rehydrate panel contains 2-5% Lysed Horse Blood.

5. Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.
6. Single Use only. All MicroScan panels and consumables should not be reused unless otherwise noted (e.g. cover trays).
7. This product contains material(s) of animal origin from dehydrated microbiology media. Observe general safety guidelines for protection when handling this product.

STORAGE

MICroSTREP *plus* panels must be stored at 2-25°C.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckmancoulter.com/techdocs

EVIDENCE OF DETERIORATION

Prolonged exposure to storage conditions other than those recommended may result in loss of potency of the antimicrobial agents. Do not use beyond the expiration date. Contact your Beckman Coulter Representative or Distributor for further assistance.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens should be collected, transported, and placed on primary isolation media according to procedures recommended in the *Manual of Clinical Microbiology*.³

MATERIALS PROVIDED

MICroSTREP *plus* Type 6 (B1016-170)

Procedural and QC Manual

MATERIALS REQUIRED, BUT NOT PROVIDED

0.5 McFarland Turbidity Standard

100 µL pipettor with disposable sterile tips

Cover Trays (B1010-56B or equivalent)

Barcode Label Paper (B1018-129)

General Laboratory Equipment

Inoculator-D Set, 60/pk (B1013-5) or 240/pk (B1013-4)

Inoculum Water, 3 mL (B1015-2)

Inoculum Saline, 3 mL (B1015-12)

Mueller-Hinton Broth with 3% lysed horse blood, 25 mL (B1015-25)

Microdilution Viewer

Quality Control Organisms (Refer to QC section of this manual)

RENOK Rehydrator/Inoculator System (B1018-14) or equivalent

Turbidity Meter

Vortex

WalkAway Tray Lids (B1018-18)

PROCEDURE

Broth Preparation

1. Allow the Mueller-Hinton Broth with 3% Lysed Horse Blood to come to room temperature before use.
2. Prior to inoculating panels, ensure that the broth has not been compromised by improper storage or handling by examining it for the following properties:
 - A. Broth should be clear and red to dark red in color.

- B. The broth may become slightly brown in color as it gets closer to expiration date but should remain clear. The brownish color will not affect panel results.
- C. Do not use any broth that is turbid and discolored.

Panel Preparation

1. Remove the panels to be used from storage. Do not use if the integrity of the packaging is compromised (unsealed, punctured or torn).
2. Cut open the pouch and remove the panel. All wells on the panel, except the locator well, should appear clear prior to inoculation. If stored in the refrigerator, remove the panel immediately from the foil pouch.
3. Panels should not be used if any of the following conditions exist:
 - A. Desiccant is not present or is broken.
 - B. Panel wells are discolored.
4. Allow panels to equilibrate to room temperature prior to rehydration. Panels may be stacked with a clean Cover Tray on top. All opened panels should be used within the same day or discarded.

Inoculum Preparation

CLSI recommends periodically checking inoculum densities by performing colony counts, refer to CLSI document M07-A10.¹⁰

Turbidity Standard Technique

THE TURBIDITY STANDARD TECHNIQUE IS THE ONLY METHOD APPROVED FOR USE WITH THIS PANEL.

1. Using a sterile applicator stick, loop or swab, touch the surface of at least 4-5 large or 5-10 small morphologically similar, well-isolated colonies or enough to create a 0.5 McFarland inoculum from a 16-20 hour noninhibitory agar plate (e.g., Trypticase Soy Agar with 5% Sheep Blood, Chocolate Agar). Do not prepare inoculum from plates with pin point colonies that exhibit autolysis. Incubation in 5-10% CO₂ is recommended.
2. Emulsify in 3 mL of Inoculum Water (autoclaved deionized water) or Inoculum Saline (autoclaved 0.9% NaCl).
3. Cap tightly and vortex the suspension for 2-3 seconds. The final turbidity should be equivalent to that of a 0.5 McFarland Turbidity Standard. An equivalent turbidity can be achieved by using a MicroScan Turbidity Meter with a range of 0.08 ± 0.02.
4. Pipet 0.1 mL (100 µL) of the standardized suspension into 25 mL of LHB broth. Cap tightly and invert 8-10 times to mix.

Extending the time between inoculation preparation and panel inoculation to greater than 15 minutes may adversely affect the viability of the streptococcus.

Panel Rehydration/Inoculation

Inoculation of the panels is performed using the RENOK Rehydrator/Inoculator System with Inoculator-D sets. If an alternate system is used, rehydrate each well with 115 µL ± 10 µL of the appropriate suspension.

To ensure viability and purity of the organism tested, a purity plate should be prepared by streaking the inoculum to an appropriate agar plate and incubating under the appropriate conditions. If two or more colony types are present on the purity plate, reisolate the colonies and retest.

1. Remove the transfer lid from the inoculator set. Pour the inoculated 25 mL broth (LHB) into the seed trough.
2. Replace the transfer lid. Gently tap the transfer lid in all four corners to ensure absence of bubbles. Allow the transfer lid to equilibrate for a minimum of 20 seconds in the seed trough.
3. Attach the RENOK Rehydrator/Inoculator to the transfer lid and draw up the inoculum. Refer to the RENOK Rehydrator/Inoculator System Operator's Manual for operating instructions.

NOTE: Due to the physical properties of the lysed horse blood supplement, a priming step is recommended to minimize the possibility of uneven filling of the panel wells. To perform the priming step, draw up the inoculum as per above. While the transfer lid is still seated on the seed trough, press the center release button to dispense the inoculum back into the seed trough. Draw up the inoculum a second time into the transfer lid.

4. Place the RENOK Rehydrator/Inoculator transfer lid apparatus on the MICroSTREP *plus* Panel and release the inoculum into the panel.

5. Visually check the volume of inoculated wells to confirm proper filling.

Incubation

1. Panels can be incubated in a WalkAway System or incubated off-line using the following steps:
 - A. To ensure even thermal distribution during incubation, stack the panels in groups of 3-5.
 - B. Place a clean Cover Tray on top of each group of panels to prevent evaporation. Cover Trays may be reused. Do not decontaminate the Cover Trays with alcohol. They may be cleaned with soap and water. Rinse well and allow to air dry.
 - C. Incubate the panels for 20-24 hours at $35^{\circ}\text{C} \pm 1^{\circ}\text{C}$ in a non- CO_2 incubator.

Reading the Panels

Panels are read manually using indirect light (such as that provided by the MicroScan Microdilution Viewer) and results can be recorded on a panel worksheet, or the panels can be read on the MicroScan WalkAway System. Refer to the appropriate operator's manual for reading panels with the MicroScan WalkAway System.

1. Following 20-24 hours incubation, remove the panels from the incubator.
2. Wipe the bottom of the panel with a lint-free tissue to remove any condensation or debris that may be present.
3. Read the panels only if the growth control well is turbid. Inadequate growth is defined as a slight haze in the well or if the broth is clear.
4. Record results on the appropriate worksheet.
5. Reading Antimicrobial Susceptibilities (MICs)
 - A. Growth in the wells appears as turbidity. Turbidity may take the form of a haze throughout the well, a button in the center of the well, or a fine granular growth throughout the well. A change in the broth color should not be interpreted as growth. Turbidity is read using indirect light.
 - B. Record MIC results as follows:
 - a. Following 20-24 hours incubation, record the MIC as the lowest antimicrobial concentration showing inhibition of growth.
 - b. When growth occurs in all concentrations of an antimicrobial, the MIC is recorded as greater than (>) the highest concentration.
 - c. When no growth occurs in any of the concentrations of antimicrobics, the MIC is recorded as less than or equal to ($\leq$) the lowest concentration.
 - d. A clear well in a series of growth wells (e.g., growth at 1, 2, and 8 $\mu\text{g/mL}$, but not at 4 $\mu\text{g/mL}$) is called a skipped well and should be ignored.
 - e. Spot contamination in isolated wells indicates contamination. The test should be repeated.

RESULTS

Interpretation of MIC Results

Susceptibility is determined by comparing the MIC of an organism to the attainable blood or urine level of the antimicrobial. The following table lists the streptococci interpretive criteria as indicated.

Streptococci Interpretive Criteria¹

Antimicrobial Agents	Abbr.	EUCAST ²			CLSI ³			OTHER ⁴		
		S	I	R	S	I	R	S	I	R
Amoxicillin-K Clavulanate - <i>S. pneumoniae</i>	Aug	-	-	-	$\leq 2/1$	4/2	$\geq 8/4$	-	-	-
Ampicillin	Am									
<i>S. pneumoniae</i>		≤ 0.5	1-2	>2	-	-	-	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤ 0.25	-	-	-	-	-
Viridans Streptococci		≤ 0.5	1-2	>2	≤ 0.25	0.5-4	≥ 8	-	-	-
Azithromycin ^{6,7}	Azi									
<i>S. pneumoniae</i>		≤ 0.25	0.5	>0.5	≤ 0.5	1	≥ 2	-	-	-
Streptococci (Group A, B, C and G)		≤ 0.25	0.5	>0.5	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤ 0.5	1	≥ 2	-	-	-
Viridans Streptococci		-	-	-	≤ 0.5	1	≥ 2	-	-	-
Cefepime	Cpe									
<i>S. pneumoniae</i>		≤ 1	2	>2	≤ 1	2	≥ 4	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤ 0.5	-	-	-	-	-
Viridans Streptococci		≤ 0.5	-	>0.5	≤ 1	2	≥ 4	-	-	-

Antimicrobial Agents	Abbr.	EUCAST ²			CLSI ³			OTHER ⁴		
		S	I	R	S	I	R	S	I	R
Cefotaxime	Cft									
<i>S. pneumoniae</i>		≤0.5	1-2	>2	-	-	-	-	-	-
<i>S. pneumoniae</i> (meningitis)		-	-	-	≤0.5	1	≥2	-	-	-
<i>S. pneumoniae</i> (non-meningitis)		-	-	-	≤1	2	≥4	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤0.5	-	-	-	-	-
Viridans Streptococci		≤0.5	-	>0.5	≤1	2	≥4	-	-	-
Ceftriaxone	Cax									
<i>S. pneumoniae</i>		≤0.5	1-2	>2	-	-	-	-	-	-
<i>S. pneumoniae</i> (meningitis) ^{8,9}		-	-	-	≤0.5	1	≥2	-	-	-
<i>S. pneumoniae</i> (non-meningitis) ^{8,9}		-	-	-	≤1	2	≥4	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤0.5	-	-	-	-	-
Viridans Streptococci		≤0.5	-	>0.5	≤1	2	≥4	-	-	-
Cefuroxime axetil (oral) - <i>S. pneumoniae</i>	Crn	-	-	-	≤1	2	≥4	-	-	-
Cefuroxime sodium (parenteral)	Crn									
<i>S. pneumoniae</i>		≤0.5	1	>1	≤0.5	1	≥2	-	-	-
Viridans Streptococci		≤0.5	-	>0.5	-	-	-	-	-	-
Chloramphenicol ⁶	C									
<i>S. pneumoniae</i>		≤8	-	>8	≤4	-	≥8	-	-	-
Streptococci (Group A, B, C and G)		≤8	-	>8	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤4	8	≥16	-	-	-
Viridans Streptococci		-	-	-	≤4	8	≥16	-	-	-
Clarithromycin ^{6,7}	Cla									
<i>S. pneumoniae</i>		≤0.25	0.5	>0.5	≤0.25	0.5	≥1	-	-	-
Streptococci (Group A, B, C and G)		≤0.25	0.5	>0.5	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤0.25	0.5	≥1	-	-	-
Viridans Streptococci		-	-	-	≤0.25	0.5	≥1	-	-	-
Clindamycin ⁶	Cd									
<i>S. pneumoniae</i>		≤0.5	-	>0.5	≤0.25	0.5	≥1	-	-	-
Streptococci (Group A, B, C and G)		≤0.5	-	>0.5	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤0.25	0.5	≥1	-	-	-
Viridans Streptococci		≤0.5	-	>0.5	≤0.25	0.5	≥1	-	-	-
Daptomycin	Dap									
Streptococci (Group A, B, C and G)		≤1	-	>1	-	-	-	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤1	-	-	-	-	-
Viridans Streptococci ⁵		-	-	-	≤1	-	-	-	-	-
Erythromycin ^{6,7}	E									
<i>S. pneumoniae</i>		≤0.25	0.5	>0.5	≤0.25	0.5	≥1	-	-	-
Streptococci (Group A, B, C and G)		≤0.25	0.5	>0.5	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤0.25	0.5	≥1	-	-	-
Viridans Streptococci		-	-	-	≤0.25	0.5	≥1	-	-	-
Levofloxacin	Lvx									
<i>S. pneumoniae</i> (EUCAST v9.0)		≤2	-	>2	≤2	4	≥8	-	-	-
Streptococci (Group A, B, C and G) (EUCAST v9.0)		≤1	2	>2	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤2	4	≥8	-	-	-
Viridans Streptococci		-	-	-	≤2	4	≥8	-	-	-
Linezolid	Lzd									
<i>S. pneumoniae</i> ⁵ (EUCAST v9.0)		≤2	4	>4	≤2	-	-	-	-	-
Streptococci (Group A, B, C and G) (EUCAST v9.0)		≤2	4	>4	-	-	-	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤2	-	-	-	-	-
Viridans Streptococci ⁵		-	-	-	≤2	-	-	-	-	-
Meropenem ⁵	Mer									
<i>S. pneumoniae</i>		-	-	-	≤0.25	0.5	≥1	-	-	-
<i>S. pneumoniae</i> (meningitis) (EUCAST v8.1)		≤0.25	0.5-1	>1	-	-	-	-	-	-
<i>S. pneumoniae</i> (non-meningitis)		≤2	-	>2	-	-	-	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤0.5	-	-	-	-	-
Viridans Streptococci ⁵		≤2	-	>2	≤0.5	-	-	-	-	-
Minocycline	Min									
<i>S. pneumoniae</i> (EUCAST v9.0)		≤0.5	1	>1	-	-	-	-	-	-
Streptococci (Group A, B, C and G) (EUCAST v9.0)		≤0.5	1	>1	-	-	-	-	-	-
Moxifloxacin	Mxf									
<i>S. pneumoniae</i>		≤0.5	-	>0.5	≤1	2	≥4	-	-	-
Streptococci (Group A, B, C and G) (EUCAST v6.0)		≤0.5	1	>1	-	-	-	-	-	-

Antimicrobial Agents	Abbr.	EUCAST ²			CLSI ³			OTHER ⁴		
		S	I	R	S	I	R	S	I	R
Penicillin	P									
<i>S. pneumoniae</i> (meningitis)		≤0.06	-	>0.06	≤0.06	-	≥0.12	-	-	-
<i>S. pneumoniae</i> (non-meningitis)		≤0.06	0.12-2	>2	≤2	4	≥8	-	-	-
<i>S. pneumoniae</i> (oral)		-	-	-	≤0.06	0.12-1	≥2	-	-	-
Streptococci (Group A, B, C and G)		≤0.25	-	>0.25	-	-	-	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤0.12	-	-	-	-	-
Viridans Streptococci		≤0.25	0.5-2	>2	≤0.12	0.25-2	≥4	-	-	-
Pristinamycin ¹⁰	Prs									
<i>S. pneumoniae</i>		-	-	-	-	-	-	≤1	-	>1
Streptococci		-	-	-	-	-	-	≤1	2	>2
Rifampin - <i>S. pneumoniae</i>	Rif	-	-	-	≤1	2	≥4	-	-	-
Tetracycline	Te									
<i>S. pneumoniae</i>		≤1	2	>2	≤1	2	≥4	-	-	-
Streptococci (Group A, B, C and G)		≤1	2	>2	-	-	-	-	-	-
Beta-Hemolytic Streptococci		-	-	-	≤2	4	≥8	-	-	-
Viridans Streptococci		-	-	-	≤2	4	≥8	-	-	-
Trimethoprim-Sulfamethoxazole	T/S									
<i>S. pneumoniae</i>		≤1/19	2/38	>2/38	≤0.5/9.5	1/19-2/38	≥4/76	-	-	-
Streptococci (Group A, B, C and G)		≤1/19	2/38	>2/38	-	-	-	-	-	-
Vancomycin	Va									
<i>S. pneumoniae</i> ⁵		≤2	-	>2	≤1	-	-	-	-	-
Streptococci (Group A, B, C and G)		≤2	-	>2	-	-	-	-	-	-
Beta-Hemolytic Streptococci ⁵		-	-	-	≤1	-	-	-	-	-
Viridans Streptococci ⁵		≤2	-	>2	≤1	-	-	-	-	-

1. The antimicrobials included in this panel are not proven to be safe and effective in treating clinical infections for all organisms tested. For reporting of antimicrobial results which have shown to be active against organism groups in vitro or in clinical infections refer to appropriate guidance documents or the pharmaceutical package insert.

2. Based on EUCAST v10.¹⁴

3. Based on Interpretive Breakpoints as indicated in CLSI Document M100, 30th ed.¹⁵

4. OTHER therapy may be based on SFM or Manufacturers Breakpoints. Please refer to subsequent footnotes for details.

5. The absence of resistant strains precludes the CLSI from defining any result categories other than "Susceptible" at this time. Strains yielding results suggestive of a "non-susceptible" category should be submitted to a reference laboratory for further testing.

6. Only systemic therapy will be reported.

7. Susceptibility and resistance to azithromycin and clarithromycin can be predicted by testing erythromycin.

8. For *S. pneumoniae* isolated from cerebral spinal fluid (CSF), report only meningitis interpretations. For *S. pneumoniae* isolated from all specimens other than CSF, report interpretations for both meningitis and non-meningitis.

9. Only results of testing with penicillin, vancomycin, and cefotaxime or ceftriaxone or meropenem should be reported routinely for CSF isolates of *S. pneumoniae*.

10. Based on Interpretive Breakpoints as indicated in 2019 Report of the Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM)¹⁶

NOTE: Interpretative criteria listed in the labeling may differ from the criteria in the LabPro Information Manager due to differences in software and panel releases.

LIMITATION OF THE PROCEDURE

1. Some fastidious organisms may require additional supplements and may not grow on the MICroSTREP *plus* panel. If there is no growth in the Growth Well, the results of the panel are not valid.
2. The MICroSTREP *plus* panels are only intended for testing streptococci (including *S. pneumoniae*).
3. A slight haziness may occur in random wells of some antimicrobics due to the incomplete solubility of some media constituents. This should not be interpreted as growth.
4. The interpretation of MIC results require trained clinical personnel who should use judgement and knowledge prior to accepting the susceptibility of an organism.
5. The MICroSTREP *plus* panel should not be incubated in a CO₂ incubator.
6. Results obtained with mixed cultures are not valid.
7. Use an alternative testing method prior to reporting results for *S. pneumoniae* with linezolid when the MIC is 2 µg/mL, if critical to patient care.
8. The ability of the MICroSTREP *plus* panel to detect resistance to moxifloxacin is unknown due the lack of sufficient resistant strains at the time of comparative testing.
9. False susceptibility and false skipped wells have been observed for *S. agalactiae* for MicroScan MICroSTREP *plus* panels read on the WalkAway System. For *S. agalactiae* testing, panels must be read manually.

QUALITY CONTROL

The acceptability of the antimicrobial agents should be checked by testing organisms with known MIC ranges. The results for the recommended American Type Culture Collection (ATCC) control organisms are listed below.*

MICroSTREP *plus* Type 6 (B1016-170)

Antimicrobial Agents	Abbr.	<i>S. pneumoniae</i> ATCC 49619 ¹ Range ²	<i>E. faecalis</i> ATCC 29212 ¹ Range ²
Amoxicillin-K Clavulanate	Aug	≤1/0.5	N/A
Ampicillin	Am	≤0.06-0.25	N/A
Azithromycin	Azi	≤0.12-0.25	N/A
Cefepime	Cpe	≤0.25	N/A
Cefotaxime	Cft	≤0.25	N/A
Ceftriaxone	Cax	≤0.25	N/A
Cefuroxime	Crm	≤0.25-1	N/A
Chloramphenicol	C	≤2-8	N/A
Clarithromycin	Cla	≤0.12	N/A
Clindamycin	Cd	≤0.06-0.12	N/A
Daptomycin	Dap	≤0.25-0.5	N/A
Erythromycin	E	≤0.06-0.12	N/A
Levofloxacin	Lvx	≤0.5-2	N/A
Linezolid	Lzd	≤1-2	N/A
Meropenem	Mer	≤0.25	N/A
Minocycline	Min	N/A	1->2
Moxifloxacin	Mxf	≤0.25	N/A
Penicillin	P	0.25-1	N/A
Pristinamycin	Prs	≤1	N/A
Rifampin	Rif	≤0.5	N/A
Tetracycline	Te	≤1	N/A
Trimethoprim-Sulfamethoxazole	T/S	≤0.25/4.75-1/19	N/A
Vancomycin	Va	≤0.25-0.5	N/A

1. American Type Culture Collection, Manassas, VA USA.

2. Expected Ranges based on the CLSI M100 document.

PERFORMANCE CHARACTERISTICS

The performance characteristics of the MICroSTREP *plus* panel was established at multiple clinical laboratories. The antimicrobial agents were tested on a MICroSTREP *plus* panel using the turbidity standard technique and were read manually. A comparison was made to a reference microdilution panel prepared according to CLSI recommended method. MIC values were considered to be in essential agreement when the MICroSTREP *plus* panel was within $\pm$ one two-fold dilution of the reference method. Categorical Agreement occurs when the MICroSTREP *plus* panel and the CLSI reference interpretive results agree (S, I or R). There may be instances when a Categorical Agreement falls below the Essential Agreement.

This may occur when a significant number of organisms' MICs cluster around a category breakpoint resulting in interpretive errors. A high corresponding Essential Agreement percentage demonstrates that a $\pm$ 1 dilution fluctuation at the breakpoint can occur within an acceptable overall performance if category errors are due in large part to minor errors.

REPRODUCIBILITY

Reproducibility studies were completed at multiple clinical sites to confirm acceptable performance with the recommended methods described in the Procedural Manual.

Percent Agreement to Reference Method (CLSI)

Antimicrobial Agent	# Tested	MIC	Categorical
Amoxicillin-K Clavulanate	296	99.0	97.3
Ampicillin	320	93.4	99.3
Azithromycin	562	99.3	98.2
Cefepime	563	96.3	96.6
Cefotaxime (meningitis) ¹	563	98.6	95.6
Cefotaxime (non-meningitis) ¹	563	98.6	98.2
Ceftriaxone (meningitis) ¹	563	97.5	94.7

*American Type Culture Collection, Manassas, VA USA

Antimicrobial Agent	# Tested	MIC	Categorical
Ceftriaxone (non-meningitis) ¹	563	97.5	97.9
Cefuroxime sodium (parenteral)	296	99.7	95.6
Cefuroxime axetil (oral) ¹	296	99.7	94.9
Chloramphenicol	562	100	99.6
Clarithromycin	562	98.6	98.0
Clindamycin	562	97.5	99.5
Daptomycin	376	96.5	97.9
Erythromycin	562	94.5	98.9
Levofloxacin	562	99.5	100
Linezolid	493	97.2	99.8
Meropenem	562	99.1	97.3
Moxifloxacin	127	99.2	99.2
Penicillin (meningitis) ¹	563	92.5	98.6
Penicillin (non-meningitis) ¹	563	92.5	97.3
Penicillin (oral) ¹	563	92.5	96.4
Rifampin	108	96.3	100
Tetracycline ²	563	92.9	99.5
Trimethoprim-Sulfamethoxazole	296	100	98.0
Vancomycin	562	99.5	100

1. "Meningitis", "non-meningitis", and "oral" refer to the *Streptococcus pneumoniae* breakpoints used in the performance calculation. All antimicrobial breakpoints are not specific to infection site or antimicrobial agent administration.

2. For tetracycline an MIC bias for tetracycline and *S. pneumoniae* with manual read methods tended to be one doubling dilution higher than the reference method.

Percent Agreement to Reference Method (EUCAST)

Antimicrobial Agent	# Tested	MIC	Categorical
Ampicillin	384	96.1	95.6
Azithromycin	257	98.1	99.2
Cefepime	366	95.9	94.3
Cefotaxime	192	97.4	93.8
Ceftriaxone	192	96.4	94.3
Cefuroxime sodium (parenteral)	387	99.2	96.4
Chloramphenicol	218	96.3	100
Clarithromycin	475	98.3	97.9
Clindamycin	346	96.5	99.1
Daptomycin	281	97.2	99.3
Erythromycin	475	94.9	99.2
Levofloxacin	215	99.1	99.1
Linezolid	398	96.7	99.7
Meropenem (meningitis) ¹	133	98.5	95.5
Meropenem (non-meningitis) ¹	91	95.6	100
Minocycline	398	98.5	97.2
Moxifloxacin	398	97.5	100
Penicillin (meningitis) ¹	496	91.7	97.2
Penicillin (non-meningitis) ¹	347	92.2	93.9
Pristinamycin (SFM)	835	98.4	97.6
Tetracycline	219	96.3	98.2
Trimethoprim-Sulfamethoxazole	218	97.7	99.1
Vancomycin	342	98.8	100

1. "Meningitis", "non-meningitis", and "oral" refer to the *Streptococcus pneumoniae* breakpoints used in the performance calculation. All antimicrobial breakpoints are not specific to infection site or antimicrobial agent administration.

WARRANTY STATEMENT

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

NOTICE TO USER


Any serious incident that has occurred in relation to MicroScan panels shall be reported to Beckman Coulter and the competent authority of the Member State in which the user and/or the patient is established.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	May 2022	Initial Release

Symbols Key

The following symbols may be used in the labeling of Beckman Coulter Microbiology products.

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1, clause 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1, clause 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1, clause 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1, clause 5.1.3
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1, clause 5.5.5
	In Vitro Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1, clause 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1, clause 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1, clause 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
	Info For USA Only	Indicates information that is required only for products intended to be sold in the US.	n/a
	Prescription Use Only	Recognized by the US FDA as an alternate to "Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner"	21CFR 801.109(b)(1)

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Biological Risks	Indicates that there are potential biological risks associated with the medical device. Note: This symbol is not to be confused with the "biohazard" sign intended to be used in the workplace.	ISO 15223-1, clause 5.4.1
	Danger	n/a	n/a
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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REFERENCES

1. Barry, A. L. 1970. The Antimicrobial Susceptibility Test, Principles and Practices. p. 26. Lea and Febiger. Philadelphia, PA.
2. Chitwood, L. A. 1969. Tube dilution antimicrobial susceptibility testing: efficacy of a microtechnique applicable to diagnostic laboratories. *Appl. Microbiol.* 17:707.
3. Jorgensen, J.H., M.A. Pfaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter, D.W. Warnock (eds), 2015. Manual of Clinical Microbiology, 11th ed. American Society for Microbiology, Washington D.C.
4. Ericsson, H. M., and J. C. Sherris. 1971. Antibiotic sensitivity testing: Report of an international collaborative study. *Acta Path. Microbiol. Scand. Suppl.* 217:1.
5. Gerlach, E. H. 1974. Microdilution I: A comparative study, p. 63. In A. Balows (ed) Current Techniques for Antibiotic Susceptibility Testing. C. C. Thomas, Springfield, IL.
6. Harwick, J. H., P. Weiss, and F. Fekety, Jr. 1968. Application of microtitration techniques to bacteriostatic and bacteriocidal antibiotic susceptibility testing. *J. Lab. Clin. Med.* 72:511.
7. Jorgensen, J.H. 1994. Detection of antimicrobial resistance in *Streptococcus pneumoniae* by use of standardized susceptibility testing methods and recently developed interpretive criteria. *Clin. Microbiol. Newsl.* 16: 97-101.
8. Jorgensen, J.H. and G.V. Doern. 1995. Practical guidelines for in vitro susceptibility testing of selected gram-positive bacteria. *Clin. Microbiol. Newsl.* 17: 81-85.
9. Jorgensen, J.H., J.M. Swenson, F.C. Tenover, M.J. Ferraro, J.A. Hindler and P.R. Murray. 1994. Development of interpretive criteria and quality control limits for broth microdilution and disk diffusion antimicrobial susceptibility testing of *Streptococcus pneumoniae*. *J. Clin. Microbiol.* 32: 2448-2459.
10. Methods for Dilution Antimicrobial Susceptibility tests for Bacteria That Grow Aerobically. 2015. Approved Standard M07-A10. Clinical and Laboratory Standards Institute, Wayne, PA.
11. Tilton, R. C. and L. Newberg. 1974. Standardization of the microdilution susceptibility test, p. 77. In A. Balows (ed) Current Techniques for Antimicrobial Susceptibility Testing. C. C. Thomas, Springfield, IL.
12. Washington, J. A., E. Warren y A. G. Karlson. 1973. Stability of barium sulfate turbidity standards. *Appl. Microbiol.* 24:1013.
13. World Health Organization 1961. Standardization of methods for conducting microbic sensitivity tests. Second Report of the Expert Committee on Antibiotics. Technical Report Series #210, p. 1. Geneva.
14. EUCAST: European Committee on Antimicrobial Susceptibility Testing, ESCMID: European Society of Clinical Microbiology and Infectious Diseases, Version 6.0, Version 8.1, Version 9.0, Version 10.0 http://www.escmid.org/research_projects/eucast/clinical_breakpoints/.
15. Performance Standards for Antimicrobial Susceptibility Testing. 2020. CLSI Document M100, 30thed.. Clinical and Laboratory Standards Institute, Wayne, PA.
16. Comité de l'antibiogramme de La Société Française de Microbiologie. Recommandations 2019 V.2.0 Mai. Société Française de Microbiologie, Paris.

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.