MicroScan

MICroSTREP *plus* 3 Procedural Manual, Multi-regional

REF B1016-89A

TABLE OF CONTENTS

English	2
Français (FR)	10
Deutsch (DE)	19
Italiano (IT)	28
Español (ES)	37
Português Portugal (PT-PT)	46
Norsk (NO)	55
Magyar (HU)	63
Русский (RU)	72
Română (RO)	82
Србија (SR)	90
Tiếng Việt (VI)	98
REFERENCES	106

MicroScan

MICroSTREP *plus* 3 Procedural Manual, Multi-regional

REF B1016-89A

INTENDED USE

For Use With MicroScan MICroSTREP *plus* Type 3 Panels. MicroScan MICroSTREP *plus* Panels are designed for use in determining antimicrobial agent susceptibility of aerobic non-enterococcal streptococci (including *Streptococcus pneumoniae*) and *Haemophilus* spp.

SUMMARY AND PRINCIPLE

The antimicrobial susceptibility tests are miniaturizations of the broth dilution susceptibility test that have been dehydrated. Various antimicrobial agents are diluted in water or buffer to concentrations bridging the range of clinical interest. Panels are rehydrated with either Mueller-Hinton broth supplemented with 2-5% lysed horse blood (LHB) (for streptococci) or Haemophilus Test Medium (HTM) broth (for *Haemophilus* spp.), after inoculation of the broth with a standardized suspension of the organism.

After incubation in a non-CO₂ incubator for 20-24 hours, the minimum inhibitory concentration (MIC) for the test organism is determined by observing the lowest antimicrobial concentration showing inhibition of growth.^{1,2,3,4,5,6,7,8,9,10,11,12,13,14}

The procedural manual content listed in the labeling (e.g., interpretive criteria, panel read times, limitations) may differ from the criteria in the LabPro Information Manager due to differences in software and panel releases.

MICroSTREP *plus* Panel Type 3 (B1016-89A)

Antimicrobial Agents	Abbr.	Dilutions µg/mL
Amoxicillin/K Clavulanate	Aug	0.25/0.12-4/2
Ampicillin	Am	0.03-8
Azithromycin	Azi	0.25-4
Cefaclor	Cfr	1-2, 8-16
Cefepime	Cpe	0.12-2
Cefotaxime	Cft	0.03-8
Cefuroxime	Crm	0.25-8
Chloramphenicol	C	1-8
Ciprofloxacin	Cp	0.06-2
Clindamycin	Cd	0.25-2
Erythromycin	E	0.25-4
Meropenem	Mer	0.06-4
Penicillin	P	0.03-8
Rifampin	Rif	0.5-2
Tetracycline	Te	1-4
Trimethoprim/Sulfamethoxazole	T/S	0.25/4.75-2/38
Vancomycin	Va	0.5-8

WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with special awareness that the inoculated panels contain potentially pathogenic organisms.
3. This material contains infectious agents and should be discarded as appropriate for biohazardous waste.
4. Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.
5. Broth used to rehydrate panel contains 2-5% Lysed Horse Blood.

STORAGE

Store MicroScan MICroSTREP *plus* Panels at 2-25°C.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckmancoulter.com/techdocs

EVIDENCE OF DETERIORATION

Prolonged exposure to storage conditions other than those recommended may result in loss of potency of the antimicrobial agents. Do not use beyond the expiration date. Contact your Beckman Coulter Representative or Distributor for further assistance.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens should be collected, transported, and placed on primary isolation media according to procedures recommended in the *Manual of Clinical Microbiology*.³

MATERIALS PROVIDED

MICroSTREP *plus* Panel Type 3 (B1016-89A)

Procedural and QC Manual

Panel Worksheet

MATERIALS REQUIRED, BUT NOT PROVIDED

0.5 McFarland Turbidity Standard

100µL pipettor with disposable sterile tips

Cover Trays (B1010-56B or equivalent)

General Laboratory Equipment

Inoculator-D Set, 60/pk (B1013-5) or 240/pk (B1013-4)

Inoculum Water, 3 mL (B1015-2)

Inoculum Saline, 3 mL (B1015-12)

Haemophilus Test Medium, 25mL (B1015-26)

Mueller-Hinton Broth with 3% lysed horse blood, 25 mL (B1015-25)

Microdilution Viewer

Quality Control Organisms (Refer to QC section of this manual)

Quality Control Report Form (B1014-248A)

RENOK Rehydrator/Inoculator System (B1018-14) or equivalent

Turbidity Meter

Vortex

PROCEDURE

Broth Preparation

1. Allow the Mueller-Hinton Broth with 3% lysed horse blood and/or the Haemophilus Test Medium broth to come to room temperature before use.
2. Prior to inoculating panels, ensure that the broth has not been compromised by improper storage or handling by examining it for the following properties:

For LHB broth:

- A. Broth should be clear and red to dark red in color.
- B. The broth may become slightly brown in color as it gets closer to expiration date but should remain clear. The brownish color will not affect panel results.
- C. Do not use any broth that is turbid and discolored.

For HTM broth:

- A. Broth should be clear and golden yellow in color.
- B. A small amount of brown/black precipitate in the bottom of the tubes is normal as the broth gets closer to the expiration date. This precipitate should go back into solution once the broth is at room temperature and has been inverted a few times.
- C. Do not use any broth that is turbid and discolored.

Panel Preparation

1. Remove the panels to be used from storage. Panels should not be used if any of the following conditions exist:
 - A. Desiccant is not present or is broken.
 - B. Panel wells are discolored.

- C. Integrity of the packaging is compromised (unsealed, punctured, or torn).
2. Cut open the pouch and remove the panel. If stored in the refrigerator, remove the panel immediately from the foil pouch. Allow panels to equilibrate to room temperature prior to rehydration. Panels may be stacked with a clean Cover Tray on top. All opened panels should be used within the same day or discarded.

Inoculum Preparation

CLSI recommends periodically checking inoculum densities by performing colony counts.

1. Turbidity Standard Technique

THE TURBIDITY STANDARD TECHNIQUE IS THE ONLY METHOD APPROVED FOR USE WITH THIS PANEL.

- A. Using a sterile applicator stick, loop or swab, touch the surface of at least 4-5 large or 5-10 small, morphologically similar well isolated colonies or enough to create a 0.5 McFarland inoculum from a 16-20 hour (streptococci) or a 18-24 hour (*Haemophilus* spp.) noninhibitory agar plate (e.g., Trypticase Soy Agar with 5% Sheep Blood, Chocolate Agar). Do not prepare inoculum from plates with pin point colonies that exhibit autolysis. Incubation in 5-10% CO₂ is recommended.
- B. Emulsify in 3 mL of Inoculum Water (autoclaved deionized water) or Inoculum Saline.
- C. Cap tightly and vortex the suspension for 2-3 seconds. The final turbidity should be equivalent to that of a 0.5 McFarland Turbidity Standard. An equivalent turbidity can be achieved by using a MicroScan Turbidity Meter with a range of 0.08 ± 0.02.
- D. Pipet 0.1 mL (100 µL) of the standardized suspension into 25 mL of either LHB broth (for streptococci) or HTM broth (for *Haemophilus* spp.). Cap tightly and invert 8-10 times to mix.

Panel Rehydration/Inoculation

Inoculation of the panels is performed using the RENOK Rehydrator/Inoculator System with Inoculator-D sets. If an alternate system is used, rehydrate each well with 115 µL ± 10 µL of the appropriate suspension.

To ensure viability and purity of the organism tested, a purity plate should be prepared by streaking the inoculum to an appropriate agar plate and incubating under the appropriate conditions. If two or more colony types are present on the purity plate, re-isolate the colonies and retest.

1. Remove the transfer lid from the inoculator set. Pour the inoculated 25 mL broth (either LHB or HTM) into the seed trough.
2. Replace the transfer lid. Gently tap the transfer lid in all four corners to ensure absence of bubbles. Allow the transfer lid to equilibrate for a minimum of 20 seconds in the seed trough.
3. Attach the RENOK Rehydrator/Inoculator to the transfer lid and draw up the inoculum. Refer to the RENOK Rehydrator/Inoculator System Operator's Manual for operating instructions.

NOTE: For Mueller-Hinton broth with lysed horse blood (streptococci) only: due to the physical properties of the lysed horse blood supplement, a priming step is recommended to minimize the possibility of uneven filling of the panel wells. To perform the priming step, draw up the inoculum as per above. While the transfer lid is still seated on the seed trough, press the center release button to dispense the inoculum back into the seed trough. Draw up the inoculum a second time into the transfer lid.

4. Place the RENOK Rehydrator/Inoculator transfer lid apparatus on the MICroSTREP *plus* Panel and release the inoculum into the panel.

Incubation

1. To ensure even thermal distribution during incubation, stack the panels in groups of 3-5.
2. Place a clean Cover Tray on top of each group of panels to prevent evaporation. Cover Trays may be reused. Do not decontaminate the Cover Trays with alcohol. They may be cleaned with soap and water. Rinse well and allow to air dry.
3. Incubate the panels for 20-24 hours at 35°C ± 1°C in a non-CO₂ incubator.

Reading the Panels

Panels are read manually using indirect light (such as that provided by the MicroScan Microdilution Viewer) and results can be recorded on a panel worksheet.

1. Following 20-24 hours incubation, remove the panels from the incubator.
2. Wipe the bottom of the panel with a lint-free tissue to remove any condensation or debris that may be present.
3. Read the panels only if the growth control well is turbid. Inadequate growth is defined as a slight haze in the well or if the broth is clear.
4. Record results on the appropriate worksheet.
5. Reading Antimicrobial Susceptibilities (MICs)

- A. Growth in the wells appears as turbidity. Turbidity may take the form of a haze throughout the well, a button in the center of the well, or a fine granular growth throughout the well. A change in the broth color should not be interpreted as growth. Turbidity is read using indirect light.
- B. Record MIC results as follows:
 - a. Following 20-24 hours incubation, record the MIC as the lowest antimicrobial concentration showing inhibition of growth.
 - b. When growth occurs in all concentrations of an antimicrobial, the MIC is recorded as greater than (>) the highest concentration.
 - c. When no growth occurs in any of the concentrations of antimicrobics, the MIC is recorded as less than or equal to ($\leq$) the lowest concentration.
 - d. A clear well in a series of growth wells (e.g., growth at 1, 2, and 8 $\mu\text{g/mL}$, but not at 4 $\mu\text{g/mL}$) is called a skipped well and should be ignored.
 - e. Spot growth in isolated wells indicates contamination. The test should be repeated.
 - f. A “trailing effect” may be observed in some drug/organism combinations. Trailing with trimethoprim/sulfamethoxazole (T/S) with the use of the RENOK Rehydrator/Inoculator System is due to the inoculum concentration. The end point should be read as the lowest concentration which when compared to the growth well shows approximately 80% reduction of growth.
 - g. Rare β -lactamase negative, ampicillin resistant (BLNAR) strains of *H. influenzae* should be considered resistant to amoxicillin/clavulanic acid, ampicillin/sulbactam, cefaclor, cefamandole, cefetamet, cefonicid, cefprozil, cefuroxime, loracarbef, and piperacillin/tazobactam despite apparent *in vitro* susceptibility of some BLNAR strains to these agents.¹¹

Typical Indicators

Large button in Growth well.

Typical growth buttons in columns 1 and 2;

Column 2 demonstrates a “skip” in growth wells, which should be ignored.

This single button in column 3 bounded by clear wells above and below indicates spot contamination and should be ignored.

MIC in column 1 = 2 $\mu\text{g/mL}$ (lowest concentration clear well).


MIC in column 2 = > 16 $\mu\text{g/mL}$ (growth in all wells).

MIC in column 3 = ≤ 0.25 $\mu\text{g/mL}$ (no growth in all wells).

MIC in column 4 = 2 $\mu\text{g/mL}$ (Trailing effect in wells 2-8).

MIC in column 5 = ≤ 0.25 $\mu\text{g/mL}$ (Trailing effect in all wells, thus MIC is $\leq$). (Typical effect of some drug/bug combinations).

ANTIMICROBIC WELLS $\mu\text{g/mL}^*$

AO_0092B3

*Diagram: 1 = GROWTH and 2 = COLUMN

RESULTS

Interpretation of MIC Results

Susceptibility is determined by comparing the MIC of an organism to the attainable blood or urine level of the antimicrobial. The following tables list the streptococci interpretive criteria as indicated in the CLSI Document M100-S26.¹⁵

1. Streptococci Interpretive Criteria

Streptococci Interpretive Criteria			
Antimicrobial Agents	Susceptible	Intermediate	Resistant
Amoxicillin/K Clavulanate			
<i>S. pneumoniae</i> (non-meningitis)	$\leq 2/1$	4/2	$\geq 8/4$
Ampicillin			
β -hemolytic streptococci ¹	≤ 0.25	-	-
Viridans streptococci ²	≤ 0.25	0.5-4	≥ 8
Azithromycin ³			
Streptococci	≤ 0.5	1	≥ 2
Cefaclor			
<i>S. pneumoniae</i>	≤ 1	2	≥ 4

Streptococci Interpretive Criteria

Antimicrobial Agents	Susceptible	Intermediate	Resistant
Cefepime			
<i>S. pneumoniae</i> (non-meningitis)	≤1	2	≥4
β-hemolytic streptococci ¹	≤0.5	-	-
Viridans streptococci ²	≤1	2	≥4
Cefotaxime			
<i>S. pneumoniae</i> (meningitis) ^{4,6}	≤0.5	1	≥2
<i>S. pneumoniae</i> (non-meningitis) ^{4,6}	≤1	2	≥4
β-hemolytic streptococci ¹	≤0.5	-	-
Viridans streptococci ²	≤1	2	≥4
Cefuroxime axetil (oral)			
<i>S. pneumoniae</i>	≤1	2	≥4
Cefuroxime sodium (parenteral)			
<i>S. pneumoniae</i>	≤0.5	1	≥2
Chloramphenicol ³			
<i>S. pneumoniae</i>	≤4	-	≥8
Streptococci other than <i>S. pneumoniae</i>	≤4	8	≥16
Clindamycin ³			
Streptococci	≤0.25	0.5	≥1
Erythromycin ^{3,5}			
Streptococci	≤0.25	0.5	≥1
Meropenem ^{3,6}			
<i>S. pneumoniae</i>	≤0.25	0.5	≥1
β-hemolytic streptococci ¹	≤0.5	-	-
Viridans streptococci ^{1,2}	≤0.5	-	-
Penicillin ⁶			
<i>S. pneumoniae</i> (meningitis)	≤0.06	-	≥0.12
<i>S. pneumoniae</i> (non-meningitis)	≤2	4	≥8
<i>S. pneumoniae</i> (oral penicillin V)	≤0.06	0.12-1	≥2
β-hemolytic streptococci ¹	≤0.12	-	-
Viridans streptococci ²	≤0.12	0.25-2	≥4
Rifampin			
<i>S. pneumoniae</i>	≤1	2	≥4
Tetracycline ⁷			
Streptococci	≤2	4	≥8
Trimethoprim/Sulfamethoxazole			
<i>S. pneumoniae</i>	≤0.5/9.5	1/19-2/38	≥4/76
Vancomycin ^{1,6}			
Streptococci	≤1	-	-

1. The absence of resistant strains precludes defining any results categories other than "susceptible." Strains yielding results suggestive of a "nonsusceptible" category should be submitted to a reference laboratory for further testing.¹⁵

2. Apply viridans streptococci breakpoints to *S. bovis* (Group D streptococci).

3. Report therapy for systemic source only.

4. For *S. pneumoniae* isolated from cerebral spinal fluid (CSF), report only meningitis interpretations. For *S. pneumoniae* isolated from all specimens other than CSF, report interpretations for both meningitis and non-meningitis.

5. Susceptibility and resistance to azithromycin and clarithromycin can be predicted by testing erythromycin.

6. Only results of testing with penicillin, cefotaxime, ceftriaxone, meropenem and vancomycin should be reported routinely for CSF isolates of *S. pneumoniae*.

7. Interpretations for *S. pneumoniae* are based on CLSI M100-S22.

NOTE: Results may be inaccurate if specified criteria are applied to organisms other than those listed.

NOTE: Streptococci susceptible to penicillin can be considered susceptible to other beta-lactam antimicrobics (e.g., cephalosporins, ampicillin, etc.). Refer to the General Comments in CLSI Tables 2G and 2H for additional information.

2. *Haemophilus* spp. Interpretive Criteria

Antimicrobial Agent	Susceptible	Intermediate	Resistant
Amoxicillin/K Clavulanate	≤4/2	-	≥8/4
Ampicillin ^{1,2}	≤1	2	≥4
Azithromycin ^{3,4}	≤4	-	-
Cefaclor	≤8	16	≥32
Cefepime ⁴	≤2	-	-
Cefotaxime ^{2,4}	≤2	-	-
Cefuroxime	≤4	8	≥16
Chloramphenicol ^{2,3}	≤2	4	≥8
Ciprofloxacin ⁴	≤1	-	-
Meropenem ^{2,3,4}	≤0.5	-	-
Rifampin	≤1	2	≥4

Antimicrobial Agent	Susceptible	Intermediate	Resistant
Tetracycline	≤2	4	≥8
Trimethoprim/Sulfamethoxazole	≤0.5/9.5	1/19-2/38	≥4/76

1. Rare β-lactamase negative, ampicillin resistant (BLNAR) strains of *H. influenzae* should be considered resistant to amoxicillin/clavulanic acid, ampicillin/sulbactam, cefaclor, cefamandole, cefetamet, cefonicid, cefprozil, cefuroxime, loracarbef, and piperacillin/tazobactam despite apparent *in vitro* susceptibility of some BLNAR strains to these agents.¹⁵
2. Only results of testing with ampicillin, cefotaxime, chloramphenicol and meropenem should be reported routinely for CSF isolates of *H. influenzae*.
3. Report therapy for systemic source only.
4. The absence of resistant strains precludes defining any results categories other than "susceptible." Strains yielding results suggestive of a "nonsusceptible" category should be submitted to a reference laboratory for further testing.¹⁵

NOTE: Results may be inaccurate if specified criteria are applied to organisms other than *Haemophilus* spp.

LIMITATION OF THE PROCEDURE

1. Some fastidious organisms may require additional supplements and may not grow on the MICroSTREP *plus* panel. If there is no growth in the Growth Well, the results of the panel are not valid.
2. The MICroSTREP *plus* Panels are only intended for testing all non-enterococcal streptococci (including *S. pneumoniae*) and *Haemophilus* spp.
3. A slight haziness may occur in random wells of some antimicrobics. This should not be interpreted as growth.
4. The interpretation of MIC results require trained clinical personnel who should use judgement and knowledge prior to accepting the susceptibility of an organism.
5. The MICroSTREP *plus* Panel should not be incubated in a CO₂ incubator.
6. Results obtained with mixed cultures are not valid.
7. β-lactamase testing is recommended by the CLSI for predicting *Haemophilus* spp. resistance to penicillin, ampicillin, and amoxicillin.¹¹ Even though a β-lactamase test is not included on this panel, this test should be performed routinely on *Haemophilus* spp., since results are utilized in the categorical interpretation of these antimicrobial agents.
8. *Haemophilus* spp. and rifampin/rifampicin should not be reported and if critical to patient care an alternative testing method should be performed. When printing is not suppressed by LabPro, the result should be manually suppressed. Consult the LabPro Operator's Manual for Instructions.

QUALITY CONTROL

The acceptability of the antimicrobial agents should be checked by testing organisms with known MIC ranges. The results for the recommended American Type Culture Collection (ATCC) control organisms are listed below.*

MICroSTREP *plus* Type 3 (B1016-89A)

Antimicrobial Agent	Abbr.	<i>S. pneumoniae</i> ¹ ATCC 49619 Range ³	<i>H. influenzae</i> ² ATCC 49247 Range	<i>H. influenzae</i> ² ATCC 49766 Range
Amoxicillin/K Clavulanate	Aug	≤0.25/0.12	2/1->4/2	-
Ampicillin	Am	0.06-0.25	2-8	-
Azithromycin	Azi	≤0.25	0.5-4 ⁴	-
Cefaclor	Cfr	≤1-2, 8	-	≤1-2, 8
Cefepime	Cpe	≤0.12-0.25	0.5-2	-
Cefotaxime	Cft	≤0.03-0.12	0.12-0.5	-
Cefuroxime	Crm	≤0.25-1	-	≤0.25-1
Chloramphenicol	C	2-8	≤1	-
Ciprofloxacin	Cp	-	≤0.06	-
Clindamycin	Cd	≤0.25	-	-
Erythromycin	E	≤0.25	-	-
Meropenem	Mer	≤0.06-0.25	-	≤0.06-0.12
Penicillin	P	0.25-1	-	-
Rifampin	Rif	≤0.5	≤0.5-1	-
Tetracycline	Te	≤1	4->4	-
Trimethoprim/Sulfamethoxazole	T/S	≤0.25/4.75-1/19	≤0.25/4.75	-
Vancomycin	Va	≤0.5	-	-

1. Based on panel rehydration with Mueller-Hinton broth with 2-5% lysed horse blood.
2. Based on panel rehydration with *Haemophilus* Test Medium broth.
3. Range = Expected value (μg/mL)
4. This range does not agree with the Acceptable Quality Control Range listed in the CLSI M100 Document.

*American Type Culture Collection, Manassas, VA USA

PERFORMANCE CHARACTERISTICS

The performance of the MicroScan MICroSTREP *plus* Panel has been established in evaluations at several clinical laboratories. Data were analyzed and interpreted according to MIC interpretive standards in NCCLS Standard M7-A6 and/or by percent agreement to a reference method.

WARRANTY STATEMENT


The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	June 2019	Creation
AB	August 2019	Added Translations
AC	March 2021	Added Do Not Report for Rifampin vs <i>Haemophilus</i> . Updated Summary and Principle section by adding information regarding differences between LabPro software and panel releases.

Symbols Key

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1; 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1; 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1; 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1; 5.3.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1; 5.5.5

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	In Vitro Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1; 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1; 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1; 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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REFERENCES

1. Barry, A. L. 1970. The Antimicrobial Susceptibility Test, Principles and Practices. p. 26. Lea and Febiger. Philadelphia, PA.
2. Chitwood, L.A. 1969. Tube dilution antimicrobial susceptibility testing: efficacy of a microtechnique applicable to diagnostic laboratories. Appl. Microbiol. 17:707.
3. Jorgensen, J.H., M.A. Pfaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter, D.W. Warnock (eds), 2015. Manual of Clinical Microbiology, 11th ed. American Society for Microbiology, Washington D.C.
4. Ericsson, H. M., and J. C. Sherris. 1971. Antibiotic sensitivity testing: Report of an international collaborative study. Acta Path. Microbiol. Scand. Suppl. 217:1.
5. Gerlach, E. H. 1974. Microdilution I: A comparative study, p. 63. In A. Balows (ed.), Current Techniques for Antibiotic Susceptibility Testing. C. C. Thomas, Springfield, IL.
6. Harwick, J. H., P. Weiss, and F. Fekety, Jr. 1968. Application of microtitration techniques to bacteriostatic and bacteriocidal antibiotic susceptibility testing. J. Lab. Clin. Med. 72:511.
7. Jorgensen, J.H. 1994. Detection of antimicrobial resistance in *Streptococcus pneumoniae* by use of standardized susceptibility testing methods and recently developed interpretive criteria. Clin. Microbiol. Newsl. 16: 97-101.
8. Jorgensen, J.H. and G.V. Doern. 1995. Practical guidelines for *in vitro* susceptibility testing of selected gram-positive bacteria. Clin. Microbiol. Newsl. 17: 81-85.
9. Jorgensen, J.H., J.M. Swenson, F.C. Tenover, M.J. Ferraro, J.A. Hindler and P.R. Murray. 1994. Development of interpretive criteria and quality control limits for broth microdilution and disk diffusion antimicrobial susceptibility testing of *Streptococcus pneumoniae* J. Clin. Microbiol. 32: 2448-2459.
10. Methods for Dilution Antimicrobial Susceptibility tests for Bacteria That Grow Aerobically. 2003. Approved Standard M7-A6. National Committee for Clinical Laboratory Standards, Wayne, PA.
11. Methods for Dilution Antimicrobial Susceptibility tests for Bacteria That Grow Aerobically. 2015. Approved Standard M07-A10. Clinical and Laboratory Standards Institute, Wayne, PA.
12. Tilton, R. C., and L. Newberg. 1974. Standardization of the microdilution susceptibility test, p. 77. In A. Balows (ed.), Current techniques for antimicrobial susceptibility testing. C. C. Thomas, Springfield, IL.
13. Washington, J. A., E. Warren, and A. G. Karlson. 1973. Stability of barium sulfate turbidity standards. Appl. Microbiol. 24:1013.
14. World Health Organization 1961. Standardization of methods for conducting microbic sensitivity tests. Second Report of the Expert Committee on Antibiotics. Technical Report Series #210, p. 1. Geneva.
15. Performance Standards for Antimicrobial Susceptibility Testing. 2016. CLSI Document M100-S26. Clinical and Laboratory Standards Institute, Wayne, PA.

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda
Alameda Rio Negro, 500, 15º andar, Torre B Alphaville Industrial
CEP 06.454-00 Barueri, São Paulo, Brasil
CNPJ: 42.160.812/0001-44 Telephone: 0800-771-8818

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.