MicroScan

Rapid Anaerobe ID Procedural Manual

REF B1017-2

TABLE OF CONTENTS

English	2
Français (FR)	8
Deutsch (DE)	15
Italiano (IT)	22
Español (ES)	29
Português Portugal (PT-PT)	36
Norsk (NO)	43
Ελληνικά (EL)	49
中文 (ZH-CN)	57
Polski (PL)	63
한국어 (KO)	70
Türkçe (TR)	76
Русский (RU)	82
Português Brasil (PT-BR)	90
Tiếng Việt (VI)	97
REFERENCES	104

MicroScan

Rapid Anaerobe ID Procedural Manual

REF B1017-2

INTENDED USE

The MicroScan Rapid Anaerobe Identification panel is designed for rapid identification of anaerobic bacteria isolated from human clinical specimens.

SUMMARY AND PRINCIPLE

Modified conventional and chromogenic tests are used for identification of anaerobic bacteria.^{1,2,3,4,5,6,7} A McFarland 3 suspension of the organism to be tested is made in autoclaved deionized water. The MicroScan Rapid Anaerobe Identification Panel, consisting of 24 dehydrated substrates, is rehydrated and inoculated with the bacterial suspension and incubated aerobically for 4 hours at 35-37°C. Because the MicroScan Rapid Anaerobe Identification Panel detects preformed enzymes, growth of the organism in the panel is not necessary. After incubation, changes in pH, utilization of certain substrates or other biochemical reactions are interpreted for identification.

WARNING AND PRECAUTIONS

1. MicroScan Rapid Anaerobe Identification panels and reagents are for *in vitro* diagnostic use only.
2. Do not allow Peptidase Reagent, Ehrlich's Reagent, Xylene, or any of the Identification substrates to come in contact with eyes, skin, or clothing.
3. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with special awareness that inoculated panels contain potentially pathogenic organisms.
4. All materials should be autoclaved prior to disposal.

STORAGE

Store panels and reagents at 2-8°C.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens, devoid of normal anaerobic flora, should be collected, transported, and placed on primary isolation media according to recommended procedures in the V.P.I., *Anaerobe Laboratory Manual*,² *Wadsworth Anaerobic Bacterial Manual*,⁸ C.D.C.'s *Laboratory Methods in Bacteriology*,⁹ and the *Manual of Clinical Microbiology*.¹⁰

NOTE: Colonies from selective media that contain antibiotics, i.e., LKV (Kanamycin/Vancomycin), or *Bacteroides* Bile Esculin (BBE) agar, Phenylethyl Alcohol (PEA), etc., should not be used to prepare inoculum for panels.

MATERIALS PROVIDED

Rapid Anaerobe Identification Panel (B1017-2)

20 Worksheets

Procedural Manual

MATERIALS REQUIRED, BUT NOT PROVIDED

Mineral Oil (B1010-40)

Cover Tray (B1010-56B)

Peptidase Reagent (B1012-30B)

0.5% N, N-Dimethyl-Alpha-Naphthylamine (B1010-45A)

0.8% Sulfanilic Acid (B1010-44A)

Ehrlich's reagent (B1015-6)

Xylene (B1015-5)

50 µL pipettor with sterile tips

Hydrogen Peroxide, 3%

Inoculum Water, 3mL autoclaved deionized water (B1015-2)

Sterile cotton-tipped applicators or sterile inoculating loops

Wooden applicator sticks

3, 4 and 5 McFarland Barium Sulfate Turbidity Standard

Vortex mixer

Incubator, non-CO₂ (35°C)

Quality Control Organisms:

C. perfringens ATCC 13124*

C. sordellii ATCC 9714

B. fragilis ATCC 25285

P. magnus ATCC 29328

PROCEDURE

Panel Preparation

1. Remove panels from refrigeration. **Panels should not be used if the desiccant is not present or is broken or if the integrity of the packaging is compromised (unsealed, punctured or torn).**
2. Cut open the pouch and remove the panel. Remove the panel immediately from the foil pouch. Allow panels to equilibrate to room temperature prior to rehydration. Panels may be stacked with clean Cover Tray on top.

Inoculum Preparation

1. Prior to inoculation of panels, organisms to be tested should be observed for aerotolerance, purity, Gram stain reaction and colony morphology.
2. Using a sterile cotton applicator or inoculating loop, suspend enough colonies from a 24-48 hour culture on non-selective media in a tube of 3mL Inoculum Water (autoclaved deionized water) to make a turbid suspension equivalent to at least a McFarland 3 standard but not greater than a McFarland 5 standard. Vortex to emulsify if necessary. Once the suspension is prepared, the panel should be inoculated within 15 minutes.**

Panel Inoculation

1. Remove panel from foil pouch and label it with the specimen number.
2. Add 50 µL of bacterial suspension to each substrates-containing well plus the two control wells (BNAC, NPC). Use of a Pasteur pipette for inoculation is not recommended.
3. After the panel has been inoculated with the 50 µL inoculum, lightly tap each side of the panel to help ensure mixing with the substrate.

Biochemical Oil Overlay

Overlay the Urea (URE) and Indole (IND) wells with 2-3 drops of mineral oil.

Incubation

1. Stack panels in groups of 3 to 5 with a cover panel over the top panel to prevent evaporation.
2. Incubate the inoculated panels for 4 hours at 35-37°C in a non-CO₂ incubator.

Reading the Panels

1. Remove the Peptidase Reagent from the refrigerator at least 30 minutes prior to reading the panels. The Peptidase Reagent may precipitate at 4°C, but it will redissolve at room temperature. Do not expose the Peptidase Reagent to elevated temperatures to accelerate the redissolving of the crystalline precipitate because high temperatures will cause the reagent to degrade.
2. The panels should be read using reflected light from a white background.

Table 1, Substrates

Substrates	Abbr.	Substrates	Abbr.
p-Nitrophenyl-β-D-Galactopyranoside	BGAL	L-Lysine-β-Naphthylamide (acid)	LYA
p-Nitrophenyl-α-D-Galactopyranoside	AGAL	Glycylglycine-β-Naphthylamide	GGLY
bis-p-Nitrophenyl-Phosphate	BPO ₄	Glycine-β-Naphthylamide	GLY
p-Nitrophenyl-N-Acetyl-β-D-Glucosaminide	NGLU	L-Proline-β-Naphthylamide	PRO
p-Nitrophenyl-α-D-Glucopyranoside	AGL	L-Arginine-β-Naphthylamide	ARG
o-Nitrophenyl-β-D-Glucopyranoside	BGL	L-Pyrrolidonyl-β-Naphthylamide	PYR
p-Nitrophenyl-Phosphate	PO ₄	L-Tryptophan-β-Naphthylamide	TRY
p-Nitrophenyl-α-L-Fucopyranoside	AFU	3-Indoxyl Phosphate	IDX
p-Nitrophenyl-α-D-Mannopyranoside	MNP	Trehalose	TRE
L-Leucine-β-Naphthylamide	LEU	Urea	URE
DL-Methionine-β-Naphthylamide	MET	Indole	IND
L-Lysine-β-Naphthylamide (alkaline)	LYB	Nitrate	NIT

*American Type Culture Collection, Manassas, VA USA

** Non-selective media include TSA, Brucella, brain heart infusion or Columbia agar, all with 5% sheep blood.

3. After 4 hours of incubation, record the results for all the substrates not requiring additional reagents: BGAL, AGAL, BPO₄, NGLU, AGL, BGL, PO₄, AFU, MNP, IDX, URE, and TRE. Compare BGAL, AGAL, BPO₄, NGLU, AGL, BGL, PO₄, AFU, and MNP to the Nitrophenyl Control (NPC) well. Any shade of yellow should be recorded as positive.
4. Add 1 drop of Peptidase Reagent to LEU, MET, LYB, LYA, GGLY, GLY, PRO, ARG, PYR, TRY, and BNAC wells.
 - A. Allow the color to develop for at least 30 seconds and then record the results. Do not wait longer than 3 minutes after adding the reagents before recording the results. The presence of any shade of pink or red should be interpreted as positive. The color on these wells will darken with time so that a negative yellow color may become a darker orange color, but this should still be interpreted as a negative reaction. Reactions should be recorded as positive **only** if there is a pink or red present. Some wells may be a darker orange compared to the control well (BNAC) but these should still be recorded as negative.
 - B. After addition of Peptidase Reagent to the β -Naphthylamide control well (BNAC), the well should be yellow in color. If there is any pink or red present, then the Peptidase Reagent should be discarded to prevent false positive interpretations.
5. Two control wells are provided to assist in interpreting the reactions.
 - A. BNAC = β -Naphthylamide control. After addition of Peptidase Reagent, this well should be yellow to orange in color.
 - B. NPC = Nitrophenyl control. This well should be colorless.
6. Add one drop of Xylene to the Indole (IND) well. With the wooden applicator stick, stir the contents of the IND well. Add 1 drop of Ehrlich's Reagent. Read after 2-3 minutes.
7. Add 1 drop each of Sulfanilic Acid to the Nitrate (NIT) well followed by 1 drop of N, N-Dimethyl- α -Naphthylamine. Read after 2-3 minutes.
8. Add 2-3 drops of 3% Hydrogen Peroxide to the BGAL well. Wait 2-5 minutes before recording the result. A positive reaction is indicated by the formation of large bubbles.

RESULTS

1. Biochemical Interpretation Chart

Substrate	Reagents	Positive	Negative
BGAL AGAL BPO ₄ NGLU AGL BGL PO ₄ ¹ AFU MNP	Compare to Nitrophenyl (NPC) control well	Any shade of yellow	Clear
LEU MET LYB LYA GGLY GLY PRO ARG PYR TRY	Add 1 drop of Peptidase reagent. Allow color reactions to develop for at least 30 seconds, but no longer than 3 minutes.	Any shade of Pink to Magenta	Yellow to Orange
BNAC	The β -Naphthylamide (BNAC) control well should be yellow.		
IDX		Blue color/Blue precipitate	Clear
TRE		Yellow	Orange/Red
URE		Pink-Red	Yellow-Orange
IND	Add 1 drop of Xylene, stir, add 1 drop of Ehrlich's reagent. Read after 2-3 minutes.	Pink	Yellow/Clear

Substrate	Reagents	Positive	Negative
NIT	Add 1 drop each 0.8% Sulfanilic Acid and 0.5% N, N-Dimethyl-Alpha-Naphthylamine. Read after 2-3 minutes.	Red	Clear to Pale Pink
Catalase	Add 2-3 drops of 3% Hydrogen Peroxide to BGAL. Wait 2-3 min. for bubbles to form.	Large Bubbles	No bubbles (or few tiny bubbles/same amount in all wells.)

1. A negative PO₄ reaction may be clear to very light yellow.

2. Principles of Identification Reactions

Nitrophenyl substrates (BGAL, AGAL, BPO₄, NGLU, AGL, BGL, PO₄, AFU, MNP): The ortho or para-nitrophenyl substrates are colorless. If the appropriate enzyme is present, the substrate is cleaved releasing the free o-or p-nitrophenyl, which is yellow.

Amino acid β-naphthylamide substrates (LEU, MET, LYB, LYA, GGLY, GLY, PRO, ARG, PYR, TRY): When the substrate is hydrolyzed by the appropriate enzyme (usually an amino acid arylamidase), β-naphthylamide is released. The β-naphthylamide is detected by the addition of p-dimethylaminocinnamaldehyde (in the Peptidase Reagent) which creates a complex that is pink to purple in color.

NOTE: The enzyme which cleaves an amino acid β-naphthylamide is an amino acid arylamidase. The term “aminopeptidase” is frequently used as a synonym for arylamidase but actually represents an enzyme with a different specificity (an aminopeptidase cleaves the N-terminal amino acid from a peptide). An arylamidase and an aminopeptidase may hydrolyze the same amino acid β-naphthylamide substrate. Thus, the panel is not necessarily measuring different and unique enzymes. A single arylamidase may hydrolyze several amino acid β-naphthylamides.

Indoxyl Phosphate (IDX): Indoxyl phosphate is cleaved by a phosphatase to release indoxyl. The indoxyl combines with oxygen to form indigo blue which is an insoluble blue or blue-green precipitate.

Trehalose (TRE): The utilization of Trehalose results in acid formation with a subsequent drop in pH which is detected by the phenol red indicator changing from orange to yellow.

Urea (URE): Urea is split by urease to form ammonia and carbon dioxide. The ammonia (in the form of ammonium carbonate) causes a rise in pH which is detected by the phenol red indicator changing from yellow to pink.

Indole (IND): The metabolism of tryptophan results in the formation of indole. Xylene extracts the indole from the broth and concentrates it. The Ehrlich's Reagent reacts with the indole to form a pink color.

Nitrate (NIT): The ability of an organism to reduce nitrate to nitrite is detected by the addition of Sulfanilic Acid and then N, N-Dimethyl-Alpha-Naphthylamine, which produce a red color in the presence of nitrite.

Catalase breaks down hydrogen peroxide to water and oxygen gas, which is liberated as bubbles.

ORGANISM IDENTIFICATION

The MicroScan Rapid Anaerobe Biotype Codebook is used for the identification of unknown test organisms. This codebook has been generated from the Beckman Coulter data base for the tests included on the identification panel.

The MicroScan Rapid Anaerobe Codebook is based on a computer analysis of the 24 tests included on the Rapid Anaerobe identification panel. The test results are transformed into a 8 digit “biotype number” for which the codebook lists a species identification and a cumulative relative probability of identification. All possible identifications are printed in descending order of highest probability to a cumulative of 99.9%. Should a biotype number occur that cannot be found in the codebook, one should first suspect a procedural error. Usually the biotype number has been added incorrectly or a reaction has been erroneously recorded. If the biotype number is correct, then the possibility of a mixed culture should be checked and the organism should be re-tested, with special emphasis on using a sufficiently turbid cell suspension. If retesting yields the same results, consult the Biotype Lookup Service on the Beckman Coulter website or consult your local distributor.

LIMITATION OF THE PROCEDURE

1. Organisms to be tested should be cultured on non-selective media such as TSA, Brucella, brain heart infusion, or Columbia agar, all with 5% sheep blood. Selective media that contain antibiotics, i.e., LKV (Kanamycin/Vancomycin), *Bacteroides* Bile Esculin (BBE) agar, Phenylethyl Alcohol (PEA) agar, etc., should not be used.
2. The turbidity of the organism suspension must be equal to at least a McFarland 3 standard but not greater than McFarland 5 standard. If the suspension of the organism to be tested is less turbid than a McFarland 3 standard, weak reactions or false negative reactions may occur.
3. The Rapid Anaerobe Identification panel is designed to identify only anaerobic bacteria. Prior to testing, each organism should be tested for aerotolerance, purity, Gram stain reaction and colonial morphology. Should more than one species identification be list for a particular biotype number, supplementary tests may be necessary for the correct identification.
4. Not all anaerobic species encountered in clinical situations are included in the Rapid Anaerobe database. The species identified on the Rapid Anaerobe panel are listed in the Rapid Anaerobe Biotype Codebook. Anaerobic identifications

should be correlated with gram stains and other supporting diagnostic tests. The user should consult current reference literature for these tests and also for current bacterial nomenclature.

5. The MicroScan Rapid Anaerobe Identification Panel has not been cleared by the FDA for use with the MicroScan Instrument Systems. This panel must be read manually.

QUALITY CONTROL

The acceptability of the MicroScan Rapid Anaerobe Identification media should be checked by testing organisms with known reactions. Negative reactions may be achieved by rehydrating a panel with 50 µL of autoclaved deionized water. All substrates should be negative as described in the interpretation chart. Culture purity should be checked by subculturing from the suspension used as the inoculum onto suitable agar media. Test results are invalid if a mixed culture is detected (See Table 2).*

Table 2
Quality Control Chart

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1. American Type Culture Collection, Manassas, VA USA.

2. For optimum results, the *C. perfringens* culture should be only 24 hours old.


PERFORMANCE CHARACTERISTICS

In clinical trials using anaerobic organisms representing a diverse set of recent clinical isolates, the Rapid Anaerobe Identification panel correctly identified 92% of the isolates.

Symbols Key

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
Ⓢ	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1; 5.4.2
🕒	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1; 5.1.4

*The results for each reaction for the recommended American Type Culture Collection (ATCC) control organisms are listed in the Quality Control Chart.

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1; 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1; 5.3.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1; 5.5.5
	<i>In Vitro</i> Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1; 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1; 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1; 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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Tabelle 1, Substrate

Substrate	Abk.	Substrate	Abk.
p-Nitrophenyl-β-D-Galactopyranosid	BGAL	L-Lysin-β-Naphthylamid (Säure)	LYA
p-Nitrophenyl-α-D-Galactopyranosid	AGAL	Glycylglycin-β-Naphthylamid	GGLY

*American Type Culture Collection, Manassas, VA, USA

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Substrate	Abk.	Substrate	Abk.
bis-p-Nitrophenylphosphat	BPO ₄	Glycin-β-Naphthylamid	GLY
p-Nitrophenyl-N-acetyl-β-D-glucosaminid	NGLU	L-Prolin-β-Naphthylamid	PRO
p-Nitrophenyl-α-D-Glucopyranosid	AGL	L-Arginin-β-Naphthylamid	ARG
o-Nitrophenyl-β-D-glucopyranosid	BGL	L-Pyrrolidonyl-β-Naphthylamid	PYR
p-Nitrophenylphosphat	PO ₄	L-Tryptophan-β-Naphthylamid	TRY
p-Nitrophenyl-α-L-Fucopyranosid	AFU	3-Indoxylphosphat	IDX
p-Nitrophenyl-α-D-Mannopyranosid	MNP	Trehalose	TRE
L-Leucin-β-Naphthylamid	LEU	Harnstoff	URE
DL-Methionin-β-Naphthylamid	MET	Indol	IND
L-Lysin-β-Naphthylamid (alkalisch)	LYB	Nitrat	NIT

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Bảng 1, Cơ chất

Cơ chất	Chữ viết tắt	Cơ chất	Chữ viết tắt
p-Nitrophenyl-β-D-Galactopyranoside	BGAL	L-Lysine-β-Naphthylamide (axit)	LYA
p-Nitrophenyl-α-D-Galactopyranoside	AGAL	Glycylglycine-β-Naphthylamide	GGLY
bis-p-Nitrophenyl-Phosphate	BPO ₄	Glycine-β-Naphthylamide	GLY
p-Nitrophenyl-N-Acetyl-β-D-Glucosaminide	NGLU	L-Proline-β-Naphthylamide	PRO
p-Nitrophenyl-α-D-Glucopyranoside	AGL	L-Arginine-β-Naphthylamide	ARG
o-Nitrophenyl-β-D-Glucopyranoside	BGL	L-Pyrrolidonyl-β-Naphthylamide	PYR

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REFERENCES

1. Appel, W. 1974. Peptidases, p. 949-977. In *Methods in Enzymatic Analysis*, vol. 2. Verlag Chemie International. Deerfield Beach, FL.
2. Holdeman, L.V., E.P. Cato, and W.E.C. Moore (ed). 1977. *Anaerobe Laboratory Manual*, 4th ed. Virginia Polytechnic Institute and State University, Blacksburg, VA.
3. MacFadden, J.F. 1980. *Biochemical Tests for Identification of Medical Bacteria*, 2nd ed. Williams and Wilkins, Baltimore, MD.
4. Murgier, M., C. Pelissier, A. Lazdunski, and C. Lazdunski. 1976. Existence, localization and regulation of the biosynthesis of aminoendopeptidase in Gram negative bacteria. *Eur. J. Biochem.* 65:517.
5. Porschen, R.K., and E.H. Spaulding. 1974. Phosphatase activity of anaerobic organisms. *Appl. Microbiol.* 27:744.
6. Watson, R.R. 1976. Substrate specificities of aminopeptidases: a specific method for microbial differentiation, p. 1-14. In J.R. Norris and D.W. Ribbons (ed.). *Methods in Microbiology*, vol. 9, Academic Press. New York, NY.
7. Westley, J.W., P.J. Anderson, V.A. Close, B. Halpern, and E.M. Lederberg. 1967. Aminopeptidase profiles of various bacteria. *Appl. Microbiol.* 15:822.
8. Sutter, V.L., D.M. Citron, and S.M. Finegold (eds.). 1980. *Wadsworth Anaerobic Bacteriology Manual*. The C.V. Mosby Company, St. Louis, MO.
9. Dowell, V.R. Jr., and T.W. Hawkins. 1979. *Laboratory Methods in Anaerobic Bacteriology*, CDC Laboratory Manual, Center for Disease Control, Atlanta, GA. DHEW Publication No. (CDC) 74-8272.
10. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (eds), 1995. *Manual of Clinical Microbiology*, 6th ed., American Society for Microbiology, Washington, D.C.

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.