MicroScan

Gram Positive ID 3 Procedural Manual

REF B1017-221

TABLE OF CONTENTS

English	2
Français (FR)	10
Deutsch (DE)	18
Italiano (IT)	26
Español (ES)	34
Português Portugal (PT-PT)	42
Ελληνικά (EL)	50
中文 (ZH-CN)	59
Magyar (HU)	66
Polski (PL)	74
한국어 (KO)	82
Türkçe (TR)	89
Русский (RU)	97
Hrvatski (HR)	106
Србија (SR)	114
Nederlands (NL)	122
Tiếng Việt (VI)	130
Македонски (MK)	138
REFERENCES	147

MicroScan

Gram Positive ID 3 Procedural Manual

REF B1017-221

INTENDED USE

For Use With MicroScan Dried Gram Positive ID Type 3 panels. MicroScan Dried Gram Positive ID Type 3 panels are designed for use in determining identification to the species level of rapidly growing aerobic and facultative gram-positive cocci, some fastidious aerobic gram-positive cocci and *Listeria monocytogenes*. Refer to Limitations of Procedure Section for use with fastidious streptococci.

SUMMARY AND PRINCIPLE

Modified conventional and chromogenic tests are used for the identification of *Micrococcaceae*, *Streptococcaceae*, and *Listeria monocytogenes*. Identification is based on detection of pH changes, substrate utilization, and growth in the presence of antimicrobial agents after 16–44 hours incubation at 35°C.^{1,2,3,4,5,6,7}

WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with special awareness that the inoculated panels contain potentially pathogenic organisms.
3. This material contains infectious agents and should be discarded as appropriate for biohazardous waste.
4. Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.
5. Growth and OPT well contain 0.005% hemin bovine crystalline.

STORAGE

Dried Gram Positive panels must be stored at 2–25°C.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

EVIDENCE OF DETERIORATION

Prolonged exposure to storage conditions other than those recommended may result in hydrolysis of identification substrates. Do not use beyond the expiration date. Contact your Beckman Coulter Representative or Distributor for further assistance.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens should be collected, transported, and placed on primary isolation media according to procedures recommended in the *Manual of Clinical Microbiology*.⁸

MATERIALS PROVIDED

Pos ID Type 3 (B1017-221)

Procedural Manual

MATERIALS REQUIRED, BUT NOT PROVIDED

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Inoculum Water with PLURONIC®, 25 mL (B1015-7)*

Microdilution Viewer

*PLURONIC® surfactants, a registered trademark of BASF Corporation, Parsippany, NJ USA

Mineral Oil, 60 mL (B1010-40)

Mineral Oil, 250 mL (B1010-40A) – ONLY for use with WalkAway *SI* and WalkAway *plus* instruments (and WalkAway instruments upgraded with the automated oil overlay feature)

Prompt® Inoculation System (B1026-10D)**

Quality Control Organisms (Refer to QC section of this manual)

Peptidase Reagent, 30 mL (B1012-30B) or 250 mL (B1015-30)

Reagent Dropper Kit (B1013-12A)

RENOK Rehydrator/Inoculator (B1018-14) or equivalent

Turbidity Meter

Vortex

Barcode Label Paper (B1018-129)

WalkAway Tray Lids (B1018-18)

PROCEDURE

Panel Preparation

1. Remove the panels to be used from storage. **Panels should not be used if any of the following conditions exist.**
 - A. Desiccant is not present or is broken.
 - B. Panel wells are discolored (e.g., PHO).
 - C. Integrity of the packaging is compromised (unsealed, punctured, or torn).
2. Cut open the pouch and remove the panel. If stored in the refrigerator, remove the panel immediately from the foil pouch. Allow panels to equilibrate to room temperature prior to rehydration. Panels may be stacked with a clean Cover Tray on top. All opened panels should be used within the same day or discarded.

Pos ID Type 3

Identification Substrates	Abbr.	Identification Substrates	Abbr.
Crystal Violet	CV	Urea	URE
Micrococcus Screen	MS	Mannitol	MAN
Nitrate	NIT	Lactose	LAC
Novobiocin	NOV	Trehalose	TRE
PNP-β-D-Glucuronide	PGR	Mannose	MNS
Indoxyl Phosphatase	IDX	Sodium Chloride 6.5%	NACL
Voges-Proskauer	VP	Sorbitol	SOR
Optochin	OPT	Arabinose	ARA
Phosphatase	PHO	Ribose	RBS
40% Bile Esculin	BE	Inulin	INU
L-Pyrrolidonyl-β-naphthylamide	PYR	Raffinose	RAF
Arginine	ARG	Bacitracin	BAC
PNP-β-D-Galactopyranoside	PGT	Pyruvate	PRV

Inoculum Preparation

CLSI recommends periodically checking your inoculum densities by performing colony counts. Refer to CLSI document M7-A7 section 10.3.1.3 for colony count recommendations. The expected results for *E. coli* ATCC 25922 should closely approximate 5×10^5 CFU/mL for the final test concentrations.⁹ The user should pay careful attention to inoculum preparation, especially with manual methods that are technique dependent such as the Prompt System or inoculum prepared without the aid of a photometric device.

Note: Log and Stationary phase techniques are not supported with MicroScan products.

1. Turbidity Standard Technique - Primary Inoculum Method

The turbidity standard technique is recommended for direct inoculation of all aerobic gram positive cocci.

- A. Using a sterile applicator stick, loop or swab, touch the surface of 4-5 large or 5-10 small morphologically similar, well-isolated colonies from an 18-24 hour non-inhibitory agar plate.
- B. Emulsify in 3 mL of Inoculum Water (autoclaved deionized water). The final turbidity should be equivalent to that of a 0.5 McFarland turbidity standard. An equivalent turbidity can be achieved by using a MicroScan Turbidity Meter with a range of 0.08 ± 0.02 . Cap tightly.
- C. Vortex the suspension for 2-3 seconds.

**3M, St. Paul, MN USA

- D. Pipet 0.1 mL (100 µL) of the standardized suspension into 25 mL of Inoculum Water with PLURONIC. Cap tightly. Invert 8-10 times to mix.

2. Prompt System

The Prompt System may be used for the more rapidly growing gram-positive cocci. Refer to the Prompt Inoculation procedural manual for the proper use of the Prompt System.

Note: Use with caution. Under-inoculation may occur with organisms that do not meet the size requirement. Under-inoculation could cause incorrect identification test results.

Panel Rehydration/Inoculation

Rehydration and inoculation is performed using the RENOK system with inoculators-D (B1013-4). Refer to RENOK Operator's Manual for use. If an alternate system is used, rehydrate with 115 + 10 µL of Inoculum Water (PLURONIC). A final well concentration of $3\text{--}7 \times 10^5$ CFU/mL must be achieved.

To ensure viability and purity of the organism tested, a purity plate should be prepared by streaking the inoculum to an appropriate agar plate and incubating under the appropriate conditions. If two or more colony types are present on the purity plate, reisolate the colonies and retest.

Biochemical Overlays

1. Using a dropper bottle, overlay the ARG and URE wells with at least 3 drops of mineral oil. (These wells are underlined on the panel.)
2. The media in the wells must be completely covered with mineral oil, but the oil should not overflow the wells.

Note: WalkAway *SI* instruments and WalkAway *plus* (and WalkAway instruments upgraded with the automated oil overlay feature) automatically add oil to the appropriate wells.

Incubation

1. Panels can be incubated in a WalkAway System or incubated off-line using the following steps:
 - A. To ensure even thermal distribution during incubation, stack the panels in groups of 3-5.
 - B. Place a clean Cover Tray on top of each group of panels to prevent evaporation. Cover Trays may be reused. Do not decontaminate the Cover Trays with alcohol. They may be cleaned with soap and water. Rinse well and allow to air dry.
 - C. Incubate the panels for 16-20 hours at 35°C in a non-CO₂ incubator.

Reading the Panels

Panels can be read manually using the MicroScan Microdilution Viewer and results recorded on a Manual Panel Worksheet (see Materials Required but not Provided section of this procedure manual) or on MicroScan instrumentation (touchSCAN-SR, autoSCAN-4, and WalkAway systems). Refer to the appropriate operator's manual for reading panels with MicroScan instrumentation.

1. Following 16-20 hours incubation, remove the panels from the incubator.
2. Wipe off the bottom of the panel with a lint-free tissue to remove any condensation or debris that may be present.
3. Read the panels only if the control well is clear and the growth control well is turbid. Turbidity may take the form of a white haze throughout the well, a white button in the center of the well, or a fine granular growth throughout the well. Inadequate growth is defined as slight whiteness in the well or the broth is clear.
4. If results are read manually, record results on the appropriate worksheet.
5. Reading Identification Substrates
 - A. Read all identification substrates with a white background except CV, MS, NOV, OPT, NACL, and BAC, which should be read against a black (indirectly lighted) background.
 - B. IDX and PHO results are recorded after 16-20 hours incubation.
 - C. Panels with negative PGT and fewer than 3 positive carbohydrates should be reincubated for an additional 24 hours before making a final reading and identification.
 - D. Only panels which have 3 carbohydrates or the PGT reaction positive should have reagents added at 16-20 hours. At 40-44 hours, add reagents regardless of number of positive tests.
 - E. Prior to addition of reagents, record all positive reactions.
 - F. Add reagents as follows:
 - a. Add 1 drop 40% Potassium Hydroxide (KOH) and 1 drop of 5% Alpha-Naphthol to the VP well. Wait at least 20 minutes for VP reaction to develop prior to reading.
 - b. Add 1 drop each of 0.8% Sulfanilic Acid and 0.5% N, N-Dimethyl-alpha-naphthylamine to the NIT well. Wait at least 5 minutes for the NIT reaction to develop prior to reading.
 - c. Add 2 drops of Peptidase reagent to PYR well. Wait 2 minutes prior to reading.

G. Refer to RESULTS section for aid in biochemical interpretation.

RESULTS

1. Biochemical interpretations

Well	Reagent	Positive	Negative
CV			
MS			
NOV			
OPT		Growth	No Growth
NACL			
BAC			
NIT	Add 1 drop of 0.8% Sulfanic Acid and 1 drop of 0.5% N, N-Dimethyl-alpha-naphthylamine. Wait 5 minutes for the reaction to develop.	Pink to Red	Clear (colorless) May be very pale pink
PGR	Any shade of yellow should be interpreted as positive.	Yellow	Clear (colorless)
PHO	Use the NOV well as a negative control.		
PGT			
IDX		Blue to Gray/Precipitate	Clear (colorless) or white precipitate
VP	Add 1 drop of 40% KOH and 1 drop of 5% Alpha Naphthol. Wait at least 20 minutes for reaction to develop.	Pink to Red	Clear (colorless) to Brown, May be cloudy or very pale pink
BE		Dark Brown to Black	Clear (colorless) to Light Brown
PYR	Add 2 drops Peptidase reagent. Wait 2 minutes for the reaction to develop.	Red/Orange to Red	Yellow to Orange/Red
ARG		Pink/Orange to Pink	Yellow to Orange
URE		Pink	Yellow to Orange
RAF			
LAC			
TRE			
MNS			
SOR	Any hint of orange should be considered negative.	Yellow	Orange to Red
ARA			
RBS			
INU			
MAN			
PRV			
HEM		Beta	Alpha or Gamma
LOC	For autoSCAN-4 and WalkAway Systems panel recognition only.		

Dried Gram Pos Identification - Streamlined ID QC ¹					
QC Organism	ATCC Number	Key Indicator Substrate	Expected Result	Out-of-Specification Result	Labile Substrate
<i>S. bovis</i>	49147	ARA	–	+	Yes

1. Key indicators for Quality Control are provided for streamlined QC per CLSI document M50-A.¹⁰ Before implementing the streamlined QC program, refer to the M50-A document for details and check with applicable regulatory agencies and inspecting groups to confirm acceptance.

2. Principles of Identification Reactions

Crystal Violet (CV): Growth in the presence of low concentrations of crystal violet is used to distinguish streptococci (positive) from staphylococci (mostly negative).

Micrococcus Screen (MS): Growth in the presence of low concentrations of bacitracin (0.05 µg/mL) is used to distinguish staphylococci (positive) from micrococci (negative).

Nitrate (NIT): Reduction of nitrate to nitrite is detected by the formation of a red color following the addition of 0.8% Sulfanilic Acid and 0.5% N, N-Dimethyl-alpha-naphthylamine. Streptococci are nitrate negative while most staphylococci are nitrate positive.

Novobiocin (NOV): Resistance to low concentrations (1.6 µg/mL) of Novobiocin is characteristic of some staphylococci.

Glycosidases (PGR, PGT): The ability of an organism to produce a specific glycosidase enzyme is detected by the splitting of the p-nitrophenyl-carbohydrate complex releasing p-nitrophenol which is yellow in color.

Indoxyl Phosphatase (IDX): Hydrolysis of indoxyl phosphate by the enzyme indoxyl phosphatase results in an insoluble blue compound. Most coagulase and DNase positive staphylococci are IDX positive.

Voges-Proskauer (VP): Acetylmethylcarbinol, produced from glucose, reacts with 40% Potassium Hydroxide and 5% Alpha Naphthol to form a red color.

Optochin (OPT): Susceptibility to optochin is characteristic of *Streptococcus pneumoniae*. Other streptococci and staphylococci are not inhibited by optochin.

Phosphatase (PHO): Alkaline phosphatase splits p-nitrophenyl phosphate into inorganic phosphate and p-nitrophenol which is yellow in color.

Bile Esculin (BE): Organisms capable of growing in 40% bile and hydrolyzing esculin are detected by the production of a black precipitate resulting from the reaction of the hydrolytic product esculetin with ferric citrate. Group D streptococci, some viridans streptococci, and some staphylococci are BE positive.

Pyrrolidonyl- β -naphthylamide (PYR): Organisms which produce pyrrolidonase split L-pyrrolidonyl- β -naphthylamide into L-pyrrolidone and β -naphthylamine which combines with Peptidase reagent (p-dimethyl-aminocinnamaldehyde) to produce a red color.

Arginine (ARG): Dehydrolization of arginine results in alkalinization of the medium which is detected by a yellow to red color change in the phenol red indicator.

Urea (URE): The enzyme urease splits urea forming ammonia. The resulting increase in pH is detected by the phenol red indicator.

Carbohydrates (RAF, LAC, TRE, MNS, SOR, ARA, RBS, INU, MAN): The fermentation of a specific carbohydrate results in acid formation. The consequent pH drop is detected by the phenol red indicator turning yellow.

6.5% NaCl (NACL): Tolerance to 6.5% sodium chloride is demonstrated by growth. Salt tolerance is used to differentiate enterococci from non-enterococci.

Bacitracin (BAC): Susceptibility to low concentrations of bacitracin is indicated by lack of growth and is characteristic of *Streptococcus pyogenes*.

Pyruvate (PRV): Utilization of pyruvate results in acid formation. The resultant drop in pH is detected by the phenol red indicator turning yellow.

Hemolysis (HEM): Streptolysin S and O which are produced by streptococci cause complete or partial lysis of the red blood cells in agar containing sheep blood.

3. Organism Identification

The MicroScan Dried Gram Positive Biotype Lookup Program is used for the identification of unknown test organisms. This program has been generated from MicroScan's own database. The Gram Positive Lookup Program provides for computer analysis of 27 tests for *Streptococcaceae* and 18 tests for *Micrococcaceae*. These tests are translated into a 9 or 6 digit biotype number respectively. Refer to the program for the method of recording results. The program lists the organism identification and the relative probabilities. All possibilities are listed in the order of the highest probability up to a cumulative total of 99.9%. The program is divided into 2 sections: *Micrococcaceae* and *Streptococcaceae*.

Should a biotype occur that results in a "Very Rare Biotype", consult the LabPro software (Utilities>System>Biotype Lookup) or the Biotype Lookup Service on the Beckman Coulter website, or consult your local distributor.

LIMITATION OF THE PROCEDURE

1. The CLSI Document M7-A7 recommends supplementation of Mueller Hinton Broth with Lysed Horse Blood when testing fastidious streptococci. The MicroScan Dried Gram Positive Panel Procedure differs from this recommendation. If there is inadequate growth in the growth well, the biochemical results are not valid and an alternate method should be used.
2. Anaerobic Bacteria - The Pos ID panel is not suitable for testing anaerobic gram positive cocci.
3. Additional tests may be necessary to determine the final identification when a low probability (<85%) identification is obtained (refer to the MicroScan Biotype Lookup Program for Aerobic Gram Positive Organisms or the Microbiology Manual).
4. The interpretation of test results require trained clinical personnel who should use judgment, knowledge and additional confirmatory tests where required prior to accepting the identification of an organism.
5. Biotype numbers should not be used for phenotype identification of strains isolated from various specimens from the same patient.
6. Vancomycin-resistant *E. faecium* may be distinguished from *E. gallinarum* by checking motility and pigment production. See the additional test table of the Microbiology Manual for more information.

EXPECTED VALUES

For expected values with biochemical substrates, please refer to the Percent Chart contained in the Microbiology Manual.

QUALITY CONTROL

The acceptability of the identification media should be checked by testing organisms with known reactions. The results for the recommended American Type Culture Collection (ATCC) control organisms are listed in the following table.*

*American Type Culture Collection, Manassas, VA USA.

Dried Gram Pos Identification Substrate Quality Control Table

	<i>S. aureus</i> ATCC 29213 ¹	<i>E. faecalis</i> ATCC 29212	<i>S. bovis</i> ATCC 49147	<i>M. luteus</i> ATCC 49732
CV	-	+	N/A	N/A
MS	+	+	+	-
NIT	+	-	-	N/A
NOV ²	-	N/A	-	N/A
PGR ²	-	-	-	N/A
IDX	+	N/A	N/A	-
VP	+	+	+	-
OPT ²	+	+	+	N/A
PHO	+	N/A	-	N/A
BE	-	+	+	N/A
PYR	N/A	+	-	N/A
ARG	+	+	-	N/A
PGT	+	+	+	-
URE	+ ³	-	-	N/A
MAN	+	+	+	-
LAC	+	N/A	+	-
TRE	N/A	+	+	-
MNS	N/A	+	+	-
NACL	N/A	+	-	N/A
SOR	N/A	+	-	N/A
ARA ²	N/A	-	-	N/A
RBS	N/A	+	-	N/A
INU	N/A	-	+	N/A
RAF	N/A	-	+	N/A
BAC ²	N/A	+	+	N/A
PRV	N/A	+	-	N/A

1. American Type Culture Collection, Manassas, VA USA

2. Refer to Additional Test Table for additional positive or negative reactions as needed for each biochemical.

3. May require 24 hour incubation.

N/A = Not Applicable

NOTE: Quality Control organisms are selected to provide positive/negative reactions for all identification substrates. As a result of this, the organism identifications may vary from that stated in the QC chart. The acceptability of product performance should be determined by comparison of test results, not identification.

Additional Test Table

Organism	Substrate	Expected Result
<i>S. saprophyticus</i> ATCC 49907	NOV	+
<i>S. xylosus</i> ATCC 49148	PGR	+
<i>S. pneumoniae</i> ATCC 49136	OPT	-
	BAC	-
<i>E. avium</i> ATCC 49464	ARA	+

PERFORMANCE CHARACTERISTICS

Identification

Performance claims for the MicroScan Dried Gram Positive Type 3 panel were established using the Dried Gram Positive Type 2 Panel in a multicenter study. Clinical isolates and stock strains were tested on the Dried Pos ID Type 3 panel and conventional tube methodologies to represent the type of bacterial population expected in a routine clinical laboratory. Refer to the MicroScan System Microbiology Information Manual (Section 4: Dried Conventional Pos Panels: Gram Positive Aerobic Bacteria) for a complete list of the organisms included in the database.

The Dried Pos ID Type 3 panel *Micrococcaceae* results were in agreement to both the species and the group level with 97.5% (157/161) of the isolates. Only 5.6% (9/161) of the isolates required additional tests to confirm a low probability species or group level identification.

The Dried Pos ID Type 3 panel *Streptococcaceae* results were in agreement to the species level with 91.3% (230/252) and they were in agreement to the group level with 92.1% (232/252). Only 14.3% (36/252) of the isolates required additional tests to confirm a low probability species identification and 4.4% (11/252) to confirm a low probability group level identification.

WARRANTY STATEMENT

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

Symbols Key

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1; 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1; 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1; 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1; 5.3.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1; 5.5.5
	<i>In Vitro</i> Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1; 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1; 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1; 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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REFERENCES

1. Facklam, R. R. 1973. Comparison of several laboratory media for presumptive identification of enterococci and group D streptococci. *Appl. Microbiol.* 26:138.
2. Facklam, R. R., and P. B. Smith. 1976. The Gram positive cocci. *Human Path.* 7:187.
3. Falk, D., and S. J. Guering. 1983. Differentiation of *Staphylococcus* and *Micrococcus* spp. with the taxo A bacitracin disk. *J. Clin. Microbiol.* 18:719.
4. Gross, K. C., M. P. Houghton, and L. B. Senterfit. 1975. Presumptive specification of *Streptococcus bovis* and other group D streptococci from human sources by using arginine and pyruvate tests. *J. Clin. Microbiol.* 1:54.
5. Kloos, W. E., and K. H. Schleifer. 1975. Simplified scheme for routine identification of human staphylococcus species. *J. Clin. Microbiol.* 1:82.
6. MacFaddin, J. F. 1980. *Biochemical Tests for Identification of Medical Bacteria*, 2nd ed., Williams and Wilkins, Baltimore Md.
7. Murray, P. R., E. J. Baron, J. H. Jorgensen, M.L. Landry, and M.A. Pfaller (eds), 2007. *Manual of Clinical Microbiology*, 9th ed. American Society for Microbiology, Washington D. C.
8. Facklam, R. R. 1971. Recognition of the group D streptococcal species of human origin by biochemical and physiological tests. *Appl. Microbiol.* 23:1131.
9. *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically*. 2006. Approved Standard M7-A7. Clinical and Laboratory Standards Institute, Wayne, PA.
10. *Quality Control for Commercial Microbial Identification Systems*. 2008. CLSI Document M50-A. Clinical and Laboratory Standards Institute CLSI, Wayne, PA

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.