MicroScan

Gram Negative ID2 Procedural Manual, LabPro ≥V5.0

REF B1017-27

TABLE OF CONTENTS

English	2
Français (FR)	11
Deutsch (DE)	21
Italiano (IT)	32
Español (ES)	42
Português Portugal (PT-PT)	52
Norsk (NO)	62
Ελληνικά (EL)	71
中文 (ZH-CN)	82
Lietuviškai (LT)	91
Magyar (HU)	101
Polski (PL)	111
Čeština (CZ)	122
Slovák (SK)	132
한국어 (KO)	142
Türkçe (TR)	151
Русский (RU)	161
Hrvatski (HR)	172
Български (BG)	182
中文 (ZH-TW)	193
Română (RO)	202
Србија (SR)	211
Português Brasil (PT-BR)	221
Nederlands (NL)	231
Tiếng Việt (VI)	241
Македонски (MK)	251
REFERENCES	261

MicroScan

Gram Negative ID2 Procedural Manual, LabPro ≥V5.0

REF B1017-27

INTENDED USE

MicroScan Neg ID Panel Type 2 is an *in vitro* diagnostic medical device intended for qualitative determination of identification to the species level of rapidly growing aerobic and facultatively anaerobic gram-negative bacilli from isolated colonies. The methods for incubation and results determination include full automation on WalkAway instrumentation, or non-automated offline incubation with either reading on autoSCAN-4 instrumentation or manual reading with the option to enter results into the LabPro data management system.

INTENDED USER

This product is intended for laboratory professional use.

SUMMARY AND PRINCIPLE

After panels are inoculated and rehydrated with a standardized suspension of organism, they are incubated for a minimum of 16 hours at 35°C. Modified conventional and chromogenic tests are used for the identification of fermentative and non-fermentative gram-negative bacilli. Identification is based on detection of pH changes, substrate utilization and growth in the presence of antimicrobial agents after 16-42 hours incubation at 35°C.¹

Neg ID Panel Type 2			
Substrates	Abbr.	Substrates	Abbr.
Glucose	GLU	Voges-Proskauer	VP
Sucrose	SUC	Citrate	CIT
Sorbitol	SOR	Malonate	MAL
Raffinose	RAF	o-Nitrophenyl-β-D-Galactopyranoside	ONPG
Rhamnose	RHA	Tartrate	TAR
Arabinose	ARA	Acetamide	ACE
Inositol	INO	Cetrimide	CET
Adonitol	ADO	OF Glucose	OF/G
Melibiose	MEL	OF Base	OF/B
Urea	URE	Decarboxylase Base	DCB
Hydrogen Sulfide	H ₂ S	Nitrate	NIT
Indole	IND	Penicillin G 4 µg/mL ¹	P ₄
Lysine	LYS	Kanamycin 4 µg/mL ¹	K ₄
Arginine	ARG	Colistin 4 µg/mL ¹	Cl ₄
Ornithine	ORN	Nitrofurantoin 64 µg/mL ¹	Fd ₆₄
Tryptophan Deaminase	TDA	Cephalothin 8 µg/mL ¹	Cf ₈
Esculin	ESC	Tobramycin 4 µg/mL ¹	To ₄

1. Antimicrobics used for identification.

CLINICAL RELEVANCE

Identification test results are used in conjunction with other laboratory results and clinical indicators to refine therapy in symptomatic patients.

WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with special awareness that the inoculated panels contain potentially pathogenic organisms.
3. This material contains infectious agents and should be discarded as appropriate for biohazardous waste.
4. Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.
5. Single Use only. All MicroScan panels and consumables should not be reused unless otherwise noted (e.g. cover trays).
6. This product contains material(s) of animal origin from dehydrated microbiology media. Observe general safety guidelines for protection when handling this product.

STORAGE

Dried Gram Negative panels must be stored at 2 - 25°C.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

EVIDENCE OF DETERIORATION

Prolonged exposure to storage conditions other than those recommended may result in hydrolysis of identification substrates. Do not use beyond the expiration date. Contact your Beckman Coulter Representative or Distributor for further assistance.

SPECIMEN COLLECTION AND PREPARATION

Appropriate specimens should be collected, transported and placed on primary isolation media according to the procedures recommended in the *Manual of Clinical Microbiology*.²

MATERIALS PROVIDED

Neg ID Type 2 (B1017-27)

Procedural and Quality Control Manual

MATERIALS REQUIRED, BUT NOT PROVIDED

5% Alpha Naphthol, 30 mL (B1010-42A)

10% Ferric Chloride, 30 mL (B1010-48A) or 250 mL (B1015-48)

0.5% N, N-Dimethyl-alpha-naphthylamine, 30 mL (B1010-45A) or 250 mL (B1015-45)

40% Potassium Hydroxide, 30 mL (B1010-43A) or 250 mL (B1015-43)

0.8% Sulfanilic Acid, 30 mL (B1010-44A) or 250 mL (B1015-44)

0.5 McFarland Turbidity Standard

100 µL pipettor

Barcode Label Paper (B1018-129)

Cover Trays (B1010-56B)

General Laboratory Equipment

Kovac's Reagent, 30 mL (B1010-41A) or 250 mL (B1015-41)

Manual Panel Worksheets (used for recording manually read results)

Microdilution viewer

Mineral Oil, 60 mL (B1010-40)

Mineral Oil, 250 mL (B1010-40A)

Oxidase reagent

Quality Control Organisms (Refer to QC section of this manual)

Quality Control Report Forms

Reagent Dropper Kit (B1013-12A)

RENOK - Rehydrator/Inoculator (B1018-14) or equivalent

Turbidity Meter

Vortex

SINGLE USE MATERIALS

0.85% Autoclaved Saline, 3 mL

Disposable sterile pipette tips

Inoculator-D Set (B1013-4)

Inoculum Water, 3 mL (B1015-2)

Inoculum Water with PLURONIC®, 25 mL (B1015-7)*

Prompt® Inoculation System (B1026-10D)**

Seal Strips (B1010-51)

WalkAway Tray Lids (B1018-18)

*PLURONIC® Surfactants. Registered trademark of BASF Corporation, Parsippany, NJ USA.

**3M, St. Paul, MN USA

PROCEDURE

Panel Preparation

1. Remove the panels to be used from storage. Do not use if the integrity of the packaging is compromised (unsealed, punctured, or torn).
2. Cut open the pouch and remove the panel. If stored in the refrigerator, remove the panel immediately from the foil pouch.
3. **Panels should not be used if any of the following conditions exist:**
 - A. Desiccant is not present or is broken.
 - B. Panel wells are discolored (e.g., DCB, various antimicrobics).
4. Allow panels to equilibrate to room temperature prior to rehydration. Panels may be stacked with a clean Cover Tray on top. All opened panels should be used within the same day or discarded.

Inoculum Preparation

CLSI recommends periodically checking your inoculum densities by doing colony counts. Refer to CLSI document M07, 11th ed. for colony count recommendations. The expected results for *E. coli* ATCC 25922 should closely approximate 5×10^5 CFU/mL for the final test concentrations.³ The user should pay careful attention to inoculum preparation especially with manual methods that are technique dependent such as the Prompt system or inoculum prepared without the aid of a photometric device.

1. Turbidity Standard Technique - *Primary Inoculum Method*

The turbidity standard technique is recommended for direct inoculation of all aerobic gram-negative bacilli.

- A. Using a sterile applicator stick, loop or swab, touch the surface of 4-5 large or 5-10 small morphologically similar, well-isolated colonies from an 18-24 hour non-inhibitory agar plate.
- B. Emulsify in 3 mL of Inoculum Water (autoclaved deionized water).
- C. Cap tightly and vortex the suspension for 2-3 seconds.
- D. The final turbidity should be equivalent to that of a 0.5 McFarland Turbidity Standard. An equivalent turbidity can be achieved by using a MicroScan Turbidity Meter with a range of 0.08 ± 0.02 .
- E. Pipet 0.1 mL (100 μ L) of the standardized suspension into 25 mL of Inoculum Water with PLURONIC. Cap tightly. Invert 8-10 times to mix.

2. Prompt System

The Prompt System may be used to inoculate gram-negative bacilli. Refer to the Prompt Inoculation procedural manual for the proper use of the Prompt System.

NOTE: Use with caution. Under-inoculation may occur with organisms that do not meet the size requirement. Under-inoculation could cause incorrect susceptibility and identification test results.

Oxidase Test

Perform an oxidase test prior to inoculating panels. Record results in the appropriate space on the worksheet or as required by the instrumentation. Oxidase results support separation of taxa and early reporting of non-fermenters (e.g., *Acinetobacter baumannii* complex). For WalkAway Systems, oxidase results are recommended to be entered prior to printing barcodes. Some non-fermenters may still hold to 42 hours, depending on growth. The recommended oxidase reagent is tetramethyl-p-phenylene-diamine-dihydrochloride.

Panel Rehydration/Inoculation

Rehydration and inoculation are performed using the RENOK system with Inoculators-D (B1013-4). Refer to the RENOK Operator's Manual for use. If an alternate system are used, rehydrate with $115 \mu\text{L} \pm 10 \mu\text{L}$ Inoculum Water with PLURONIC. A final well concentration of $3\text{--}7 \times 10^5$ CFU/mL should be achieved. To ensure viability and purity of the organism tested, a purity plate should be prepared by streaking the inoculum to an appropriate agar plate and incubating under the appropriate conditions. If two or more colony types are present on the purity plate, reisolate the colonies and retest.

Biochemical Overlays

1. Use a dropper bottle to overlay the GLU, URE, H₂S, LYS, ARG, ORN, and DCB wells with 3 drops of mineral oil. These wells are underlined on the panel.
2. The media in the wells must be completely covered with mineral oil, but the oil should not overflow the wells.

NOTE: WalkAway Systems (and WalkAway Systems upgraded with the automated oil overlay feature) automatically add oil to the appropriate wells.

Seal Strip

For oxidase positive gram-negative organisms only, place a seal strip over the CIT, MAL, ONPG, TAR, ACE, CET, OF/G, OF/B, and DCB wells. The quarter-inch locator hole in the tape should be aligned over the DCB well. For WalkAway Systems, the lid is used in lieu of a seal strip.

Incubation

1. Panels can be incubated in a WalkAway System or incubated off-line using the following steps:
 - A. To ensure even thermal distribution during incubation, stack the panels in groups of 3-5.
 - B. Place a clean Cover Tray on top of each group of panels to prevent evaporation. Cover Trays may be reused. Do not decontaminate Cover Trays with alcohol. They may be cleaned with soap and water. Rinse well and allow to air dry.
 - C. Incubate the panels for a minimum of 16 hours at 35°C in a non-CO₂ incubator.

Reading the Panels

Panels can be read manually and results recorded on a Manual Panel Worksheet or on MicroScan instrumentation (autoSCAN-4, and WalkAway Systems). Refer to the Operator's Guide for reading panels with MicroScan instrumentation.

1. Following 16-20 hours of incubation, remove the panels from the incubator.
2. Wipe off the bottom of the panel with a lint-free tissue to remove any condensation or debris that may be present.
3. Read the panels only if the control well is clear and the growth well is turbid. Do not read the panel if the control well is turbid, or if there is no growth in the growth well. Turbidity may take the form of a white haze throughout the well, a white button in the center of the well, or a fine granular growth throughout the well. Inadequate or no growth is defined as a slight whiteness in the well or the broth is clear.
4. If results are read manually, record the results on the appropriate worksheet.

5. Reading Identification Substrates

- A. Read all identification substrates with a white background except for CET and the antimicrobials used for identification (Cl₄, Cf₈, P₄, K₄, Fd₆₄, To₄) which should be read against a black (indirectly lighted) background.
- B. The following substrates are always read at 16-24 hours: GLU, SUC, SOR, RAF, RHA, ARA, INO, ADO, MEL, P₄, K₄, Cl₄, Fd₆₄, Cf₈, To₄.
- C. If any of the following is true:
 - a. GLU or SUC or SOR positive, or
 - b. ARG positive and OF/G positive, or
 - c. ARG positive and CET positive, or
 - d. LYS positive, or
 - e. ORN positive, or
 - f. OXI negative and OF/G positive and Fd₆₄ positive and MAL positive.

Then add reagents and read the remaining substrates at 16-24 hours of incubation. If none of the conditions above is true, reincubate panels for an additional 24 hours before adding reagents and reading the remaining substrates.

- D. If, after reincubation, the organism still appears nonreactive, check to ensure that the organism has grown in the phenol red carbohydrates (particularly if there is poor growth or no growth in the Mueller-Hinton growth control). If there is no growth in the phenol red carbohydrates, one should suspect a fastidious organism such as a halophilic *Vibrio* species. A Gram stain should be repeated to confirm that it is gram-negative.
 - a. Should a clinically significant disease due to halophilic *Vibrio* species be suspected, retest by emulsifying several colonies into 3.0 mL of sterile 0.85% saline and also perform additional testing including growth in salt, as recommended by the software. Final turbidity should be equivalent to a 0.5 McFarland Turbidity Standard. An equivalent turbidity can be achieved by using a MicroScan Turbidity Meter with a range of 0.08 ± 0.02. Rehydrate the panel with 25 mL of uninoculated Inoculum Water with PLURONIC using a RENOK. Using a sterile transfer pipette, add one drop (45-50 µL) of the saline suspension to each identification well including the growth, Cl₄, Cf₈, To₄, P₄, K₄, and Fd₆₄ wells. All other panel processing steps should be followed.
- E. Prior to reading, add the reagents as follows:
 - a. Add 1 drop of 40% Potassium Hydroxide (KOH) and 1 drop of 5% Alpha Naphthol to the VP well. Wait at least 20 minutes for the VP reaction to develop.
 - b. Add 1 drop of 10% Ferric Chloride to the TDA well. The color will develop immediately.
 - c. Add 3 drops of MicroScan Kovac's reagent to the IND well. The color will develop immediately.

NOTE: The MicroScan Kovac's reagent is a special formula designed for the use with MicroScan Instrumentation. Only one drop need be used if panels are read manually.
 - d. Add 1 drop of 0.8% Sulfanilic Acid then 1 drop of 0.5% N, N-Dimethyl-alpha-naphthylamine to the NIT well. Wait at least 5 minutes prior to reading for the NIT reaction to develop.
- F. Refer to the RESULTS section for aid in biochemical interpretation.

RESULTS

1. Biochemical Interpretations

Well	Reagent	Positive	Negative
GLU		Strong Yellow only	Orange to Red
		Some non-fermenters may produce a golden color which should be considered negative.	
SUC		Yellow to Yellow/Orange	Orange to Red
SOR			
RAF			
RHA			
ARA			
INO			
ADO			
MEL			
URE		Magenta to Pink	Yellow, Orange or Light Pink
H ₂ S		Black Precipitate or Black Button	No Blackening
IND	Add 3 drops (or one drop if panel is read visually) of MicroScan Kovac's Reagent. Color develops immediately.	Pink to Red	Pale Yellow to Orange
		<i>Providencia</i> sp. may produce a pale pink ring which should be considered positive.	
Compare LYS, ARG, and ORN with DCB. These wells must be overlaid with oil for results to be valid. Purple is considered positive when comparing to a gray DCB well. Gray is considered positive when comparing to a yellow DCB well. For non-fermenters, a positive test must be significantly more purple than the base control. If the base control is purple, the basal medium has been alkalized and the LYS, ARG and ORN should all be recorded as negative. <i>Achromobacter</i> species and <i>Ochrobactrum anthropi</i> tend to do this.			
LYS		Purple to Gray	Fermenters: Yellow
ARG			Non-fermenters:
ORN		Purple	Colorless to Gray
TDA	Add 1 drop of 10% Ferric Chloride. Color develops immediately.	Any Shade of Brown	Yellow to Orange
ESC		Light Brown to Black	Beige to Colorless
VP	Add 1 drop of 40% KOH and 1 drop of 5% Alpha Naphthol. Wait at least 20 minutes for the reaction to develop.	Red	Colorless
		Some non-fermenters may produce a pale pink color after 18 hours incubation which should be considered negative.	
ONPG		Yellow	Colorless
		Compare any ONPG well which looks clear to the cetrimide well. If the ONPG shows any yellow as compared to the cetrimide well, record as positive.	
CIT		Blue to Blue-Green	Green to Yellow
MAL		Any shade of blue is considered positive	
TAR			
ACE			
CET		Growth	No Growth
OF/G		Yellow	If OF base is Green: Green to Blue
		Yellow to Green	If OF base is Blue to Blue Green: Blue
		Compare OF/G well with OF base control.	
NIT	Add 1 drop of 0.8% Sulfanilic Acid then 1 drop of 0.5% N, N-Dimethyl-alpha-naphthylamine. Wait at least 5 minutes for the reaction to develop.	Red	Colorless to Pale Pink

Well	Reagent	Positive	Negative
P ₄		Growth (Resistant)	No Growth
K ₄			(Susceptible)
Cl ₄			
Fd ₆₄			
To ₄			
Cf ₈			
LOC	For autoSCAN-4 and WalkAway System panel recognition only.		

2. Principles of Identification Reactions

Carbohydrate Fermentation (GLU, SUC, SOR, RAF, RHA, ARA, INO, ADO, MEL): The fermentation of a specific carbohydrate results in acid formation and a pH drop is detected by a phenol red indicator.

Urea (URE): The enzyme urease splits urea forming ammonia. The resulting increase in pH is detected by the phenol red indicator.

Hydrogen Sulfide (H₂S): Hydrogen sulfide gas is produced from sodium thiosulfate and reacts with ferric ions in the medium to produce a black precipitate.

Indole (IND): The metabolism of tryptophan results in the formation of indole, which is detected by the addition of Kovac's reagent. If indole is present, a red color develops.

Lysine, Arginine, Ornithine (LYS, ARG, ORN): Decarboxylation of these amino acids results in the formation of basic amines which are detected by the bromocresol purple indicator.

Tryptophan Deaminase (TDA): Bacteria capable of deaminating tryptophan produce indolepyruvic acid, which reacts with ferric ammonium citrate in the medium after the addition of ferric chloride to produce a brown color.

Esculin Hydrolysis (ESC): Hydrolysis of esculin is detected by ferric ammonium citrate in the medium, which reacts with the hydrolytic products to form a black precipitate.

Voges-Proskauer (VP): Acetoin is produced from sodium pyruvate and is indicated by the formation of a red color after the addition of 40% Potassium Hydroxide and 5% Alpha Naphthol.

Galactosidase (ONPG): β-galactosidase hydrolyzes ortho-nitrophenyl-β-D-galactopyranoside, which liberates the yellow-colored ortho-nitrophenol.

Citrate, Malonate, Acetamide, Tartrate (CIT, MAL, ACE, TAR): The utilization of these substrates as the sole source of carbon for metabolism results in a rise in pH that is detected by bromthymol blue indicator.

Oxidation-Fermentation (OF/G): The oxidation of glucose results in acid formation and a pH drop is detected by the bromthymol blue indicator. The OF/G is compared to OF/B (base) to determine if any acid has been produced.

Nitrate (NIT): The ability of an organism to reduce nitrate to nitrite is detected by the addition of N, N-Dimethyl-α-naphthylamine and Sulfanilic Acid, which produce a red color in the presence of nitrite.

Cetrimide (CET): Tolerance to cetrimide is demonstrated by growth in Mueller-Hinton broth supplemented with cetrimide.

Penicillin, Kanamycin, Colistin, Cephalothin, Nitrofurantoin, Tobramycin (P₄, K₄, Cl₄, Cf₈, Fd₆₄, To₄): Resistance to specific concentrations of these antimicrobial agents is demonstrated by growth.

3. Organism Identification

The Biotype Lookup Program in the LabPro Information Manager and on the Beckman Coulter website is used for the identification of unknown test organisms. The program lists the organism identification and the relative probabilities, in the order of the highest probability, up to a cumulative total of 99.9%.

Should a biotype number occur that results in a "Very Rare Biotype", consult the LabPro Software (Utilities>System>Biotype Lookup), the Biotype Lookup Program on the Beckman Coulter website, or your Beckman Coulter Representative or Distributor.

LIMITATION OF THE PROCEDURE

1. *Haemophilus* and other fastidious gram-negative bacteria will not grow in Mueller-Hinton broth that is not supplemented with growth enrichment factors.
2. The procedures described herein are suitable only for facultative anaerobic and aerobic gram-negative bacilli.
3. Additional tests may be necessary to determine the final identification when a low probability (<85%) identification is obtained (refer to the MicroScan Biotype Lookup Program or the Microbiology Manual).
4. The interpretation of test results requires trained clinical personnel who should use judgment, knowledge and additional confirmatory tests where required prior to accepting the identification of an organism.
5. Biotype numbers should not be used for phenotype identification of strains isolated from various specimens from the same patient.

6. Public health officials should be notified about all isolates presumptively identified as *B. anthracis*, *Y. pestis*, *B. mallei*, *B. pseudomallei* or *F. tularensis*. Confirmation of isolates of these bacteria may require specialized testing only available in reference or public health laboratories.

QUALITY CONTROL

The acceptability of the identification media should be checked by testing organisms with known reactions. The results for the recommended American Type Culture Collection (ATCC) control organisms are listed in the following table.*

	<i>E. coli</i> ATCC 25922 ¹	<i>P. aeruginosa</i> ATCC 27853	<i>K. oxytoca</i> ATCC 49131	<i>P. vulgaris</i> ATCC 49132
GLU	+	-	+	+
SUC	-	-	+	+
SOR	V	-	+	-
RAF	-	-	+	-
RHA	+	-	+	-
ARA	+	-	+	-
INO	-	-	+	-
ADO	-	-	+	-
MEL	V	-	+	-
URE ²	-	V	V	V
H ₂ S	-	-	-	+
IND	+	-	+	+
LYS	+	-	+	-
ARG	-	+	-	-
ORN	+	-	-	-
TDA	-	-	-	+
ESC	-	-	+	V
VP	-	- ³	+	V
CIT	-	+	+	V
MAL	-	+	+	-
ONPG	+	-	+	-
TAR ²	- ⁴	V	V	V
ACE	-	+ ⁵	-	V
CET	-	+	-	-
OF/G ²	+	V	+	+
OF/B	Blue, Blue-Green, Green	Blue, Blue-Green	Blue, Blue-Green, Green	Blue, Blue-Green, Green
P ₄ ²	+	+	+	V
NIT ²	+	V	+	+
DCB	Gray to Yellow	Gray to Yellow-Gray	Yellow	N/A
Cl ₄	-	-	-	+
Fd ₆₄	-	+	-	N/A
Cf ₆ ²	V	+	V	V
K ₄ ²	V	+	V	N/A
To ₄ ²	-	-	-	N/A

1. American Type Culture Collection, Manassas, VA USA

2. Refer to the Additional Test Table for additional positive or negative reactions as needed for each biochemical.

3. *P. aeruginosa* ATCC 27853 may produce a pale pink color in VP after 18 hours of incubation which should be interpreted as negative.

4. A weak positive reaction may be obtained with *E. coli* ATCC 25922 if panels are incubated longer than 18 hours.

5. May require 48 hours incubation.

N/A = Not Applicable

NOTE: Quality Control organisms are selected to provide positive/negative reactions for all identification substrates. As a result, the organism identifications may vary from that stated in the QC chart. The acceptability of product performance should be determined by comparison of test results, not identification.

Additional Test Table

Organism	Substrate	Expected Result
<i>P. stuartii</i> ATCC 49809	URE	+
<i>P. putida</i> ATCC 49128	TAR	+
	NIT	-
	K ₄	-
<i>S. halotis</i> ATCC 49138	OF/G	-
<i>E. faecalis</i> ATCC 29212	P ₄	-

*American Type Culture Collection, Manassas, VA USA

Organism	Substrate	Expected Result
<i>S. aureus</i> ATCC 29213	Cf ₈	-
<i>R. insidiosa</i> ATCC 49129	To ₄	+

Dried Gram Neg Identification - Streamlined ID QC ¹					
QC Organism	ATCC Number	Key Indicator Substrate	Expected Result	Out-of- Specification Result	Labile Substrate
<i>E. coli</i>	25922	ONPG	+	-	Yes
<i>E. coli</i>	25922	ARA	+	-	Yes

1. Key indicators for Quality Control are provided for streamlined QC per CLSI document M50-A.⁴ Before implementing the streamlined QC program, refer to the M50-A document for details and check with applicable regulatory agencies and inspecting groups to confirm acceptance.

PERFORMANCE CHARACTERISTICS

Performance claims for the MicroScan Dried Gram Negative Identification panel were established in a multicenter study. Clinical isolates and stock strains were tested on the Dried Neg ID Type 2 panel and conventional tube methodologies to represent the type of bacterial population expected in a routine clinical laboratory. Refer to the Dried (Conventional) Neg Panels: Gram-Negative Aerobic Bacilli section of the Microbiology Information Manual for a complete list of organisms included in the database.

Percent Agreement to Reference Method

# Tested	# Correct Identification	% Correct Identification
609	593	97.4

REPRODUCIBILITY

Reproducibility is established by testing multiple replicates under different conditions (e.g., panel lots, days, instruments, users). This may include QC strains, clinical isolates or both. Overall agreement was ≥95%.

WARRANTY STATEMENT

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

NOTICE TO USER

Any serious incident that has occurred in relation to MicroScan panels shall be reported to Beckman Coulter and the competent authority of the Member State in which the user and/or the patient is established.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	May 2022	Initial Release
AB	December 2022	Updated: Reproducibility

Symbols Key

The following symbols may be used in the labeling of Beckman Coulter Microbiology products.

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1, clause 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1, clause 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1, clause 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1, clause 5.1.3
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1, clause 5.5.5
	<i>In Vitro</i> Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1, clause 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1, clause 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1, clause 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
	Info For USA Only	Indicates information that is required only for products intended to be sold in the US.	n/a
	Prescription Use Only	Recognized by the US FDA as an alternate to "Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner"	21CFR 801.109(b)(1)
	Biological Risks	Indicates that there are potential biological risks associated with the medical device. Note: This symbol is not to be confused with the "biohazard" sign intended to be used in the workplace.	ISO 15223-1, clause 5.4.1
	Danger	n/a	n/a
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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REFERENCES

1. MacFaddin, J. F. 2000. Biochemical Tests for Identification of Medical Bacteria, 3rd ed., Williams and Wilkens, Baltimore, MD.
2. Jorgensen, J.H., M.A. Pfaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter, D.W. Warnock (eds), 2015. Manual of Clinical Microbiology, 11th ed. American Society for Microbiology, Washington D. C.
3. CLSI. *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically* 11th ed. CLSI standard M07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
4. Quality Control for Commercial Microbial Identification Systems. 2008. CLSI Document M50-A. Clinical and Laboratory Standards Institute, Wayne, PA

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.