MicroScan

Rapid Yeast ID Procedural Manual

REF B1017-70

TABLE OF CONTENTS

English	2
Français (FR)	8
Deutsch (DE)	15
Italiano (IT)	22
Español (ES)	29
Português Portugal (PT-PT)	36
Norsk (NO)	43
Ελληνικά (EL)	49
中文 (ZH-CN)	56
Lietuviškai (LT)	62
Magyar (HU)	68
Polski (PL)	75
한국어 (KO)	82
Türkçe (TR)	88
Русский (RU)	94
Български (BG)	101
Română (RO)	108
Србија (SR)	114
Português Brasil (PT-BR)	120
Nederlands (NL)	127
Tiếng Việt (VI)	134
REFERENCES	140

MicroScan

Rapid Yeast ID Procedural Manual

REF B1017-70

INTENDED USE

The Rapid Yeast Identification Panel is designed for rapid identification of yeast and yeast-like species isolated from clinical specimens.

SUMMARY AND PRINCIPLE

Chromogenic and modified conventional tests are used for the identification of yeast species isolated from clinical specimens. The Rapid Yeast Identification Panel is a 96 well microdilution panel that utilizes 27 dehydrated substrates (Refer to Table 1-1). A heavy suspension of the yeast in autoclaved water is used to rehydrate the substrates. An identification is obtained after incubation of the panel at 35-37°C for 4 hours.

WARNING AND PRECAUTIONS

1. Rapid Yeast Identification Panels are for *in vitro* diagnostic use only.
2. Do not allow Peptidase or Sodium Hydroxide reagent or any of the identification substrates to come in contact with eyes, skin, or clothing.
3. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with awareness that inoculated panels contain potentially pathogenic organisms.
4. All materials should be autoclaved prior to disposal.

STORAGE

Store panels and Peptidase reagent at 2-8°C.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION AND PREPARATION

Refer to the American Society for Microbiology Manual of Clinical Microbiology¹ or other mycology reference manuals for recommended procedures for the transport and isolation of yeasts. Rapid chromogenic tests are based upon the presence of preformed enzymes in the unknown organism.^{2,3,4} Thus, the plating medium on which the organism is grown will influence the test results because of the induction of different enzymes by different media. The Rapid Yeast Identification Panel is based upon the use of Sabouraud Dextrose Agar (either standard or Emmons' modification). The presence of antibiotics in the growth medium is not a critical variable, however media containing blood, Yeast Extract-Malt Extract (YM) Agar and Inhibitory Mold Agar (IMA) should not be used.

MATERIALS PROVIDED

Rapid Yeast Identification Panels (B1017-70)

Worksheet

Procedural Manual

Quality Control Chart

MATERIALS REQUIRED, BUT NOT PROVIDED

Inoculum Water (B1015-2)

0.05N Sodium Hydroxide (B1015-3)

Peptidase Reagent (B1012-30B)

Yeast Turbidity standard (B1015-18)

Cover Trays (B1010-56B)

MicroScan Yeast Biotype Codebook (B1014-1A)

Inoculating loops

Sabouraud Dextrose Agar plates

Vortex mixer

50 µL pipettor with tips

Incubator (35-37°C), non-CO₂

Quality control organisms:

Candida albicans ATCC 66027*

Candida kefyr ATCC 66028

Candida tropicalis ATCC 66029

Cryptococcus albidus ATCC 66030

Cryptococcus neoformans ATCC 66031

Candida glabrata ATCC 66032

Cryptococcus uniguttulatus ATCC 66033

PROCEDURE

Preparation of Panels

1. Remove the panels to be used from storage. **Panels should not be used if the desiccant is not present or is broken or if the integrity of the packaging is compromised (unsealed, punctured, or torn).**
2. Cut open the pouch and remove the panel. Remove the panel immediately from the foil pouch. Allow panels to equilibrate to room temperature prior to rehydration. Panels may be stacked with clean Cover Tray on top.

Inoculum Preparation

1. Inoculate the Rapid Yeast Identification Panel with an actively growing organism. Colonies from a primary isolation Sabouraud Dextrose Agar plate may be used if sufficient growth is present. If necessary, subculture the yeast isolate onto a Sabouraud Dextrose Agar plate and incubate until sufficient growth occurs. Incubation for 48 hours at 30°C or 48 to 72 hours at 25°C (room temperature) is recommended.

NOTE: Cultures incubated for more than 48 hours at 30°C should be used only if additional incubation is required to obtain adequate growth.

2. Using a sterile swab or inoculating loop, remove growth from the agar plate and emulsify in 3 mL of Inoculum Water (autoclaved deionized water in a 13 x 84 mm or 13 x 100 mm tube). Each suspension must be compared and be equivalent to the MicroScan Yeast Turbidity Standard. The comparison to the test inoculum suspension and the Yeast Turbidity Standard should be done using an adequate light source and by comparing the tubes to a white card with contrasting black lines.

NOTE: Take care to insure that no agar medium is removed with the organism.

3. Mix or vortex the suspension until it is homogeneous. Some yeast strains do not readily disperse in water. If vortexing does not result in a homogeneous suspension, let the suspension sit for 10-15 minutes and revortex. If the suspension still has clumps, allow the clumps to settle to the bottom of the tube and use the supernatant to inoculate the yeast panel. The turbidity of the supernatant must match the MicroScan Yeast Turbidity Standard. If it does not, transfer more growth to the Inoculum Water and repeat the above procedure.

NOTE: DO NOT pipette clumps of cells into the yeast panel.

Panel Inoculation

1. Label the panel with the specimen number.
2. Add 50 µL of yeast suspension to each substrate-containing well plus the control wells BNAC and NPC. Use of a Pasteur pipette for inoculation is not recommended.

Table 1-1, Substrates

Substrates	Abbr.	Substrates	Abbr.
Hydroxyproline-β-Naphthylamide	HPR	L-Histidine-β-Naphthylamide	HIS
L-Isoleucine-β-Naphthylamide	ILE	Sucrose	SUC1
L-Proline-β-Naphthylamide	PRO	Sucrose	SUC2 ¹
L-Tyrosine-β-Naphthylamide	TYR	Trehalose	TRE
Glycine β-Naphthylamide	GLY	p-Nitrophenyl-α-D-Glucopyranoside	AGL1
Glycylglycine-β-Naphthylamide	GGLY	p-Nitrophenyl-α-D-Glucopyranoside	AGL2 ¹
Glycyl-L-Arginine-4-Methoxy-β-Naphthylamide	GLAR	p-Nitrophenyl-β-D-Glucopyranoside	BGL
Glycyl-L-Proline-4-Methoxy-β-Naphthylamide	GLPR	o-Nitrophenyl-β-D-Galactopyranoside	BGAL
L-Arginyl-L-Arginine-β-Naphthylamide	AARG	p-Nitrophenyl-β-D-Fucopyranoside	BDF
L-Lysyl-L-Alanine-4-Methoxy-β-Naphthylamide	LYAL	p-Nitrophenyl-α-D-Galactopyranoside	AGAL
L-Alanine-4-Methoxy-β-Naphthylamide	ALA	p-Nitrophenyl-N-Acetyl-β-D-Glucosamine	NAG
L-Seryl-L-Tyrosine-β-Naphthylamide	STY	p-Nitrophenyl-β-D-Cellobiose	CELL
Urea	URE	p-Nitrophenyl-N-Acetyl-β-D-Galactosaminide	NGAL
3-Indoxyl Phosphate	IDX		

1. Different formulations of the same Substrates are used.

*American Type Culture Collection, Manassas, VA USA

- After the panel has been inoculated with the 50 μ L inoculum, lightly tap each side of the panel to help ensure mixing with the substrate.

NOTE: If the panel is to be read on the WalkAway or autoSCAN-4 the LOC well must also be inoculated.

Incubation

- Stack panels in groups of 3 to 5 with a cover tray over the top panel to prevent evaporation.
- Incubate the inoculated panels for 4 hours at 35-37°C in a non-CO₂ incubator.

Reading the Panels

- Remove the Peptidase reagent from the refrigerator at least 30 minutes prior to reading the panels. The Peptidase reagent may precipitate at 4°C, but it will redissolve at room temperature. Do not expose the Peptidase reagent to elevated temperatures to accelerate the redissolving of the crystalline precipitate because high temperatures will cause the reagent to degrade.
- Read the panels using reflected light from a white background.
- After 4 hours of incubation, add the Peptidase and Sodium Hydroxide (NaOH) reagents to the appropriate wells.

- Add **1 drop** of Peptidase reagent to the BNAC, HPR, ILE, PRO, TYR, GLY, GGLY, GLAR, GLPR, AARG, LYAL, ALA, STY and HIS wells. Allow the color to develop for 30 seconds before recording the results. Results must be recorded within 3 minutes of reagent addition. The presence of ANY shade of pink or red in the solution throughout the entire well or any portion of the well is interpreted as positive. Colors should be compared to the BNAC well. Some wells may appear a darker orange than the BNAC well. These should be recorded as negative. Reactions are positive only if there is pink or red present in the solution. The color in these wells will darken with time. A negative yellow color may become a darker orange color. This should be interpreted as a negative reaction.

NOTE: Following addition of Peptidase to the BNAC well, it should be yellow. If there is any pink or red present, the Peptidase reagent should not be used.

- Add **1 drop** of 0.05N NaOH to AGL2, BDF, AGAL, NPC, NAG, CELL and NGAL wells. Wait at least 5 seconds, but no longer than 5 minutes, before recording results. Any shade of yellow should be recorded as a positive. Compare the substrate-containing wells to the NPC well.

NOTE: These wells may appear darker than the NPC well. This should be interpreted as negative, but a yellow color must be present to be considered positive.

- Two control wells are provided to assist in interpreting the reactions.
 - BNAC = β -naphthylamide control. After addition of Peptidase reagent, this well will be yellow in color.
 - NPC = nitrophenyl control. This well should be colorless.
- Refer to the RESULTS section for aid in biochemical interpretation.

NOTE: The WalkAway reads wells that do not require reagents first to prevent altered results due to reagent fumes that may be trapped under the tray lid, therefore, biochemical reactions cannot be verified manually.

RESULTS

Substrate Interpretations

Well	Reagent	Positive	Negative
HPR ILE PRO TYR GLY GGLY GLAR GLPR AARG LYAL ALA STY HIS	Add 1 drop of Peptidase reagent. Allow color reaction to develop for at least 30 seconds, but no longer than 3 min.	Any shade of Pink, Red to Magenta in the solution ¹	Yellow to Orange
SUC1 SUC2 TRE		Any shade of Brown or Yellow	Purple
AGL1 BGL BGAL	Compare to the NPC control well.	Any shade of yellow ²	Clear
URE		Pink to Magenta	Clear/Beige to Yellow

Well	Reagent	Positive	Negative
IDX		Any shade of Blue	Clear
AGL2 BDF AGAL NAG CELL NGAL	Add 1 drop of 0.05 N NaOH. Wait at least 5 seconds, but no longer than 5 minutes before recording result. Compare to the NPC control well.	Any shade of Yellow ²	Clear

1. A positive color may not be present throughout the entire well.

2. A pink color may be present with some yeasts. This should be considered a positive reaction.

PRINCIPLES OF IDENTIFICATION REACTIONS

Amino acid β -naphthylamide substrates (HPR, ILE, PRO, TYR, GLY, GGLY, GLAR, GLPR, AARG, LYAL, ALA, STY, HIS) - When the substrate is hydrolyzed by the appropriate enzyme (usually an amino acid arylamidase), β -naphthylamine is released. The β -naphthylamine is detected by the addition of p-dimethylaminocinnamaldehyde (in the Peptidase reagent) which creates a complex that is pink to purple in color.

NOTE: The enzyme which cleaves an amino acid β -naphthylamide is an amino acid arylamidase. The term “aminopeptidase” is frequently used as a synonym for arylamidase but actually represents an enzyme with different specificity (an aminopeptidase cleaves the N-terminal amino acid from a peptide). An arylamidase and an aminopeptidase may hydrolyze the same amino acid β -naphthylamide substrate. Thus, the panel is not necessarily measuring different and unique enzymes. A single arylamidase may hydrolyze several amino acid β -naphthylamides.

Carbohydrates (SUC1, SUC2, TRE) - The utilization of sucrose and trehalose results in a pH drop which causes the chlorophenol red to change from purple to yellow. These carbohydrate reactions are selective and may not conform to conventional reactions.

Nitrophenyl Substrates (AGL1, BGL, BGAL, AGL2, BDF, AGAL, NAG, CELL, NGAL) - If the appropriate enzyme is present, the substrate is cleaved releasing ortho- or para-nitrophenol. At alkaline pH, these compounds are yellow. If the reaction occurs at an acidic pH, NaOH must be added after incubation before the results can be read. The wells may be either colorless or pale yellow before the sodium hydroxide is added.

Indoxyl Phosphatase (IDX) - Indoxyl phosphate is cleaved by a phosphatase to release indoxyl. The indoxyl combines with oxygen to form indigo blue which is an insoluble blue or blue-gray precipitate.

Urea (URE) - Urea is split by urease to form ammonia and carbon dioxide. The ammonia (in the form of ammonium carbonate) causes a rise in pH which is detected by phenol red changing from yellow to red.

ORGANISM IDENTIFICATION

The MicroScan Rapid Yeast Biotype Codebook is used for the identification of unknown test organisms. This codebook has been generated from MicroScan's data base for the tests included on the identification panel. The Yeast Codebook is based on a computer analysis of the 27 tests on the Rapid Yeast Identification Panel. The test results are transformed into a 9 digit biotype number for which the codebook lists a species identification and a cumulative relative probability of identification. All possible identifications are printed in descending order of highest probability to a cumulative total of 99.9%.

Should a biotype number occur that cannot be found in the codebook, one should first suspect a procedural error; for example, the biotype number has been added incorrectly or a reaction has been erroneously recorded. If the biotype number is correct, then the possibility of a mixed culture should be checked and the organism should be re-tested, with special emphasis on using a sufficiently turbid cell suspension. If retesting yields the same results, consult the Biotype Lookup Service on the Beckman Coulter website or consult your local distributor.

LIMITATION OF THE PROCEDURE

1. Rapid Yeast Identification Panels are designed for identification of yeast and yeast-like organisms (e.g., *Prototheca*). Care must be taken with organisms with excessive hyphal production.
2. Should more than one species identification be listed for a particular biotype number, supplemental tests may be necessary. Accepted identification procedures of yeast and yeast-like organisms include colonial morphology, growth at elevated temperature, pigment formation and color on phenoloxidase agar. Microscopic morphology on media designed to elicit typical morphology can also aid in identification of these organisms. Refer to appropriate mycology reference procedures.^{5,6,7,8,9,10,11,12,13}
3. Variation in chromogenic results may occur due to excessive over or under inoculation of the system. Each inoculum must be compared to the turbidity standard to achieve the correct concentration of inoculum for optimum performance.

QUALITY CONTROL

The acceptability of the Identification substrates should be checked by testing organisms of known reactions. The results for each reaction for the recommended MicroScan control organisms are found in the Quality Control Chart located in this manual.

Quality Control Chart Rapid Yeast Identification Substrates

	<i>Candida albicans</i> ATCC 66027 ¹	<i>Candida kefyr</i> ATCC 66028	<i>Candida tropicalis</i> ATCC 66029	<i>Cryptococcus albidus</i> ATCC 66030	<i>Cryptococcus neoformans</i> ATCC 66031	<i>Candida glabrata</i> ATCC 66032	<i>Cryptococcus uniguttulatus</i> ATCC 66033
HPR	+	V	V	V	-	V	V
ILE ²	V	V	V	V	-	V	V
PRO	+	V	V	V	-	V	V
TYR	V	V	+	V	-	V	V
GLY	V	-	+	V	V	V	V
GGLY	V	V	+	V	-	V	V
GLAR	+	V	V	V	-	V	V
GLPR ³	V	V	+	V	V	V	V
AARG	-	V	V	+	V	V	V
LYAL	-	V	V	V	V	V	+
ALA	V	V	V	V	-	+	V
STY	-	V	V	V	V	V	+
URE	-	V	V	V	+	V	V
IDX	V	V	+ ⁴	V	-	V	V
HIS	V	V	V	V	-	+	V
SUC1	-	+	V	V	V	V	V
SUC2	-	+	V	V	V	V	V
TRE	-	V	V	V	V	+	V
AGL1	-	V	+	V	V	V	V
BGL	-	+	V	V	V	V	V
BGAL	-	+	V	V	V	V	V
AGL2	V	-	V	+	V	V	V
BDF	-	+	V	V	V	V	V
AGAL ⁵	-	V	V	V	V	V	V
NAG	V	-	V	V	+	V	V
CELL	-	V	V	+	V	V	V
NGAL	+	-	V	V	V	V	V

Shaded QC specifications in the table above will be implemented in LabPro software versions ≥1.90.

1. American Type Culture Collection, Manassas, VA USA
2. A positive result may be obtained with *Rhodotorula rubra*, ATCC 66034.
3. A negative result may be obtained with *Candida lusitanae*, ATCC 66035.
4. Instrument reads may require visual verification.
5. A positive result may be obtained with *Cryptococcus laurentii*, ATCC 66036.



Note: Quality Control organisms are selected to provide positive/negative reactions for all identification substrates. As a result of this, the organism identification may vary from that stated in the QC chart. The acceptability of product performance should be determined by comparison of test results, not identification.

WARRANTY STATEMENT

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

Symbols Key

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1; 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1; 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1; 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1; 5.3.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1; 5.5.5
	In Vitro Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1; 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1; 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1; 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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Tabelle 1-1, Substrate

Substrate	Abk.	Substrate	Abk.
Hydroxyprolin-β-Naphthylamid	HPR	L-Histidin-β-Naphthylamid	HIS
L-Isoleucin-β-Naphthylamid	ILE	Saccharose	SUC1
L-Prolin-β-Naphthylamid	PRO	Saccharose	SUC2 ¹
L-Tyrosin-β-Naphthylamid	TYR	Trehalose	TRE
Glycin-β-Naphthylamid	GLY	p-Nitrophenyl-α-D-Glucopyranosid	AGL1
Glycylglycin-β-Naphthylamid	GGLY	p-Nitrophenyl-α-D-Glucopyranosid	AGL2 ¹
Glycyl-L-Arginin-4-Methoxy-β-Naphthylamid	GLAR	p-Nitrophenyl-β-D-Glucopyranosid	BGL
Glycyl-L-Prolin-4-Methoxy-β-Naphthylamid	GLPR	o-Nitrophenyl-β-D-Galactopyranosid	BGAL

*American Type Culture Collection, Manassas, VA, USA

Substrate	Abk.	Substrate	Abk.
L-Arginyl-L-Arginin- β -Naphthylamid	AARG	p-Nitrophenyl- β -D-Fucopyranosid	BDF
L-Lysin-L-Alanin-4-Methoxy- β -Naphthylamid	LYAL	p-Nitrophenyl- α -D-Galactopyranosid	AGAL
L-Alanin-4-Methoxy- β -Naphthylamid	ALA	p-Nitrophenyl-N-Acetyl- β -D-Glucosamin	NAG
L-Seryl-L-Tyrosin- β -Naphthylamid	STY	p-Nitrophenyl- β -D-Cellobiose	CELL
Harnstoff	URE	p-Nitrophenyl-N-Acetyl- β -D-Galactosaminid	NGAL
3-Indoxylphosphat	IDX		

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Bảng 1-1, Cơ chất

Cơ chất	Chữ viết tắt	Cơ chất	Chữ viết tắt
Hydroxyproline-β-Naphthylamide	HPR	L-Histidine-β-Naphthylamide	HIS
L-Isoleucine-β-Naphthylamide	ILE	Sucrose	SUC1
L-Proline-β-Naphthylamide	PRO	Sucrose	SUC2 ¹
L-Tyrosine-β-Naphthylamide	TYR	Trehalose	TRE
Glycine β-Naphthylamide	GLY	p-Nitrophenyl-α-D-Glucopyranoside	AGL1
Glycylglycine-β-Naphthylamide	GGLY	p-Nitrophenyl-α-D-Glucopyranoside	AGL2 ¹
Glycyl-L-Arginine-4-Methoxy-β-Naphthylamide	GLAR	p-Nitrophenyl-β-D-Glucopyranoside	BGL
Glycyl-L-Proline-4-Methoxy-β-Naphthylamide	GLPR	o-Nitrophenyl-β-D-Galactopyranoside	BGAL
L-Arginyl-L-Arginine-β-Naphthylamide	AARG	p-Nitrophenyl-β-D-Fucopyranoside	BDF
L-Lysyl-L-Alanine-4-Methoxy-β-Naphthylamide	LYAL	p-Nitrophenyl-α-D-Galactopyranoside	AGAL
L-Alanine-4-Methoxy-β-Naphthylamide	ALA	p-Nitrophenyl-N-Acetyl-β-D-Glucosamine	NAG
L-Seryl-L-Tyrosine-β-Naphthylamide	STY	p-Nitrophenyl-β-D-Cellobiose	CELL
Urê	URE	p-Nitrophenyl-N-Acetyl-β-D-Galactosaminide	NGAL
3-Indoxyl Phosphate	IDX		

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REFERENCES

1. Merz, W.G. and G.D. Roberts. 1996. Detection and Recovery of Fungi from Clinical Specimens. P.709-722. In P. R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Yolken (eds.), Manual of Clinical Microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
2. Bobey, D.G., and G.M. Ederer. 1981. Rapid detection of yeast enzymes by using 4-methylumberyl substrates. J. Clin. Microbiol. 13:393.
3. Lee, K.L., M.E. Reza, R.R. Watson, and C.C. Campbell. 1975, Identification of yeast phase of pathogenic fungi by the specificity of their aminopeptidase(S). Sabouraudia 13:132.
4. Watson, R.R. 1976. Substrates Specificities of Aminopeptidase: A Specific Method of Microbial Differentiation, p. 1-14. In J.R. Norris, and D.W. Ribbons (ed.), Methods in Microbiology, Vol. 9, Academic Press, Inc., New York, NY.
5. Callaway, C.J., and L.D. Haley. 1978. Laboratory Methods in Medical Mycology, 4th ed. Centers for Disease control, Atlanta, GA.
6. Emmons, C.W., C.H. Binford, J.P. Utz, and K.J. Kwon-Chung. 1977. Medical Mycology, 3rd ed. Lea and Febiger, Philadelphia, PA.
7. Koneman, E.W., G.D. Roberts, and S.F. Wright. 1985. Practical Laboratory Mycology, 3rd ed. The Williams & Wilkins Co., Baltimore, MD.
8. Kreger-van Rij, N.J.W (ed). 1984. The Yeast-A Taxonomic Study, 3rd ed. Elsevier Science Publishers B.V., Amsterdam.
9. LaRone, D.H. 1987. Medically Important Fungi. A Guide to Identification. 2nd ed. Elsevier Science Publishing Co., New York, NY.
10. McGinnis, M.R. 1980. Laboratory Handbook of Medical Mycology. Academic Press New York.
11. Pore, R.S. 1985, Prototheca Taxonomy. Mycopathol. 90:129-139.
12. Roberts, G.D. 1977. Laboratory identification of fungi, p. 147-160. In E.W. Koneman, and M.S. Britt (ed.), Clinical Microbiology. Health Education Resources, Inc., Bethesda, MD.
13. Warren, N.G., and K. C Hunter. 1996. Candida, Cryptococcus, and Other Yeasts of Medical Importance. P. 723-738. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Yolken (eds.), Manual of Clinical Microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.