
An Indirect Immunofluorescence Test for the Detection of P. jirovecii in Human Clinical Specimens / Test d'immunofluorescence indirecte pour la détection de P. jirovecii dans les échantillons cliniques humains / Prueba indirecta de inmunofluorescencia para la detección de P. jirovecii en muestras clínicas humanas / Indirekter Immunofluoreszenztest (IFFT) zur Bestimmung von P. jirovecii in Humanproben / Test di immunofluorescenza indiretta per la rivelazione di P. jirovecii in campioni clinici umani / Teste de imunofluorescência Indirecta para Detecção de P. jirovecii em Amostras Clínicas de Origem Humana / En indirekte immunofluorescens test for opdagelse af P. jirovecii i kliniske humanprøver / En indirekt immunofluorescenstest för detektion av P. jirovecii i humana kliniska prov / Μια εξέταση έμμεσου ανοσοφθορισμού για την ανίχνευση της P. jirovecii σε ανθρώπινα κλινικά δείγματα

FOR PROFESSIONAL USE ONLY
USAGE RESERVE AUX PROFESSIONNELS
PARA USO PROFESIONAL ÚNICAMENTE
AUSSCHLIESSLICH FÜR BERUFLICHE ZWECKE
SOLO PER USO PROFESSIONALE
EXCLUSIVAMENTE PARA USO PROFISSIONAL
KUN TIL FAGLIGT BRUG
ENDAST FÖR PROFESSIONELL ANVÄNDNING
ΓΙΑ ΕΠΑΓΓΕΛΜΑΤΙΚΗ ΧΡΗΣΗ ΜΟΝΟ

Axis-Shield Diagnostics Limited
The Technology Park, Dundee DD2 1XA
United Kingdom
Tel. : +44 (0) 1382 422000
Fax. : +44 (0) 1382 422088

DISTR.

Bio-Rad
3, boulevard Raymond Poincaré
92430 Marnes-la-Coquette France
Tel. : +33 (0) 1 47 95 60 00
Fax. : +33 (0) 1 47 41 91 33

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1. Stringer JR, Beard CB, Miller RF, Wakefield AE: A New Name (*Pneumocystis jiroveci*) for *Pneumocystis* from Humans. *Emerg Infect Dis* **8(9)**:891-896, 2002.
2. Kovacs JA et al. Diagnosis of *Pneumocystis carinii* Pneumonia: Improved Detection in Sputum with use of Monoclonal Antibodies. *The New Eng J Med*, March 10 1988, **318** No. 10, 499-593.
3. Posner D, Khan FA. Respiratory Infections in AIDS. A Comprehensive Care Approach. *J Respir Dis*, **June 1987**, 83-94.
4. Broaddus CM et al. Bronchoalveolar Lavage and Transbronchial Biopsy for the Diagnosis of Pulmonary Infections in the Acquired Immunodeficiency Syndrome. *Ann Intern Med*, **102**, 747-752, 1985.
5. Mills J. *Pneumocystis carinii* and *Toxoplasma Gondii* Infections in Patients with AIDS. *Rev Infect Dis*, **8**, 1001-1011.
6. Cregan P et al. Comparison of Four Methods for the Rapid Detection of *Pneumocystis carinii* in Respiratory Systems. *J Clin Microb*, **28** No. 11, 2432-2436, 1990
7. Dournon E et al. Diagnosis of *Pneumocystis carinii* pneumonia by non-experts. *The Lancet*, **Jan 14 1989**, 107-108.
8. Midgley J et al. Increased Sensitivity of Immunofluorescence for the Detection of *Pneumocystis carinii*. *The Lancet*, **Dec 23/30 1989**, 1523.
9. Magee JG et al. *Pneumocystis carinii* pneumonia: Detection of Parasites by Immunofluorescence Based on a Monoclonal Antibody. *Medical Laboratory Sciences*, **48**, 235-237, 1991.
10. Midgley J et al. Monoclonal Immunofluorescence Compared with Silver Stain for Investigating *Pneumocystis carinii* pneumonia. *J Clin Path*, **44**, 75-76, 1991.
11. Gill VJ et al. Detection of *Pneumocystis carinii* by Fluorescent-Antibody Stain Using a Combination of Three Monoclonal Antibodies. *J Clin Microbiol*. **25 No. 10**, 1837-1840, 1987.
12. Datry A., Rozenbaum W., Rosenheim M., Danis M., Duflo B. et Gentilini M. - Parasitoses Opportunistes et SIDA. Méthodes de diagnostic. Feuil. Biol., 1985, **XXVI**, 41-46.
13. Petithory J.C., Derouin F., Ardoin F. - Les parasitoses dans le SIDA. Rev. Franç. Lab., Suppl : 15 ème Congrès National des Biologistes des Hôpitaux Généraux, Port Barcarès, Sept. 1986, 21-36.
14. Medrano F, Montes-Cano M, Conde M, de la Horra C, Respaldiza N, Gasch A, Perez-Lozano MJ, Varela JM, Calderon EJ: *Pneumocystis jirovecii* in general population. *Emerg Infect Dis* **11(2)**:245-50,2005.
15. US Department of Health and Human Services. Biosafety in Microbiological and Biomedical Laboratories, Fifth Edition. Washington, DC: US Government Printing Office, January, 2007.
16. Koch M & Heizmann W. Problems in the Detection of *Pneumocystis carinii* by Indirect Immunofluorescence. *Eur. J. Clin. Microbiol. Infect Dis*, **1**, 58-59, 1990

For in vitro diagnostic use / Pour diagnostic in vitro / Para uso diagnóstico in vitro
In Vitro Diagnosticum / Per uso diagnostico in vitro
Para utilizar no diagnóstico *in vitro* / Kun for in vitro diagnostik.
För diagnostisk användning in vitro / Για διαγνωστική χρήση in vitro

Catalogue number / Numero catalogue / Numero de catalogo / Bestellnummer
Numero di catalogo / Número de catálogo Bestillingsnummer / Katalognummer
Αριθμός καταλόγου

Lot / Lot / Lote / Ch.-B. / Lotto
Lote / Parti / Lot / Παρτίδα

60 tests / 60 determinations / 60 pruebas / 60 Bestimmungen
60 tests / 60 testes / 60 tests / 60 tester / 60 προσδιορισμοί

Caution / Avertissez / Adverta / Verwarnen / Cautela
Cuidado! / Advarsel / lakttag försiktighet / Προσοχή

See instructions for use / Voir les consignes d'utilisation
Ver las instrucciones de uso / Gebrauchsinformation beachten
Vedere le istruzioni per l'uso / Ver as instruções de utilização / Se
brugsvejledningen / Se instruktioner för användning / Δείτε τις Οδηγίες χρήσης

Use by / Utiliser avant / Utilizar antes de / Verwendbar bis / Scadenza
Utilizar antes de / Brug før / Använd före / Ημερομ. Λήξης

Store at 2-8°C / Conserver a 2-8°C / Conservar a 2-8°C / Lagerung bei 2-8°C
Conservare a 2-8°C / Armazenar entre 2 e 8°C / Opbevar ved 2-8°C
Förvara vid 2-8°C / Φυλάσσετε στους 2 -8° C

Manufactured by / Fabrique par / Fabricado por / Hergestellt von / Prodotto da
Fabricado por / Fremstillet af / Tillverkad av / Κατασκευάζεται από την

Keep away from light / Conserver à l'abri de la lumière / Mantener alejado de la luz /
Vor Licht schützen / Tenere al riparo dalla luce / Manter ao abrigo da luz / Må ikke
udsættes for lys / Förvara i väg från ljus / Φυλάξτε το μακριά από το φως

Axis-Shield Diagnostics Limited

The Technology Park, Dundee DD2 1XA, United Kingdom.

Tel: +44 (0) 1382 422000, Fax: +44 (0) 1382 422088.

E-mail: shield@uk.axis-shield.com

Website: www.axis-shield.com

Bio-Rad

3, boulevard Raymond Poincaré

92430 Marnes-la-Coquette, France

Tel: +33 (0) 1 47 95 60 00, Fax: +33 (0) 1 47 41 91 33

Website: www.bio-rad.com