MICROBIOLOGY

MicroScan RENOK Rehydrator/Inoculator System Instrument Guide

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9020-7657, Rev. BD

US FMR MR

See each language section for information about copyrights, limitations, and trademarks.

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Revision history

This revision history applies to the Multi-Language document assembly for the RENOK Rehydrator/Inoculator System. See each language for individual revision histories. When subsequent changes affect the information in this document, a new document will be released to the Beckman Coulter website.

9020-7657		
Revision	Date:	Summary of changes:
Rev. BD	2022-06	Minor formatting changes for fr-FR, it-IT, and el-GR.
Rev. BC	2021-09	Additional topics for IVDR compliance.
Rev. BB	2019-06	Minor safety and formatting updates

2022-06

9020-7657, Rev. BD

LC: Multiple Languages (EN, de-DE, es-ES, fr-FR, it-IT, pt-PT, el-GR)

Original instructions (English)

All translations in this multi-language guide are based on Eng. Ref.: 9020-7657

Made in USA
Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 USA
www.beckmancoulter.com

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare, Ireland
Tel: 353 (0) 65 6831100

RENOK Rehydrator/Inoculator System

Intended Purpose

Intended use

The RENOK is a manual pipettor that simultaneously rehydrates and inoculates MicroScan panels (conventional dried panels or designated rapid panels) used in conjunction with disposable MicroScan inoculator sets that consist of a transfer lid (to hold and dispense the inoculum) and a seed trough (to contain the inoculum).

Intended user

This product is intended for laboratory professional use.

Clinical relevance

Not applicable, product is an IVD accessory.

See [RENOK Accessories \(inoculator sets\)](#) (page 7) for additional accessory Intended Purpose information.

Precautions and limitations

Only trained personnel should operate the RENOK Rehydrator/Inoculator.

WARNING

To prevent personal injury, damage to the instrument, or unexpected system performance, observe all safety precautions and warnings.

Always wear Personal Protective Equipment (PPE) as defined by your laboratory policies and follow universal, safe laboratory practices when handling specimens, plates, or any surfaces or instruments that can contact biological specimens.

Consider all materials that come in contact with patient samples as potential biohazards and treat according to local biohazard handling and disposal procedures.

Follow your local or regional regulations for disposing of biohazardous specimens or consumables.


Observe all warnings and cautions on the reagent labels, the safety data sheets (SDS), and hazard warnings when handling any chemicals or flammable liquid materials.

Info for USA only-California Proposition 65: This product may contain chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

For additional MicroScan instrument safety information, including the current hazard symbol glossary, compliance standards, and other related information, see the *Microbiology Instrument Safety Guide*.


Introduction

The RENOK is a manual pipettor that simultaneously rehydrates and inoculates MicroScan panels. It is used in conjunction with disposable MicroScan inoculator sets. An inoculator set consists of a transfer lid (to hold and dispense the inoculum) and a seed trough (to contain the inoculum).

(1)	RENOK	(5)	Partitioned Seed Trough
(2)	Standard Seed Trough	(6)	Spare Rubber Seals
(3)	Transfer Lid	(7)	RENOK Instrument Guide
(4)	Stand		

The following figure shows the parts of the RENOK.

(1)	Release Button	(3)	Pick-Up Levers (2)
(2)	Center Lever	(4)	Base


How RENOK works

The RENOK operates by creating a vacuum when a lever is lifted.

- A chamber within the unit draws in a precisely-controlled amount of air.
- The resulting vacuum draws inoculum into the inoculator transfer lid. No liquid enters the RENOK so contamination of the apparatus is more easily controlled than with direct contact methods.
- The inoculum-filled transfer lid is then placed over a MicroScan panel.
- The inoculum is released, which rehydrates the panel and at the same time inoculates it with the organism to be evaluated.
- The panel is then ready for incubation.

Operation

The following figure shows the proper way to hold the RENOK.

NOTE


When operating the RENOK, follow these usage instructions as well as good ergonomic practices, (for example, breaks, stretching, rotating work activities) to prevent repetitive motion injuries.

Inoculation procedure


Select a level working surface on which to perform the following procedures.

To inoculate a panel:

1. Select and inoculate a seed trough (Standard or Partitioned) according to the appropriate panel procedural manual.
2. If necessary, rock the seed trough gently to distribute the inoculum fluid evenly throughout the seed trough.

3. Place the transfer lid over the seed trough. **Do not attach RENOK to the transfer lid at this time.**
4. Gently tap the transfer lid in all four corners to remove any bubbles. Allow the transfer lid to equilibrate for a minimum of **20 seconds** in the seed trough.
5. Pick up the RENOK from its stand by depressing pick-up levers using thumb and fingers, and place RENOK on top of the transfer lid.

**NOTE**


If a Partitioned Seed Trough is used, make note of its orientation to the transfer lid to ensure the proper inoculation of the corresponding MicroScan panel. See the panel procedural manual for more information.

6. Release the RENOK pick-up levers.

NOTE

For use with Mueller-Hinton broth with lysed horse blood, see the MICroSTREP plus or MICroFAST procedural manual for more information.

7. Fully lift the center lever of the RENOK to draw inoculum into the transfer lid and then release hand pressure. The mechanical catch holds the center lever in place.


8. Pick up the RENOK with the transfer lid attached, holding it level (shaking or rocking can cause spillage), then position the transfer lid onto a MicroScan panel.

NOTE

If inoculating from a Partitioned Seed Trough, ensure that the RENOK and transfer lid are positioned on the MicroScan panel to dispense inoculum according to the appropriate panel procedural manual.


9. Depress the center release button of the RENOK to inoculate panel. There should be a smooth, even release of the inoculum without hesitation or sticking.

10. Dispose of the used transfer lid by positioning the RENOK, with the transfer lid attached, over a biohazard waste container.

Depress the pick-up levers using thumb and fingers to release the transfer lid from the RENOK.

You can also dispose of the seed trough at this time, or you can release the transfer lid onto the seed trough and dispose of the transfer lid and seed trough simultaneously.

11. Place the RENOK on its stand for protection when not in use.

Maintenance

These maintenance procedures describe how to check the dispense volume of the RENOK, replace the rubber seal, and clean the RENOK.

Required tools

The following tools are required to maintain the RENOK.

- Several transfer lids, standard seed troughs, and cover trays
- Gravimetric scale (at least 0.1 gram accuracy)
- Small Phillips screwdriver
- Several vials of Inoculum water with PLURONIC® or Prompt™

Checking the dispense volume

Check the fill volume at least once per month, and whenever low fills are suspected.

To check the dispense volume:

1. Place one cover tray in the gravimetric scale, weigh the tray, and "zero" the scale.
2. With the cover tray on the gravimetric scale, follow the dispensing procedures presented in *Operating Instructions*, but inoculate the tray with PLURONIC.
3. Record the final weight to determine the appropriate inoculation volume. The proper weight of the inoculated cover tray should be within a 10.1 to 12.0 gram threshold.

If the weight is outside the threshold, retest the RENOK using another transfer lid and trough, preferably from another box.

If the weight is still outside the threshold, try replacing the rubber seal.

If the weight is still outside the threshold after you replace the rubber seal, contact your Beckman Coulter representative or distributor.

NOTE


Make sure to follow the dispensing techniques described in this guide to ensure proper fill volumes. Follow these techniques before contacting your Beckman Coulter representative or distributor.

It is still possible that the fill volume was inconsistent due to technique problems, a faulty transfer lid, or a marginal rubber seal. Assure proper technique and replace the rubber seal if needed.


Replacing the rubber seal

Replace the rubber seal if cleaning is required or if you observe dispense volume problems. Complete the following steps to remove part of the RENOK unit body to expose the rubber seal.


1. Use a Phillips screwdriver to remove the two screws from the face of the RENOK and the four screws from the bottom.

2. Pull the body up and away from the base.

3. Remove the rubber seal from the RENOK base and discard in a biohazard waste container. Install a new rubber seal and reassemble the RENOK in reverse order of the preceding instructions.

Cleaning

Clean the RENOK only when contamination occurs. If it is contaminated, the RENOK can be autoclaved.

- Before autoclaving, apply a small amount of stopcock grease to the rubber seal at the base of the unit.
- Use non-fluid autoclave settings.
- **Do not clean by immersion.**

Spare parts list

- Rubber Seal (2)
- Screw Kit (1 set)

Contact your Beckman Coulter representative for spare part ordering information.

Specifications

Nominal dispense volume per well:	115 µL
Well-to-well repeatability, single panel:	± 10 µL
Total dispense volume:	11.04 mL ± 0.96 mL

RENOK Accessories (inoculator sets)

D-Inoculator Sets (for Dried Conventional panels) and R-Inoculator Sets (for Rapid panels) are IVD accessories and include disposable seed troughs, and trough lids (vacuum covers with pipette plates). See [Introduction \(page 3\)](#) for component descriptions.

Intended purpose

All inoculator sets are intended for use only with the MicroScan RENOK by laboratory professionals and are single-use materials. Only use once per inoculum.

Summary and principles

- MicroScan RENOK Disposable D-Inoculator Sets (for Dried Conventional panels) are an in-vitro diagnostic accessory which aids the RENOK in rehydrating and inoculating MicroScan Panels for quantitative and qualitative AST and ID bacterial identification. Refer to MicroScan panel instructions for methods of incubation and results determination which may include manual, semi-automated, or full automation on MicroScan instrumentation.
- MicroScan RENOK Disposable R-Inoculator Sets (for Rapid panels) are an in-vitro diagnostic accessory which aids the RENOK in rehydrating and inoculating MicroScan Panels for qualitative bacterial identification. Refer to MicroScan Rapid Neg ID and Rapid Pos ID panel instructions for methods of incubation and results determination which may include manual, semi-automated, or full automation on MicroScan instrumentation.

See [Inoculation procedure \(page 4\)](#) for information about using inoculator sets with the RENOK.

Inoculator set details

D-Inoculator Sets (for Dried Conventional panels)
• Inoculator set, dried conventional, 240 count (B1013-4) • Inoculator set, dried conventional, 240 count, US (B1013-40) • Inoculator set, dried conventional, 60 count, US (B1013-5)
R-Inoculator Sets (for Rapid panels)
Inoculator set, rapid, 240 count (B1013-6)

Disposal of inoculator sets

Consider trough and lid materials that come in contact with patient samples as potential biohazards and treat according to local biohazard handling and disposal procedures.

- Always follow universal, safe laboratory practices when handling specimens, troughs, lids, panels, or any surfaces or instruments that can contact biological specimens.
- Follow your local or regional regulations for disposing of biohazardous specimens or consumables.

See [Inoculation procedure \(page 4\)](#), Step 10 for more information.

Ordering inoculator sets

Contact your Beckman Coulter representative for inoculator set ordering information. Refer to the part numbers listed above.

Instrument disposal

The Instrument disposal process includes decontaminating and e-waste recycling or disposing of mechanical or electronic components. Bench-top instruments that have reached end-of-life, can be decontaminated and disposed by the customer according to laboratory, local, and regional e-waste regulations.

Instruments must be decontaminated first. This minimizes the risk of biohazards to the environment and to the personnel. See the “Maintenance” topic in this document for more information.

Warranty statement

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents.

The customer is responsible for routine preventative maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.

The RENOK Rehydrator/Inoculator System is guaranteed to be free of mechanical defects for a period of one year. During the warranty period, Beckman Coulter will replace a unit which is found to exhibit defects in materials, workmanship or performance when used according to the procedural manual.

Expected lifetime - RENOK Instruments

The RENOK Rehydrator/Inoculator System has been manufactured and marketed in the U.S. (since the early 1980s) and worldwide (since 1990). Based on historical usage, Beckman Coulter claims that the service life of the product can be expected to exceed 7 years under the defined laboratory conditions identified in the product labeling.

As long as the instrument follows the recommended preventive maintenance schedule, the product will perform within the stated product claims.

Notice to user

Any serious incident that has occurred in relation to this software or device shall be reported to the manufacturer and the appropriate U.S. and country authority (for IVD products) or appropriate authority of the EU Member State (for IVDR products) in which the user and/or the patient is established.

Copyrights

No part of this work covered by the copyrights herein may be reproduced or copied in any form or by any means (graphic, electronic, or mechanical; including photocopying, recording, typing, or information storage and retrieval systems) without the written permission of the publisher.

Trademarks

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

PLURONIC® surfactants. Registered Trademark of BASF Corporation, Parsippany, NJ.

Revision history

9020-7657		
Revision	Date:	Summary of changes:
Rev. BC	2021-08	Additional topics for IVDR compliance.
Rev. BB	2019-06	Minor safety and formatting updates

Future document updates

techdocs.beckmancoulter.com

RENOK Rehydrator/Inokulator System

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