MicroScan

Prompt Procedural Manual

REF B1026-10D

TABLE OF CONTENTS

English	2
Français (FR)	10
Deutsch (DE)	18
Italiano (IT)	27
Español (ES)	35
Português Portugal (PT-PT)	43
Norsk (NO)	51
Ελληνικά (EL)	59
中文 (ZH-CN)	68
Lietuviškai (LT)	75
Magyar (HU)	83
Polski (PL)	92
Čeština (CZ)	101
Slovenčina (SK)	109
한국어 (KO)	117
Türkçe (TR)	125
Русский (RU)	133
Hrvatski (HR)	143
Български (BG)	151
中文 (ZH-TW)	160
Română (RO)	167
Србија (SR)	175
Português Brasil (PT-BR)	183
Nederlands (NL)	191
Tiếng Việt (VI)	200
Македонски (MK)	208
REFERENCES	217

MicroScan

Prompt Procedural Manual

REF B1026-10D

INTENDED USE

The Prompt® Inoculation System-D is used to standardize inocula for microdilution antimicrobial susceptibility and bacterial identification tests.*

INTENDED USER

This product is intended for laboratory professional use.

CLINICAL RELEVANCE

Not applicable, product is an IVD accessory.

SUMMARY AND PRINCIPLE

The Prompt Inoculation System-D is an *in vitro* diagnostic accessory which aids in the inoculum preparation of MicroScan panels for quantitative and qualitative AST and bacterial identification. Refer to MicroScan panel instructions for methods of inoculation, incubation and results determination which may include manual, semi-automated, or full automation on MicroScan instrumentation.

The microdilution procedure for antimicrobial susceptibility testing has provided the clinical microbiologist with a reliable method for obtaining meaningful susceptibility test results. This procedure is used to determine the minimum inhibitory concentration (MIC) of antimicrobial agents and has rapidly gained broad acceptance in the clinical laboratory.¹ Accuracy and reproducibility in the MIC procedure depend on use of defined materials and methods.

One of the important requirements in the MIC procedure is control of the bacterial population of the inoculum within defined limits. This step may be accomplished in two ways:

1. Manual adjustment of the inoculum to match a 0.5 McFarland turbidity standard² followed by appropriate dilution or
2. Incubation to stationary phase in broth culture³ followed by appropriate dilution.

The Prompt Inoculation System-D provides the microbiologist with a method for obtaining standardized inocula while eliminating the need for incubation and turbidity adjustment.^{4,5,6,7} When used according to directions the PROMPT Inoculation System-D will provide the proper inoculum for the microdilution MIC and Biochemical tests for most clinical isolates. MIC results obtained with inocula prepared with the PROMPT Inoculation System-D show greater than 97% agreement with MIC results obtained with inocula prepared according to the CLSI procedure.² Biochemical results obtained with inocula prepared with the PROMPT Inoculation System-D show greater than 97% agreement with Biochemical results obtained with inocula prepared using the turbidity method.

The Prompt Inoculation System-D consists of an inoculation wand and a bottle of diluent. The wand is a polypropylene rod with a breakaway collar that serves as a wiping mechanism. The rod is attached to a stopper. At the tip of the wand is a groove designed to hold a specific number of bacteria. Thirty (30) mL of diluent are provided in the plastic bottle. The wand is touched to several bacterial colonies on a primary isolation plate, wiped, then placed in the plastic bottle. The bacteria are suspended by shaking the bottle. The bacterial suspension is stable for four hours. When the Prompt Inoculation System-D is used within four hours after preparation, 97% of MICs are within $\pm$ one (1) dilution of the MICs determined by the conventional method detailed by CLSI.

The Prompt Inoculation System-D facilitates the MIC inoculum preparation by eliminating:

1) the incubation period, and 2) the need to adjust the inoculum concentration.

WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Colonies must be selected from a FRESH culture plate (<24 hours).
3. The bacterial suspension should be used within four hours of preparation.
4. DO NOT FLAME the plastic inoculation wand under ANY circumstances.
5. To avoid contamination, keep fingers above the ridge of the stopper and on the outer surface of the collar.

*Registered trademark of 3M Company, St. Paul, MN USA

6. Keep wands covered when not in use.
7. Do not use if colonies are smaller than the wand tip.
8. Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures, with special awareness that the inoculated panels contain potentially pathogenic organisms.
9. Single Use Only. All MicroScan consumables should not be re-used unless otherwise noted (e.g. cover trays).

NOTE: After inoculation, the Prompt inoculation vial, wand and collar should be considered potential pathogen carriers. Handle and discard the product accordingly.

STORAGE

Prompt units should be stored at 2-27°C. Avoid direct light. Do not use units beyond the expiration date shown on the package. Inoculated units may be stored for up to four hours at room temperature.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckmancoulter.com/techdocs

EVIDENCE OF DETERIORATION

Prompt bottles should not be used if the bottle or cap is cracked. Inoculation bottles should contain a colorless, clear transparent and odorless liquid. Inoculation wands should be inspected for damage to the tip, wand or collar. Do not use if damage is noted. Contact your Beckman Coulter representative or distributor for further assistance.

SPECIMEN COLLECTION AND PREPARATION

Specimens received in the laboratory should be treated in the usual manner for the preparation of a primary culture plate. Three isolated colonies larger than the wand tip from an eighteen to twenty-four (18-24) hour culture are needed for use with the Prompt Inoculation System-D.

MATERIALS PROVIDED

Prompt Inoculation System-D catalog number B1026-10D contains:*

62 Prompt inoculation wands

60 Prompt inoculation bottles containing 30 mL ($\pm$ 1.0 mL) of stabilized, aqueous PLURONIC® surfactants.**

Prompt Inoculation System-D can be used for preparation of inocula for microdilution MIC testing procedures that use the inoculum suspended in a non-nutritive solution to rehydrate the tray. To use, follow the instructions below:

COMPOSITION

Prompt Inoculation System-D is a proprietary blend of the following components:

Pluronic P-104

Polyvinylpyrrolidone

Potassium phosphate

NOTE: This product does not contain materials of animal origin.

MATERIALS REQUIRED, BUT NOT PROVIDED

Culture media, inoculator tray, quality control organisms and laboratory equipment as required.

*Registered trademark of 3M Company, St. Paul, MN USA

** PLURONIC® surfactants, a registered trademark of BASF Corporation, Parsippany, NJ USA

PROCEDURE

A. Preparation of bacterial suspension

1. Remove the required number of Prompt inoculation bottles from the box.
2. Remove an inoculation wand from the box.
3. Holding the wand tip *perpendicular* to the agar surface, touch 3 isolated colonies as large as or larger than the tip. *Do not* penetrate the agar. *Do not* scrape or drag the tip across the colonies (Figure 1-1).


Figure 1-1

NOTE:

For very small, pin-point colonies, continue incubation of the primary plate until they reach the diameter of the wand tip. If the colony diameter is not likely to reach this size (for example, some streptococci), an alternate method for inoculum preparation should be used.

4. Holding the wand by the handle with one hand, grasp the collar with the other hand and **pull** firmly to break the connection between the collar and wand shaft. **Do Not twist or bend the collar. Keep the tip pointed toward the floor, away from yourself and others** (Figure 1-2).



Figure 1-2

A. Preparation of bacterial suspension

5. Slide the collar slowly down and off the wand shaft and discard the collar (Figure 1-3).

6. While holding the inoculation wand with one hand, pick up a Prompt inoculation bottle.

A0_0821B

Figure 1-3

7. Bend the cap of the bottle sideways until it snaps off (Figure 1-4).

A0_0822B

Figure 1-4

8. Place the inoculation wand into the bottle and press down with a twisting motion to assure a tight seal (Figure 1-5).

A0_0823B

Figure 1-5

A. Preparation of bacterial suspension

9. Shake the bottle vigorously 8 to 10 times to release the bacteria from the wand tip. If the organism is not released from the wand, let the solution sit for 5 minutes and shake again (Figure 1-6).

10. The bacterial suspension should be used within four hours of preparation. If not used immediately after preparation, shake vigorously to resuspend the bacteria just prior to use.


Figure 1-6

B. Addition of Suspension to Seed Tray

1. Remove the inoculation wand from the bottle and discard.
2. Pour the suspension into the seed tray by gently squeezing the bottle with a pumping action (Figure 1-7).
3. Proceed with the appropriate panel procedure.


Figure 1-7

NOTE:

The inoculated Prompt bottle, wand and collar should be considered potential pathogen carriers. Handle and discard appropriately. (See Warning and Precautions for disposal precautions).

QUALITY CONTROL

Refer to appropriate MicroScan panel Procedural Manual.

RESULTS

When used according to directions the PROMPT Inoculation System-D will provide the proper inoculum for the microdilution MIC and Biochemical tests for most clinical isolates. MIC results obtained with inocula prepared with the PROMPT Inoculation System-D show greater than 97% agreement with MIC results obtained with inocula prepared according to the CLSI procedure. Biochemical results obtained with inocula prepared with the PROMPT Inoculation System-D show greater than 97% agreement with Biochemical results obtained with inocula prepared using the turbidity method.²

PERFORMANCE CHARACTERISTICS

Microdilution susceptibility results obtained with inocula prepared by the PROMPT Inoculation System-D were compared to the results obtained with inocula prepared according to the CLSI procedure.² Nine hundred sixteen isolates were tested on MicroScan dried conventional panels. Testing was performed at two sites, using both MicroScan instrumentation and manual verification.

Inocula were prepared according to the Prompt Inoculation System-D and the CLSI procedure. The mean inoculum density for 100 organisms counted immediately after inoculum preparation was 6.9×10^5 cfu/mL. The overall percent of results within +/- 1 dilution was 97%. The MICs for *P. aeruginosa* obtained with aminoglycosides (e.g., amikacin,

gentamicin, tobramycin) with Prompt are 97.5% essential agreement (MICs within +/- 1 dilution of MICs obtained with the CLSI method). Additional drug-specific performance data is on file.

The penicillin and ampicillin MICs for some β -lactamase-producing staphylococci are extremely inoculum dependent. *S. aureus* ATCC 29213 is one of these organisms. Penicillin and ampicillin results should be considered to be in control if this quality control strain shows any penicillin MIC ≥ 0.25 $\mu\text{g/mL}$ and any ampicillin MIC ≥ 0.5 $\mu\text{g/mL}$ (i.e., no upper limit should be applied to these control ranges). When discrepancies for penicillin and ampicillin were seen between the Prompt system and standard-method tests, the organisms involved were β -lactamase producers.

All observed discrepant MICs would be interpreted as resistant (MIC 0.25 $\mu\text{g/mL}$ for penicillin). The Prompt system MICs were invariably higher than the standard MICs, better indicating the resistant nature of the bacteria. This phenomenon has also been noted by Flejzor and Bokkenheuser.⁸ A similar phenomenon was also noted with methicillin-resistant staphylococci and the penicillinase-resistant penicillins.

Staphylococci colony counts may be elevated with the Prompt method and may be greater than the CLSI expected range. Elevated colony counts may adversely affect antibiotic results that are affected by inoculum. Refer to the Limitations section for antimicrobials which are known to be inoculum dependent.

Identification test results obtained with inocula prepared by the the PROMPT Inoculation System-D were compared to the results obtained with inocula prepared using the turbidity method. Six hundred isolates were tested on MicroScan dried conventional panels. Testing was performed at three sites, using both MicroScan instrumentation and manual verification.

Biochemical results obtained with inocula prepared with the PROMPT Inoculation System-D show greater than 97% agreement with Biochemical results obtained with inocula prepared using the turbidity method.

REPRODUCIBILITY

Reproducibility is established by testing multiple replicates under different conditions (e.g., panel lots, days, instruments, users). This may include QC strains, clinical isolates or both. Overall agreement was $\geq 95\%$.

LIMITATIONS

1. The Prompt Inoculation System-D should not be used when the colony size is smaller than the wand tip. Examples of organisms that may not meet the size requirement are *Streptococcus* spp. other than *S. bovis* and *S. agalactiae* (Group B). Picking colonies which are too small will result in under-inoculation and may cause a resistant organism to appear to be susceptible. An alternate method of inoculum preparation should be used for colonies smaller than the wand tip.
2. Some mucoid strains of organisms such as *Klebsiella* spp. or *Pseudomonas* spp. may not adhere to the Prompt wand when attempting to pick the colony. This will be visually apparent. An alternate method of inoculum preparation should be used for such organisms.
3. Refer to appropriate panel instructions for use for additional Prompt System antimicrobial performance limitations.

EXPECTED RESULTS

When used according to directions the Prompt Inoculation System-D will provide the proper inoculum for the microdilution MIC test for most clinical isolates. MIC results obtained with inocula prepared with the Prompt Inoculation System-D show greater than 97% agreement with MIC results obtained with inocula prepared according to the CLSI procedure.²

WARRANTY STATEMENT

The system is covered by and subject to the provisions of the warranty included in your contractual agreement for the system or its reagents. The customer is responsible for routine preventive maintenance procedures. Repairs arising from the failure to perform these maintenance procedures at the indicated time intervals are made at the discretion of Beckman Coulter, and at the customer's expense.


NOTICE TO USER

Any serious incident that has occurred in relation to MicroScan consumables shall be reported to Beckman Coulter and the competent authority of the Member State in which the user and/or the patient is established.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	August 2021	Initial Release

Symbols Key

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1; 5.4.2
	Use-by-date	Indicates the date after which the medical device is not to be used.	ISO 15223-1; 5.1.4
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1; 5.1.5
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1, clause 5.1.6
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.	ISO 15223-1, clause 5.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1; 5.3.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1, clause 5.1.2
	Contains sufficient for <n> tests	Indicates the total number of IVD tests that can be performed with the IVD. (Typically included on reagent kits).	ISO 15223-1; 5.5.5
	<i>In Vitro</i> Diagnostic Medical Device	Indicates a medical device that is intended to be used as an <i>in vitro</i> diagnostic medical device.	ISO 15223-1, clause 5.5.1
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1; 5.3.7
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1; 5.4.3
	CE Mark	European Conformity Regulatory Mark	n/a
	Contents	n/a	n/a
	Contents (package)	Antimicrobial Agent (abbreviation)	n/a
	Contents (package)	Identification Substrate (abbreviation)	n/a
	Safety Data Sheet	Indicates a safety data sheet.	n/a
	Made in USA	n/a	n/a

Symbols Key			
Symbol	Symbol Title	Description of Standard	Standard
	For Export Use Only	n/a	n/a
	Reconstitution Volume	n/a	n/a
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	ISO 15223-1; 5.3.1
	This way up	Indicates the correct upright position	ISO 7000: 0623
	Info For USA Only	Indicates information that is required only for products intended to be sold in the US.	n/a
	Prescription Use Only	Recognized by the US FDA as an alternate to "Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner"	21CFR 801.109(b)(1)
	Stabilized Aqueous PLURONIC® surfactants	n/a	n/a
	Inoculation Wands	n/a	n/a
ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied. Part 1: General Requirements.			

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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REFERENCES

1. Gerlach, E.H. in *Current Techniques for Antibiotic Susceptibility Testing*, A. Balows, ed. (Charles C. Thomas, Publisher, Springfield, Illinois) 1974. pp. 63-76.
2. *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically*. 2018. Approved Standard M07-A11. Clinical and Laboratory Standards Institute, Wayne, PA.
3. Barry, A.L., F. Garcia, and L.D. Thrupp. 1970. An improved single-disc method for testing the antibiotic susceptibility of rapidly-growing pathogens, *Am. J. Clin. Path.* 53:49-158.
4. Wicks, J.H., R.L. Nelson, and G.E. Krejcarek. 1981; A rapid inoculum standardization system: a novel device for standardization of inocula in antimicrobial susceptibility testing. *J. Clin. Microbiol.* 17:1114-1119.
5. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
6. Lund, M.E. and R.W. Hawkinson. 1983. Evaluation of the Prompt Inoculation System for preparation of standardized bacterial inocula. *J. Clin. Microbiol.* 18:84-91.
7. Barry, A.L., R.E. Badal, and R. W. Hawkinson. 1983. Influence of inoculum growth phase on microdilution susceptibility tests. *J. Clin. Microbiol.* 18:645-651.
8. Flejzor, B. and V.D. Bokkenheuser. 1985. Performance of the Prompt System in identification and antimicrobial susceptibility testing of clinical isolates. *J. Clin. Microbiol.* 21:267-268.

Beckman Coulter Ireland Inc.,
Lismeehan, O'Callaghan's Mills
Co. Clare, Ireland +(353) (0) 65 683 1100

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.